

IPS e.max[®] CAD

Technical Guide, May 2022

IPS e.max CAD THE ORIGINAL



IPS e.max CAD



IPS e.max CAD is a tried-and-tested lithium disilicate glass-ceramic block (LS₂) for the fabrication of fixed anterior and posterior restorations.

IPS e.max CAD can be processed in an authorized CAD/CAM machine^[1] in the intermediate crystalline state (≥ 130 MPa). After wet-processing of the block, the restoration is crystallized in a ceramic furnace.^[2]

Proven & reliable

- Since 2005 on your side
- The world's best-selling glass ceramic^[3]
- 97.2% survival rate^[4]

High strength, esthetic &

530 MPa^[7]

Maximum flexibility

- Comprehensive range of restoration types
- Adhesive, self-adhesive or conventional cementation^[8-10]

Trust builds confidence

10 years guarantee

[1] e.g. PrograMill, CEREC/inLab, PlanMill. The complete list is available at www.ivoclar.com.

[2] e.g. Programat CS6, CS4, CS3, CS2, CS, CEREC SpeedFire.

[3] Based on international sales figures.

[4] The survival rate of monolithic IPS e.max CAD posterior crowns over a period of 10 years was evaluated using the Kaplan-Meier method. The failure rate refers to technical failures such as fractures and chipping, R&D Ivoclar, Schaan.

[5] Reich S. et al., Clin Oral Invest, 2014, p. 2171-2178.

[6] Ojeda G et al., Case Reports in Dentistry, 2017, p. 1-6.

[7] Mean biaxial flexural strength, result following more than 10 years of continuous quality testing, R&D Ivoclar, Schaan.

[8] Aslan Y. U et al., Eur. J. Prosthodont. Restor. Dent. 2019 (27), p. 131-140.

[9] Lyann S. K et al., J. Adhes. Dent. 2018 (20), p. 261-268.

[10] Schmitz JH, Beani M, J. Adv. Prosthodont. 2016 (115), p. 678-683.

Solutions

IPS e.max CAD

IPS e.max CAD

Monolithic Solutions



IPS e.max CAD

Abutment Solutions



Types of restoration

Monolithic Solutions

- Veneers
- Inlays
- Onlays (e.g. occlusal veneers, partial crowns)
- Minimally invasive crowns in the anterior and posterior region
- Crowns in the anterior and posterior region
- 3-unit bridges in the anterior and posterior region (up to the second premolar as the terminal abutment)



Abutment Solutions

- Hybrid abutments in the anterior and posterior region as single-tooth restorations
- Hybrid abutment crowns in the anterior and posterior region as single-tooth restorations

Authorized Chairsides^[1]



Manufacturer of CAD-CAM units	System		IPS e.max® CAD				IPS e.max® ZirCAD		IPS Empress® CAD	Tetric® CAD	Tello® CAD	
			Inlays, onlays, veneers and crowns	Three-unit bridges ⁽²⁾	Implant-supported hybrid abutments	Implant-supported hybrid abutment crowns	Crowns	Three-unit bridges	Inlays, onlays, veneers and crowns	Inlays, onlays, veneers and crowns	Temporary crowns and bridges ⁽³⁾	Implant-supported hybrid abutment crowns
	CEREC® MC ⁽⁴⁾	1 block up to 20 mm										
	CEREC® MC X ^(4,5)	1 block up to 55 mm										
	CEREC® MC XL ⁽⁴⁾	1 block up to 85 mm										
	CEREC® Primemill	1 block up to 70 mm										
PLANMECA	PlanMill® 30 S	1 block up to 85 mm										
	PlanMill® 40	1 block up to 55 mm										
	PlanMill® 40 S	1 block up to 85 mm										
	PrograMill® One	Material changer for up to 5 blocks and 45 mm										
	ceramill® motion dts	1 block up to 55 mm										

⁽¹⁾ The range of products may vary from country to country.

⁽²⁾ Up to the second premolar as the terminal abutment

⁽³⁾ With up to two connected pontics

⁽⁴⁾ Wet Milling is possible with new machines since 2013, or with older machines after retrofitting. Dry Milling is possible with construction years 2016 and higher. No retrofitting of older machines.

⁽⁵⁾ B45 (IPS e.max ZirCAD) and B55 (Tello CAD) blocks are possible as of CEREC SW 5.13 (Jan. 2021)

Authorized Labside^[1]



Manufacturer of CAD-CAM units	System	IPS e.max® CAD				IPS e.max® ZirCAD		IPS Empress® CAD	Tetric® CAD	Telio® CAD	
		Inlays, onlays, veneers and crowns	Three-unit bridges ^[2]	Implant-supported hybrid abutments	Implant-supported hybrid abutment crowns	Crowns	Bridges	Inlays, onlays, veneers and crowns	Inlays, onlays, veneers and crowns	Temporary crowns and bridges ^[3]	Implant-supported hybrid abutment crowns
ivoclar	PrograMill® One										
	PrograMill® Dry										
	PrograMill® PM3 / PM5										
	PrograMill® PM7										
	Zenotec® mini										
	Zenotec® select										
	Zenotec® select hybrid										
AMANN GIRRBACH	ceramill® mikro ic										
	ceramill® motion 2 / motion 3										
	ceramill® matik										
Dentsply Sirona	inLab® MC XL										
	inLab® MC X5										

^[1] The range of products may vary from country to country

^[2] Up to the second premolar as the terminal abutment





^[3] With up to two connected pontics

blocks 98,5 mm discs

A male dentist in a white lab coat is shown in profile, smiling and holding a white dental model of a lower jaw with teeth. He is presenting it to a female patient with short blonde hair, who is also smiling and holding a small, single tooth model between her fingers. They are in a dental office with a green chair and a computer monitor in the background.

IPS e.max CAD Monolithic Solutions

The translucency concept

Degree of translucency	Block sizes/ shades	Processing technique				Types of restorations						
		Polishing technique 	Staining technique 	Cut-back technique 	Layering technique 	Occlusal veneer [a]	Thin veneer [a]	Veneer	Inlay, onlay	Partial crown	Crown	3-unit bridge
HT High Translucency	I12, C14, B40, B40 L available in 20 shades ^[a]	✓	✓	✓		✓	✓ ^[d]	✓ ^[d]	✓ ^[d]	✓ ^[d]	✓ ^[d]	
MT Medium Translucency	C14, available in 7 shades	✓	✓	✓		✓	✓ ^[d]	✓ ^[d]	✓ ^[d]	✓ ^[d]	✓ ^[d]	
LT Low Translucency	I12, C14, C16, B32 in 20 shades ^[a]	✓	✓	✓				✓ ^[d]		✓ ^[d]	✓ ^[d]	✓ ^[b]
M0 Medium Opacity	C14, in M0 0 – 4 ^[a]				✓						✓ ^[c]	
I Impulse	C14 in O1 and O2	✓	✓	✓		✓	✓	✓				

^[a] The cut-back technique must not be used when fabricating thin anterior veneers and occlusal veneers.












^[b] Only up to the second premolar as the terminal abutment

^[c] Up to the second premolar

^[d] Max. 2 units if the IPS Speed Tray is used

^[e] The shade range may vary depending on the translucency level/block size or the CAD/CAM machine.

Block concept

		Translucency levels	Shades
HT The minimally invasive		 High translucency Thin veneers, occlusal veneers, veneers, inlays, onlays, partial crowns I12, C14	20 (16 A-D, 4 Bleach)
MT The bright		 Medium translucency Thin veneers, occlusal veneers, veneers, partial crowns, crowns C14	7 (BL2, BL3, BL4, A1, A2, A3, B1)
LT The versatile		 Low translucency Veneers, partial crowns, crowns, bridges, hybrid abutment crowns I12, C14, C16*, A14 & A16 (S&L)*. B32	20 (16 A-D, 4 Bleach)
MO The classical		 Average opacity frameworks on slightly discoloured preparations C14, A14 (S&L)	5 (MO 0, MO 1, MO 2, MO 3, MO 4)
 The opalescent		 Impulse Thin veneers, occlusal veneers, veneers C14	2 (Opal 1, Opal 2)

Block selection



IPS e.max®
Shade Navigation App

Only 5 steps to find the right shade and translucency level

recommendation




*10 shades: BL2, A1, A2, A3, A3.5, B1, B2, C1, C2, D2

Block selection

The IPS e.max Shade Navigation App provides the correct shade and level of translucency for your IPS e.max restoration. All important factors that influence the overall shade design are taken into account.





5 steps to find the correct shade and translucency level



Preparation guidelines

Basic preparation guidelines for all-ceramic restorations

- No angles or edges
- Shoulder preparation with rounded inner edges and/or pronounced chamfer preparation
- The indicated dimensions reflect the minimum layer thickness for IPS e.max CAD restorations.
- The thickness of the preparation edge, particularly for anterior teeth, should be at least 1.0 mm in order to permit optimum milling during CAD/CAM processing.



Preparation guidelines

Minimally invasive posterior crown (adhesive cementation mandatory)	Minimally invasive anterior crown (adhesive cementation mandatory)	Inlay	Onlay	Thin veneer
Veneer	Posterior crown/bridge abutment in the premolar region	Anterior crown/bridge abutment in the anterior region	Occlusal veneer (table top)	Partial crown

Dimensions in mm

Minimum layer thicknesses

Cementation	Mandatory adhesive cementation				Optional adhesive, self-adhesive or conventional cementation			
Types of restorations	Thin veneer	Inlay	Onlay (e.g. occlusal veneer, partial crown)	Minimally invasive crown in the anterior and posterior region	Crown		Bridge	
					Anterior region	Posterior region	Anterior region	Posterior region
Minimum layer thickness IPS e.max CAD – Polishing technique								
Minimum layer thickness IPS e.max CAD – Staining technique								
Incisal/occlusal	0.5	1.0 Depth of the fissures	1.0	1.0	1.5	1.5	1.5	1.5
Circular	0.4	1.0 Isthmus width	1.0	1.0	1.2	1.5	1.2	1.5
Connector dimension	–	–	–	–	–	–	16 mm ² In general: height ≥ width	

The design of bridge connectors should be extended in a vertical direction rather than in a horizontal direction.

Minimum layer thicknesses

Cementation	Mandatory adhesive cementation				Optional adhesive, self-adhesive or conventional cementation				
Types of restorations	Veneer	Inlay	Onlay	Partial crown	Crown			Bridge	
					Anterior region	Premolar region	Molar region	Anterior region	Premolar region
Minimum layer thickness IPS e.max CAD – Cut-back technique									
Incisal/occlusal	0.4	–	–	1.3	0.8	1.0	1.3	0.8	1.0
Circular	0.6	–	–	1.5	1.2	1.5	1.5	1.2	1.5
Connector dimension	–	–	–	–	–	–	–	16 mm ² In general: height ≥ width	
Minimum layer thickness IPS e.max CAD – Layering technique									
Incisal/occlusal	–	–	–	–	0.8	1.0	–	–	–
Circular	–	–	–	–	0.8	0.8	–	–	–
Design type	–	–	–	–	Supporting the tooth shape		–	–	–

Dimensions in mm

Minimum layer thicknesses


IPS e.max CAD is the high-strength component of the restoration and must, therefore, always make up at least 50% of the total layer thickness of the restoration. The total layer thickness of the restoration (depending on the type of restoration) consists of:

Total layer thickness of the restoration	0.8	1.0	1.2	1.5	1.8	2.0	2.5	3.0
Minimum thickness of the IPS e.max CAD framework	0.4	0.5	0.6	0.8	1.0	1.1	1.3	1.6
Maximum layer thickness of the IPS e.max Ceram veneer	0.4	0.5	0.6	0.7	0.8	0.9	1.2	1.4

Dimensions in mm


Example 1:

Minimum layer thicknesses for anterior and posterior crowns in the cut-back technique



Example 2:

Minimum layer thicknesses for anterior and posterior crowns in the layering technique



Model and tooth preparation

Fabricate a working model with removable segments as usual. The manufacturer's instructions of the different CAD/CAM systems must be observed regarding the type of plaster to be used.



Inlay and Onlay




Veneer and
anterior crown



Posterior crown

Processing techniques















Efficiency		Esthetics					
Technique	Polishing technique (self-glaze)	Staining technique Staining technique on the "blue" restoration			Staining technique on the "tooth-coloured" restoration	Cut-back technique	Layering technique
Description	  Polishing of the "blue" restoration, followed by crystallization without individual characterization and glaze.	Variant A   Glazing with IPS e.max CAD Crystall./Glaze Spray on the blue restoration followed by Speed Crystallization (Speed Crystallization and Glaze firing in one step).	Variant B   Staining and glazing with IPS e.max CAD Crystall./Glaze Paste/FLUO on the blue restoration followed by Crystallization and the Stain/Glaze firing in one step.	Variant C   Staining and glazing with IPS e.max CAD Crystall./Glaze Spray on the blue restoration followed by Crystallization and the Stain/Glaze firing in one step.	  Crystallization without the application of materials. Stain/Glaze firing on the tooth-coloured restoration either with IPS e.max CAD Crystall./ or IPS Ivocolor materials	  Crystallization without the application of materials. Incisal firing with IPS e.max Ceram materials. Stain/Glaze firing with IPS Ivocolor materials.	  Crystallization without the application of materials. Dentin/Incisal firing with IPS e.max Ceram materials. Stain/Glaze firing with IPS Ivocolor materials.
Total firing time	<div><div>[1]</div><div><div><div>14:55 min</div><div>or</div><div>25:00 min</div></div></div><div>11:10 min^[2]</div></div>	<div><div><div>14:55 min</div></div></div> <div>11:10 min^[2]</div>	<div><div><div>25:00 min</div></div></div> <div>16:50 min^[3]</div>	<div><div><div>25:00 min</div></div></div> <div>16:50 min^[3]</div>	<div><div><div>37:00 min</div></div></div>	<div><div><div>63:00 min</div></div></div>	<div><div><div>87:00 min</div></div></div>

[1] The respective crystallization and firing parameters are included in the Ivoclar furnaces. Firing times, R&D Ivoclar, Schaan.

[2] Programat CS6, Superspeed crystallization, 11:10 minutes, IPS e.max CAD HT, MT, LT, IPS e.max CAD Crystall./ Glaze Spray or Self Glaze technique, max. two restorations, R&D Ivoclar, Schaan.

[3] Programat CS6, Speed crystallization, 16:50 minutes, IPS e.max CAD HT, MT, LT, IPS e.max CAD Crystall./ Glaze Spray or Paste and Shades/Stains, max. three restorations, R&D Ivoclar, Schaan.

Processing techniques

		Efficiency		Esthetics			
Technique		Staining technique Staining technique on the "blue" restoration			Staining technique on the "tooth-coloured" restoration	Cut-back technique	Layering technique
Description		Variant A	Variant B	Variant C			
	  Polishing of the "blue" restoration, followed by crystallization without individual characterization and glaze.	  Glazing with IPS e.max CAD Crystall./Glaze Spray on the blue restoration followed by Speed Crystallization (Speed Crystallization and Glaze firing in one step).	  Staining and glazing with IPS e.max CAD Crystall./Glaze Paste/FLUO on the blue restoration followed by Crystallization and the Stain/Glaze firing in one step.	  Staining and glazing with IPS e.max CAD Crystall./Glaze Spray on the blue restoration followed by Crystallization and the Stain/Glaze firing in one step.	  Crystallization without the application of materials. Stain/Glaze firing on the tooth-coloured restoration either with IPS e.max CAD Crystall./ or IPS Ivocolor materials	  Crystallization without the application of materials. Incisal firing with IPS e.max Ceram materials. Stain/Glaze firing with IPS Ivocolor materials.	  Crystallization without the application of materials. Dentin/Incisal firing with IPS e.max Ceram materials. Stain/Glaze firing with IPS Ivocolor materials.
Total firing time	<div>[1]</div> <div><div><div>14:55 min</div><div>or</div><div>25:00 min</div></div><div>11:10 min^[2]</div></div>	<div><div>14:55 min</div><div>11:10 min^[2]</div></div>	<div><div>25:00 min</div><div>16:50 min^[3]</div></div>	<div><div>25:00 min</div><div>16:50 min^[3]</div></div>	<div><div>37:00 min</div></div>	<div><div>63:00 min</div></div>	<div><div>87:00 min</div></div>

[1] The respective crystallization and firing parameters are included in the Ivoclar furnaces. Firing times, R&D Ivoclar, Schaan.

[2] Programat CS6, Superspeed crystallization, 11:10 minutes, IPS e.max CAD HT, MT, LT, IPS e.max CAD Crystall./ Glaze Spray or Self Glaze technique, max. two restorations, R&D Ivoclar, Schaan.

[3] Programat CS6, Speed crystallization, 16:50 minutes, IPS e.max CAD HT, MT, LT, IPS e.max CAD Crystall./ Glaze Spray or Paste and Shades/Stains, max. three restorations, R&D Ivoclar, Schaan.

Polishing technique – Finishing and polishing




Smooth out the attachment point of the holder and finish the restoration using suitable grinding instruments.


Pre-polishing with e.g. OptraGloss PP. Prevent the glass-ceramic from overheating.

High-gloss polishing with e.g. OptraGloss HP. Overheating of the glass-ceramic must be prevented.

Polishing technique – Speed crystallization or crystallization



Select the largest possible
IPS e.max CAD Crystallization Pin.




Fill the inside of the crown with
IPS Object Fix Putty or Flow.

Polishing technique – Speed crystallization or crystallization



Press the selected IPS e.max CAD Crystallization Pin deeply into the IPS Object Fix Putty or Flow.








Place the restoration in the centre of the IPS e.max CAD Speed Crystallization Tray or IPS Speed Tray and crystallize using the stipulated firing parameters.

Polishing technique – Speed crystallization or crystallization



Completed, polished and crystallized IPS e.max CAD restoration.

Processing techniques

Efficiency					Esthetics		
Technique	Polishing technique (self-glaze)	Staining technique Staining technique on the "blue" restoration			Staining technique on the "tooth-coloured" restoration	Cut-back technique	Layering technique
Description	  <p>Polishing of the "blue" restoration, followed by crystallization without individual characterization and glaze.</p>	<p>Variant A</p>   <p>Glazing with IPS e.max CAD Crystall./Glaze Spray on the blue restoration followed by Speed Crystallization (Speed Crystallization and Glaze firing in one step).</p>	<p>Variant B</p>   <p>Staining and glazing with IPS e.max CAD Crystall./Glaze Paste/FLUO on the blue restoration followed by Crystallization and the Stain/Glaze firing in one step.</p>	<p>Variant C</p>   <p>Staining and glazing with IPS e.max CAD Crystall./Glaze Spray on the blue restoration followed by Crystallization and the Stain/Glaze firing in one step.</p>	  <p>Crystallization without the application of materials. Stain/Glaze firing on the tooth-coloured restoration either with IPS e.max CAD Crystall./ or IPS Ivocolor materials</p>	  <p>Crystallization without the application of materials. Incisal firing with IPS e.max Ceram materials. Stain/Glaze firing with IPS Ivocolor materials.</p>	  <p>Crystallization without the application of materials. Dentin/Incisal firing with IPS e.max Ceram materials. Stain/Glaze firing with IPS Ivocolor materials.</p>
Total firing time	<p>[1]</p>  <p>11:10 min^[2]</p>	 <p>11:10 min^[2]</p>	 <p>16:50 min^[3]</p>	 <p>16:50 min^[3]</p>	 <p>37:00 min</p>	 <p>63:00 min</p>	 <p>87:00 min</p>

[1] The respective crystallization and firing parameters are included in the Ivoclar furnaces. Firing times, R&D Ivoclar, Schaan.

[2] Programat CS6, Superspeed crystallization, 11:10 minutes, IPS e.max CAD HT, MT, LT, IPS e.max CAD Crystall./ Glaze Spray or Self Glaze technique, max. two restorations, R&D Ivoclar, Schaan.

[3] Programat CS6, Speed crystallization, 16:50 minutes, IPS e.max CAD HT, MT, LT, IPS e.max CAD Crystall./ Glaze Spray or Paste and Shades/Stains, max. three restorations, R&D Ivoclar, Schaan.

Staining technique on the "blue" restoration – Finishing

Thin veneers



Carefully separate the thin veneer from the block using a diamond-coated separating disk.




Smooth out attachment point and finish the restoration.



Finish margins with suitable instruments.

Staining technique – Finishing

Inlays, onlays, crowns




Smooth out attachment point.




Finish the outer surfaces ...

Staining technique – Finishing



... and margins using suitable grinding instruments.



Surface-grind the functional areas in particular with a fine diamond instrument to smooth out the surface structure created by the CAD/CAM procedure.

Staining technique – Finishing



Fitted IPS e.max CAD LT crown ...



... and IPS e.max CAD HT inlay and onlay on the model.

Staining technique – Finishing

Bridges




Bridge after the CAD/CAM process.




Smooth out the attachment point and take proximal contacts into account.

Staining technique – Finishing



Do not use a diamond disc to finish the interdental areas, as this will initiate predetermined breaking points.




Finish interdental areas with fine diamond instruments (grain size 40 – 50 μm), pay attention to the connector dimensions.

Staining technique – Finishing























Check the occlusion, articulation and proximal contact points.



Finish the outer surface, particularly the functional areas of the restoration with a fine diamond instrument to smooth out the surface structure created by the CAD/CAM process.

Processing techniques

Efficiency		Esthetics					
Technique	Polishing technique (self-glaze)	Staining technique on the "blue" restoration		Staining technique on the "tooth-coloured" restoration	Cut-back technique	Layering technique	
Description	  Polishing of the "blue" restoration, followed by crystallization without individual characterization and glaze.	Variant A   Glazing with IPS e.max CAD Crystall./Glaze Spray on the blue restoration followed by Speed Crystallization (Speed Crystallization and Glaze firing in one step).	Variant B   Staining and glazing with IPS e.max CAD Crystall./Glaze Paste/FLUO on the blue restoration followed by Crystallization and the Stain/Glaze firing in one step.	Variant C   Staining and glazing with IPS e.max CAD Crystall./Glaze Spray on the blue restoration followed by Crystallization and the Stain/Glaze firing in one step.	  Crystallization without the application of materials. Stain/Glaze firing on the tooth-coloured restoration either with IPS e.max CAD Crystall./ or IPS Ivocolor materials	  Crystallization without the application of materials. Incisal firing with IPS e.max Ceram materials. Stain/Glaze firing with IPS Ivocolor materials.	  Crystallization without the application of materials. Dentin/Incisal firing with IPS e.max Ceram materials. Stain/Glaze firing with IPS Ivocolor materials.
	Total firing time	 11:10 min ^[2]	 11:10 min ^[2]	 16:50 min ^[3]	 16:50 min ^[3]	 37:00 min	 63:00 min

[1] The respective crystallization and firing parameters are included in the Ivoclar furnaces. Firing times, R&D Ivoclar, Schaan.

[2] Programat CS6, Superspeed crystallization, 11:10 minutes, IPS e.max CAD HT, MT, LT, IPS e.max CAD Crystall./ Glaze Spray or Self Glaze technique, max. two restorations, R&D Ivoclar, Schaan.

[3] Programat CS6, Speed crystallization, 16:50 minutes, IPS e.max CAD HT, MT, LT, IPS e.max CAD Crystall./ Glaze Spray or Paste and Shades/Stains, max. three restorations, R&D Ivoclar, Schaan.

Staining technique (Variant A) – Speed Crystallization and Glaze firing in one step with glazing spray



Glazing with IPS e.max CAD Crystall./ Glaze spray on the "blue" restoration followed by Crystallization.

Only IPS e.max CAD Crystall./Glaze Spray must be used for this processing technique. Crystallization and Glaze firing are performed in one step. The requirements for Speed Crystallization must be observed.

Required materials

- **IPS e.max CAD Crystall./Glaze Spray** is a ready-to-use glazing spray.
- **IPS e.max CAD Speed Crystallization Tray** is a special firing tray for the Speed Crystallization of IPS e.max CAD.




Staining technique (Variant A) – Requirements for Speed Crystallization

IPS e.max CAD	Indication	max. 2 units with IPS e.max CAD Crystall./Glaze Spray
HT High Translucency	Thin veneers, veneers, inlays, onlays, partial crowns, crowns	✓
MT Medium Translucency	Thin veneers, veneers, inlays, onlays, partial crowns, crowns	✓
LT Low Translucency	Veneers, inlays, onlays, crowns	✓
	Hybrid abutment crowns, 3-unit bridges	—
MO Medium Opacity	Hybrid abutments, crowns (framework)	—
I Impulse	Thin veneers, veneers	—

Staining technique (Variant A) – Speed Crystallization (Speed Crystallization and Glaze firing in one step)



Secure the restoration on a suitable IPS e.max CAD Crystallization Pin as described in Variant B and C and apply IPS e.max CAD Crystall./Glaze Spray.



Position the IPS Speed Tray with max. two IPS e.max CAD HT/MT/LT in the furnace and fire using the parameters for the Speed Crystallization firing.



Remove the cool restoration from the set IPS Object Fix Putty or Flow.


Staining technique (Variant A) – Speed Crystallization (Speed Crystallization and Glaze firing in one step)



Remove residue with ultrasound in a water bath ...



... or with the steam jet.

















Do **not** remove residue with Al_2O_3 or glass polishing beads.

Staining technique (Variant A) – Speed Crystallization (Speed Crystallization and Glaze firing in one step)



Completed IPS e.max CAD LT
crown after Speed Crystallization.

Processing techniques


Efficiency		Esthetics					
Technique	Polishing technique (self-glaze)	Staining technique on the "blue" restoration			Staining technique on the "tooth-coloured" restoration	Cut-back technique	Layering technique
Description	  Polishing of the "blue" restoration, followed by crystallization without individual characterization and glaze.	Variant A   Glazing with IPS e.max CAD Crystall./Glaze Spray on the blue restoration followed by Speed Crystallization (Speed Crystallization and Glaze firing in one step).	Variant B   Staining and glazing with IPS e.max CAD Crystall./Glaze Paste/FLUO on the blue restoration followed by Crystallization and the Stain/Glaze firing in one step.	Variant C   Staining and glazing with IPS e.max CAD Crystall./Glaze Spray on the blue restoration followed by Crystallization and the Stain/Glaze firing in one step.	  Crystallization without the application of materials. Stain/Glaze firing on the tooth-coloured restoration either with IPS e.max CAD Crystall./ or IPS Ivocolor materials	  Crystallization without the application of materials. Incisal firing with IPS e.max Ceram materials. Stain/Glaze firing with IPS Ivocolor materials.	  Crystallization without the application of materials. Dentin/Incisal firing with IPS e.max Ceram materials. Stain/Glaze firing with IPS Ivocolor materials.
Total firing time	<div><div></div><div>14:55 min</div></div> or <div><div></div><div>25:00 min</div></div> 11:10 min ^[2]	<div><div></div><div>14:55 min</div></div> 11:10 min ^[2]	<div><div></div><div>25:00 min</div></div> 16:50 min ^[3]	<div><div></div><div>25:00 min</div></div> 16:50 min ^[3]	<div><div></div><div>37:00 min</div></div>	<div><div></div><div>63:00 min</div></div>	<div><div></div><div>87:00 min</div></div>

[1] The respective crystallization and firing parameters are included in the Ivoclar furnaces. Firing times, R&D Ivoclar, Schaan.

[2] Programat CS6, Superspeed crystallization, 11:10 minutes, IPS e.max CAD HT, MT, LT, IPS e.max CAD Crystall./ Glaze Spray or Self Glaze technique, max. two restorations, R&D Ivoclar, Schaan.

[3] Programat CS6, Speed crystallization, 16:50 minutes, IPS e.max CAD HT, MT, LT, IPS e.max CAD Crystall./ Glaze Spray or Paste and Shades/Stains, max. three restorations, R&D Ivoclar, Schaan.

Staining technique (Variant B) – Crystallization and Stain/Glaze Firing in one step with glazing paste



Staining and glazing with IPS e.max CAD Crystall./ Shades, Stain and Glaze Paste/FLUO on the "blue" restoration followed by Crystallization.

The following paragraphs will explain the steps for staining and glazing with IPS e.max CAD Crystall./Shades, Stains and Glaze Paste/FLUO. In this processing technique, Crystallization and the Stain/Glaze firing are performed in one step. Characterizations are applied using IPS e.max CAD Crystall./Shades and Stains.

Required materials

- **IPS e.max CAD Crystall./Shades** are ready-to-use "Dentin" stains in syringes.
- **IPS e.max CAD Crystall./Stains** are ready-to-use intensive stains in syringes.
- **IPS e.max CAD Crystall./Glaze Paste/FLUO** is a ready-to-use glazing paste.
- **IPS e.max CAD Crystall./Glaze Liquid** is a special liquid for mixing with Shades, Stains and Glaze.



Staining technique (Variant B) - Preparation

with
IPS e.max CAD Crystallization Pin



**Partial crowns, anterior crowns,
posterior crowns, bridges**

Use either IPS Object Fix Putty or Flow to secure the restoration on the pin.

without (optional)
IPS e.max CAD Crystallization Pin



Thin veneers, veneers, inlays, onlays

To apply Glaze, Shades, and Stains, secure the restoration

- with an OptraStick,
- with diamond tweezers, or
- directly on the die.

Staining technique (Variant B) - Preparation

Thin veneers – veneers (optional)



Fill the inside of the veneer with IPS Object Fix Flow.



Place the restoration on the IPS e.max CAD Crystallization Pin.



Adapt IPS Object Fix Putty to the pin and restoration margin. Remove any contamination from the outer surface.

Staining technique (Variant B) - Preparation


Crowns – i n l a y s, o n l a y s (o p t i o n a l)



IPS Object Fix Flow and IPS Object Fix Putty.




Select the largest possible IPS e.max CAD Crystallization Pin.



Fill the inside of the restoration with IPS Object Fix Putty or Flow.

Staining technique (Variant B) - Preparation



Press the selected IPS e.max CAD Crystallization Pin deeply into the IPS Object Fix Putty or Flow.




Smooth out displaced IPS Object Fix Putty or Flow with a plastic spatula from the margin towards the support pin so that the pin is secured in the paste and the restoration wall is exactly supported.



Clean off any possible residue adhering to the outer surface of the crown with a brush dampened with water and then dry.

Staining technique (Variant B) - Preparation


Bridge



Select the largest possible IPS e.max CAD Crystallization Pin.



Fill the inner aspects of the crowns of the cleaned bridge with IPS Object Fix Putty or Glow and press the suitable IPS e.max CAD Crystallization Pin into it.



Immediately place the bridge on the IPS e.max CAD Crystallization Tray.

Staining technique (Variant B) - Preparation



Smooth out displaced IPS Object Fix Putty or Flow with a plastic spatula from the margin towards the support pin so that the pin is secured in the paste and the crown margins are exactly supported.



Clean off any possible residue adhering to the outer surface with a brush dampened with water and then dry.

Staining technique (Variant B) – Crystallisation and Stain/Glaze firing in one

Thin veneers, veneers, inlays, onlays
step




Individualized characterization with IPS e.max CAD Crystall./Shades, Stains and Glaze Paste/FLUQ



Secure inlays, onlays on an OptraStick for the application of IPS e.max CAD Crystall./ Glaze Paste/FLUQ, Shades and Stains.


Staining technique (Variant B) – Crystallisation and Stain/Glaze firing in one

Crowns, bridges
step



Apply IPS e.max CAD Crystall./Glaze Paste/FLUO evenly on the blue restoration.

Staining technique (Variant B) – Crystallisation and Stain/Glaze firing in one step



Apply the mixed IPS e.max CAD Crystall./Shades and Stains directly into the unfired IPS e.max CAD Crystall./Glaze Paste/FLUQ

Staining technique (Variant B) – Characterization Guide

Characterization Guide

With the IPS e.max CAD Crystall./Shades and IPS e.max CAD Crystall./Stains, it is possible to characterize the restoration whilst it is in the “blue” state prior to Crystallization. The following Shades and Stains are available for characterization:

IPS e.max CAD Crystall./Shades




IPS e.max CAD Crystall./Stains



Staining technique (Variant B) – Characterization Guide

Depending on the individual patient situation, the characterizations may be applied as follows (example: shade A2):



Minor characterizations on the buccal surface using IPS e.max CAD Crystall./Shade Incisal and Stains.




Occlusal characterizations using IPS e.max CAD Crystall./Shades and Stains.

- Cusp inclinations: Shade Incisal I1
- Fissures: Stains mahogany
- Cusps, marginal ridges: Stains white/cream
- Enhancing the chroma: Stains sunset/copper

Staining technique (Variant B) – Wrong application




Example of IPS e.max CAD Crystall./Shades and Stains applied too thickly



Too thick a layer of IPS e.max CAD Crystall./Shades and Stains

Staining technique (Variant B) – Positioning

- Place the restoration in the centre of the IPS e.max CAD Crystallization Tray.
- Place veneers, inlays and onlays directly on the IPS e.max CAD Crystallization Pin with a small amount of IPS Object Fix Flow.
- A maximum of 6 units can be positioned on the firing tray and crystallized with IPS e.max CAD Crystall./Glaze Paste/FLUO.



Place the restoration in the centre of the IPS e.max CAD Crystallization Tray and fire using the stipulated firing parameters.

Staining technique (Variant B) – Cleaning



Remove the cool restoration from the set IPS Object Fix Putty or Flow.




Remove residue with ultrasound in a water bath ...

Staining technique (Variant B) – Cleaning



... or with the steam jet.



Do **not** remove residue with Al_2O_3 or glass polishing beads.

Staining technique (Variant B) – Finished restoration



IPS e.max CAD LT crown after crystallization.

















Onlay and inlay made of aus IPS e.max CAD HT after Crystallization.

Staining technique (Variant B) – Finished restoration



IPS e.max CAD LT bridge after Crystallization.

Processing techniques

Efficiency		Esthetics					
Technique	Polishing technique (self-glaze)	Staining technique Staining technique on the "blue" restoration			Staining technique on the "tooth-coloured" restoration	Cut-back technique	Layering technique
Description	  Polishing of the "blue" restoration, followed by crystallization without individual characterization and glaze.	Variant A   Glazing with IPS e.max CAD Crystall./Glaze Spray on the blue restoration followed by Speed Crystallization (Speed Crystallization and Glaze firing in one step).	Variant B   Staining and glazing with IPS e.max CAD Crystall./Glaze Paste/FLUO on the blue restoration followed by Crystallization and the Stain/Glaze firing in one step.	Variant C   Staining and glazing with IPS e.max CAD Crystall./Glaze Spray on the blue restoration followed by Crystallization and the Stain/Glaze firing in one step.	  Crystallization without the application of materials. Stain/Glaze firing on the tooth-coloured restoration either with IPS e.max CAD Crystall./ or IPS Ivocolor materials	  Crystallization without the application of materials. Incisal firing with IPS e.max Ceram materials. Stain/Glaze firing with IPS Ivocolor materials.	  Crystallization without the application of materials. Dentin/Incisal firing with IPS e.max Ceram materials. Stain/Glaze firing with IPS Ivocolor materials.
Total firing time	<div><div>[1]</div><div><div><div>14:55 min</div>or<div>25:00 min</div></div></div><div>11:10 min^[2]</div></div>	<div><div><div>14:55 min</div></div><div>11:10 min^[2]</div></div>	<div><div><div>25:00 min</div></div><div>16:50 min^[3]</div></div>	<div><div><div>25:00 min</div></div><div>16:50 min^[3]</div></div>	<div><div><div>37:00 min</div></div></div>	<div><div><div>63:00 min</div></div></div>	<div><div><div>87:00 min</div></div></div>


[1] The respective crystallization and firing parameters are included in the Ivoclar furnaces. Firing times, R&D Ivoclar, Schaan.

[2] Programat CS6, Superspeed crystallization, 11:10 minutes, IPS e.max CAD HT, MT, LT, IPS e.max CAD Crystall./ Glaze Spray or Self Glaze technique, max. two restorations, R&D Ivoclar, Schaan.

[3] Programat CS6, Speed crystallization, 16:50 minutes, IPS e.max CAD HT, MT, LT, IPS e.max CAD Crystall./ Glaze Spray or Paste and Shades/Stains, max. three restorations, R&D Ivoclar, Schaan.

Staining technique (Variant C) – Crystallization and Stain/Glaze Firing in one step with glazing

draw



Staining and glazing with IPS e.max CAD Crystall./ Shades, Stains and Glaze Spray on the "blue" restoration followed by Crystallization.

In this processing technique, IPS e.max CAD Crystall./Glaze Spray is used instead of IPS e.max CAD Crystall./Glaze Paste/FLUO. Crystallization and the Stain/Glaze firing are performed in one step. Characterizations are applied using IPS e.max CAD Crystall./Shades and Stains.

The spray should only be used if IPS Object Fix auxiliary firing paste is adapted up to the margins.

Required materials

- **IPS e.max CAD Crystall./Shades** are ready-to-use "Dentin" stains in syringes.
- **IPS e.max CAD Crystall./Stains** are ready-to-use intensive stains in syringes.
- **IPS e.max CAD Crystall./Glaze Spray** is a ready-to-use glazing spray.



Staining technique (Variant C) – Preparation


with
IPS e.max CAD Crystallization Pin



**Thin veneers, veneers,
inlays, onlays, partial crowns, anterior crowns,
posterior crowns**

Adapt IPS Object Fix Putty or Flow exactly up to the restoration margin so that the glazing spray does not touch the inner areas.


without
IPS e.max CAD Crystallization Pin



Staining technique (Variant C) – Crystallization and Stain/Glaze firing in one step



Extrude IPS e.max CAD Crystall./Shades and Stains from the syringe and mix thoroughly. If required, thin with IPS e.max CAD Crystall./Glaze Liquid.




Apply mixed Shades and Stains directly on the blue restoration.


Staining technique (Variant C) – Crystallization and Stain/Glaze firing in one step



Scuturati bine
flaconul de CAD
Crystall Glaze
Spray



Tineti lucrarea folosind IPS
e.max CAD Crystallization Pin.



Spreiati IPS e.max CAD Crystall./
Glaze Spray direct pe IPS e.max
CAD Crystall nearsa./Shades and
Stains. Spreiati din toate
partile, rotind simultan.

Staining technique (Variant C) – Crystallization and Stain/Glaze firing in one step




Shake the spray can again between individual bursts.




Spray an even layer onto the restoration.

Staining technique (Variant C) – Crystallization and Stain/Glaze firing in one step



Allow the IPS e.max CAD Crystall./Glaze Spray to dry briefly until a whitish layer has formed. If required, spray the restoration again to achieve an even Glaze Spray layer on the IPS e.max CAD restoration.







Place the restoration in the centre of the IPS e.max CAD Crystallization Tray and fire using the stipulated firing parameters.

Staining technique (Variant C) – Example of incorrect Glaze Spray application



Example of incorrect Glaze Spray application

Problem/Cause	Before Firing Application of the Glaze Spray	After Firing Detailed view of the surface
Problem: Not enough Glaze Spray on the restoration Possible cause: <ul style="list-style-type: none"> – Distance between the spray can and the restoration too far – Spraying too short – Spray is not shaken sufficiently – Spray can held not upright during spraying 	 Insufficient application of IPS e.max CAD Crystall./Glaze Spray	 Insufficient gloss or incomplete glaze layer
Problem: Too much Glaze Spray on the restoration Possible cause: <ul style="list-style-type: none"> – Distance between the spray can and the restoration too small – Too much Glaze Spray applied 	 Too much IPS e.max CAD Crystall./Glaze Spray was applied.	 Loss of texture and too glossy surface

Adjustments with IPS e.max CAD Crystall./Add-On

Optional:

- For minor shape adjustments (e.g. proximal contact points), IPS e.max CAD Crystall./Add-On is available.
- The adjustments may be made with both the Crystallization firing or a separate Corrective firing.



Mixing IPS e.max CAD Crystall./Add-On with IPS e.max CAD Crystall./Add-On Liquid to an easy-tocontour consistency.
























Application of the mixed Add-On on the blue restoration before crystallization.



Application of the mixed Add-On on the crystallized restoration.

Processing techniques


Efficiency					Esthetics		
Technique	Polishing technique (self-glaze)	Staining technique Staining technique on the "blue" restoration			Cut-back technique	Layering technique	
Description	  Polishing of the "blue" restoration, followed by crystallization without individual characterization and glaze.	Variant A   Glazing with IPS e.max CAD Crystall./Glaze Spray on the blue restoration followed by Speed Crystallization (Speed Crystallization and Glaze firing in one step).	Variant B   Staining and glazing with IPS e.max CAD Crystall./Glaze Paste/FLUO on the blue restoration followed by Crystallization and the Stain/Glaze firing in one step.	Variant C   Staining and glazing with IPS e.max CAD Crystall./Glaze Spray on the blue restoration followed by Crystallization and the Stain/Glaze firing in one step.	  Crystallization without the application of materials. Stain/Glaze firing on the tooth-coloured restoration either with IPS e.max CAD Crystall./ or IPS Ivocolor materials	  Crystallization without the application of materials. Incisal firing with IPS e.max Ceram materials. Stain/Glaze firing with IPS Ivocolor materials.	  Crystallization without the application of materials. Dentin/Incisal firing with IPS e.max Ceram materials. Stain/Glaze firing with IPS Ivocolor materials.
Total firing time	<div>[1]</div>  11:10 min ^[2]	 11:10 min ^[2]	 16:50 min ^[3]	 16:50 min ^[3]	 37:00 min	 63:00 min	 87:00 min

[1] The respective crystallization and firing parameters are included in the Ivoclar furnaces. Firing times, R&D Ivoclar, Schaan.

[2] Programat CS6, Superspeed crystallization, 11:10 minutes, IPS e.max CAD HT, MT, LT, IPS e.max CAD Crystall./ Glaze Spray or Self Glaze technique, max. two restorations, R&D Ivoclar, Schaan.

[3] Programat CS6, Speed crystallization, 16:50 minutes, IPS e.max CAD HT, MT, LT, IPS e.max CAD Crystall./ Glaze Spray or Paste and Shades/Stains, max. three restorations, R&D Ivoclar, Schaan.

Staining technique on the "tooth-coloured" restoration – Crystallization and Separate Stain/Glaze



Crystallization without application of any materials; separate Stain/Glaze firing either with IPS e.max CAD Crystall./ or IPS Ivocolor materials.

In this processing technique, the IPS e.max CAD restorations are crystallized in a first step without applying any Stains and Glaze materials. Subsequently, the Stain and Glaze firing of the tooth-coloured restoration is conducted.

Required materials

- **IPS e.max CAD Crystall./Shades** are ready-to-use "Dentin" stains in syringes.
- **IPS e.max CAD Crystall./Stains** are ready-to-use intensive stains in syringes.
- **IPS e.max CAD Crystall./Glaze Paste/FLUO** is a ready-to-use glazing paste.
- **IPS e.max CAD Crystall./Glaze Liquid** is a special liquid for mixing with Shades, Stains and Glaze.



or

- **IPS Ivocolor Essences** are intensively shaded stains in powder form.
- **IPS Ivocolor Shades** are ready-to-use stains in jars.
- **IPS Ivocolor Glaze Paste/FLUO, Glaze Powder/FLUO** are glazing materials in paste and powder forms.
- **IPS Ivocolor Mixing Liquids** (allround, longlife) to mix the materials in powder form (Essences, Glaze), as well as to thin paste materials (Shades, Glaze).
- **IPS Ivocolor Essence Fluid** to mix the Essences in powder form to a pasty consistency.



Staining technique – Crystallization



Fill the entire cavity with IPS Object Fix Putty or Flow and place the restoration on the IPS e.max CAD Crystallization Tray.

Remove the crystallization tray from the furnace once the Crystallization program has been completed and allow it to cool.


Staining technique – Cleaning



Remove residue with
ultrasound in a water bath
...




... or with the steam jet.




Do **not** remove residue with
 Al_2O_3 or glass polishing
beads.


Staining technique – Die fabrication with IPS Natural Die Material




Coat the inner surfaces of the ceramic restoration with IPS Natural Die Material Separator and allow it to react for a short time.




Apply IPS Natural Die Material in the restoration and insert the die holder.




Cure in a customary light curing device.




A die made of IPS Natural Die Material provides the optimum basis for true-to-nature all-ceramic restorations.



Staining technique – Stain and Glaze firing with IPS Ivocolor or IPS e.max CAD Crystall./materials



Application of IPS Ivocolor Glaze/FLUO before Glaze firing.
























Application of IPS Ivocolor Essences and Shades before the Stain and Characterization firing.



Completed IPS e.max CAD LT crown after Glaze firing.

Processing techniques


Efficiency						Esthetics	
Technique	Polishing technique (self-glaze)	Staining technique Staining technique on the "blue" restoration			Staining technique on the "tooth-coloured" restoration	Cut-back technique	Layering technique
Description	  <p>Polishing of the "blue" restoration, followed by crystallization without individual characterization and glaze.</p>	<p>Variant A</p>   <p>Glazing with IPS e.max CAD Crystall./Glaze Spray on the blue restoration followed by Speed Crystallization (Speed Crystallization and Glaze firing in one step).</p>	<p>Variant B</p>   <p>Staining and glazing with IPS e.max CAD Crystall./Glaze Paste/FLUO on the blue restoration followed by Crystallization and the Stain/Glaze firing in one step.</p>	<p>Variant C</p>   <p>Staining and glazing with IPS e.max CAD Crystall./Glaze Spray on the blue restoration followed by Crystallization and the Stain/Glaze firing in one step.</p>	  <p>Crystallization without the application of materials. Stain/Glaze firing on the tooth-coloured restoration either with IPS e.max CAD Crystall./ or IPS Ivocolor materials</p>	  <p>Crystallization without the application of materials. Incisal firing with IPS e.max Ceram materials. Stain/Glaze firing with IPS Ivocolor materials.</p>	  <p>Crystallization without the application of materials. Dentin/Incisal firing with IPS e.max Ceram materials. Stain/Glaze firing with IPS Ivocolor materials.</p>
Total firing time	<p>[1]</p>  <p>11:10 min^[2]</p>	 <p>11:10 min^[2]</p>	 <p>16:50 min^[3]</p>	 <p>16:50 min^[3]</p>			

[1] The respective crystallization and firing parameters are included in the Ivoclar furnaces. Firing times, R&D Ivoclar, Schaan.

[2] Programat CS6, Superspeed crystallization, 11:10 minutes, IPS e.max CAD HT, MT, LT, IPS e.max CAD Crystall./ Glaze Spray or Self Glaze technique, max. two restorations, R&D Ivoclar, Schaan.

[3] Programat CS6, Speed crystallization, 16:50 minutes, IPS e.max CAD HT, MT, LT, IPS e.max CAD Crystall./ Glaze Spray or Paste and Shades/Stains, max. three restorations, R&D Ivoclar, Schaan.

Cut-back technique – Crystallization – Incisal firing – Stain/Glaze




Crystallization without the application of materials. Incisal firing with IPS e.max Ceram materials. Stain/Glaze firing with IPS Ivocolor materials.

In the cut-back technique, IPS e.max Ceram Impulse and Incisal materials are layered in the incisal and/or occlusal area of the crystallized, reduced IPS e.max CAD restoration. Subsequently, the Stain and Glaze firing is conducted. The minimal amount of build-up material required means that highly esthetic restorations can be fabricated within only a few working steps.

Cut-back technique – Finishing and Preparation for Crystallization




Fit the milled restoration on the model.
Tooth 11: Anterior crown with cut-back
Tooth 21: Full-contour veneer




Make sure that the minimum layer thickness
of the
restoration is maintained during finishing.

Cut-back technique – Finishing and Preparation for Crystallization



Finish the restoration surface with suitable grinding instruments.




Refrain from designing extreme morphologies with undercuts for mamelons.

Cut-back technique – Finishing and Preparation for Crystallization




Veneer and anterior crown with cut-back on dies.



Bridge with cut-back.

Cut-back technique – Crystallization



Do not place the IPS e.max CAD restoration on metal pins and do not use a honey-combed firing tray.




Slightly overfill the restoration cavity with IPS Object Fix Putty or Flow.

Place the restoration on the IPS e.max CAD Crystallization Tray.

Cut-back technique – Crystallization




Fill the inner aspects of the crowns of the cleaned bridge with IPS Object Fix Putty or Glow and press the suitable IPS e.max CAD Crystallization Pin into it.



Immediately place the bridge on the IPS e.max CAD Crystallization Tray.

Remove the crystallization tray from the furnace once the Crystallization program has been completed and allow the IPS e.max CAD restoration to cool to room temperature.

Cut-back technique – Preparing for veneering



Do **not** remove residue with Al_2O_3 or glass polishing beads.



Remove residue with ultrasound in a water bath ...

Cut-back technique – Preparing for veneering



... or with the steam jet.



Reduced IPS e.max CAD frameworks prepared for veneering.


Cut-back technique – Die fabrication with IPS Natural Die Material




Coat the inner surfaces of the ceramic restoration with IPS Natural Die Material Separator and allow it to react for a short time.



Apply IPS Natural Die Material in the restoration and insert the die holder.



Cure in a customary light curing device.




Cut-back technique – Die fabrication with IPS Natural Die Material





A die made of IPS Natural Die Material provides the optimum basis for true-to-nature all-ceramic restorations.

Cut-back technique – Preparation for veneering



Before veneering, clean the framework under running water or with the steam jet.




Do **not** blast the framework with Al_2O_3 or glass polishing beads.

Cut-back technique – Veneering with IPS e.max Ceram

Firing tray and pins

Use a honey-combed firing tray (Programat® firing tray) and the corresponding support pins to fire the restorations (do not use IPS e.max CAD Crystallization Tray or IPS e.max CAD Crystallization Pins). Round the top edges of the support pin to prevent the restoration from sticking to the pin. Another method of reducing this risk is to cover the pins with platinum foil or a small amount of IPS Object Fix Putty or Flow. Regularly clean the support pins. Do not use contaminated pins.



Cut-back technique – Wash firing – **Variant A with IPS e.max Ceram**




Apply the wash using Incisal and/or Impulse materials ...



... and fire using the stipulated firing parameters.

Cut-back technique – Wash firing – **Variant B with IPS Ivocolor**



Apply the wash using Glaze, Shades, and Essence ...




... and fire using the stipulated firing parameters.

Cut-back technique – Wash firing – **Variant B with IPS Ivocolor**



Finish the full-contour areas of the IPS e.max CAD bridge with diamond grinding instruments and create a true-tonature shape and surface texture.



Apply IPS Ivocolor Glaze on the entire bridge and individually characterize the bridge using IPS Ivocolor Shades and Essences.



The wash must be fired before the actual layering procedure is started. Place the honey-comb firing tray in the furnace and conduct the Wash firing (foundation) using the stipulated firing parameters.

Cut-back technique - Incisal firing



Design the incisal edge using Impulse and Transpa materials.



Complete the layering procedure with Incisal and Transpa materials.



Fire with the firing parameters for the Incisal firing.

Cut-back technique - Incisal firing



Complete the anatomical shape using IPS e.max Ceram Incisal and Transpa materials. Do not veneer the connectors or separate the interdental space.



Place the honey-comb firing tray in the furnace and conduct the Incisal firing using the stipulated firing parameters.



Restoration after the Incisal firing.

Cut-back technique – Finishing and preparing for the Stain and Glaze



Finish the restoration with diamond grinding instruments and give it a true-to-nature shape and surface structure.

Cut-back technique – Stain and Glaze firing
























IPS e.max CAD LT veneer and anterior crown after Glaze firing.



Completed IPS e.max CAD LT bridge after Glaze firing.

Processing techniques


Efficiency							Esthetics
Technique	Polishing technique (self-glaze)	Staining technique Staining technique on the "blue" restoration			Staining technique on the "tooth-coloured" restoration	Cut-back technique	Layering technique
Description	  Polishing of the "blue" restoration, followed by crystallization without individual characterization and glaze.	<p>Variant A</p>   Glazing with IPS e.max CAD Crystall./Glaze Spray on the blue restoration followed by Speed Crystallization (Speed Crystallization and Glaze firing in one step).	<p>Variant B</p>   Staining and glazing with IPS e.max CAD Crystall./Glaze Paste/FLUO on the blue restoration followed by Crystallization and the Stain/Glaze firing in one step.	<p>Variant C</p>   Staining and glazing with IPS e.max CAD Crystall./Glaze Spray on the blue restoration followed by Crystallization and the Stain/Glaze firing in one step.	  Crystallization without the application of materials. Stain/Glaze firing on the tooth-coloured restoration either with IPS e.max CAD Crystall./ or IPS Ivocolor materials	  Crystallization without the application of materials. Incisal firing with IPS e.max Ceram materials. Stain/Glaze firing with IPS Ivocolor materials.	  Crystallization without the application of materials. Dentin/Incisal firing with IPS e.max Ceram materials. Stain/Glaze firing with IPS Ivocolor materials.
Total firing time	<div>[1]</div>  11:10 min ^[2]	 11:10 min ^[2]	 16:50 min ^[3]	 16:50 min ^[3]	 37:00 min	 63:00 min	 87:00 min

[1] The respective crystallization and firing parameters are included in the Ivoclar furnaces. Firing times, R&D Ivoclar, Schaan.

[2] Programat CS6, Superspeed crystallization, 11:10 minutes, IPS e.max CAD HT, MT, LT, IPS e.max CAD Crystall./ Glaze Spray or Self Glaze technique, max. two restorations, R&D Ivoclar, Schaan.

[3] Programat CS6, Speed crystallization, 16:50 minutes, IPS e.max CAD HT, MT, LT, IPS e.max CAD Crystall./ Glaze Spray or Paste and Shades/Stains, max. three restorations, R&D Ivoclar, Schaan.

Layering technique – Crystallization – Dentin-/Incisal firing – Stain/Glaze firing




Crystallization without the application of materials. Dentin/Incisal firing with IPS e.max Ceram materials. Stain/Glaze firing with IPS Ivocolor materials.

In the layering technique, the IPS e.max Ceram layering materials are fired onto the framework made from IPS e.max CAD MO. This enables very individualized design possibilities. The opacity of the IPS e.max CAD MO framework permits the fabrication of highly esthetic restorations on discoloured dies as well as metal core build-ups and Ti abutments.

Layering technique – Finishing and preparation for Crystallization




Place the milled framework on the model and check fit.




Finish the framework with suitable grinding instruments.

Layering technique – Finishing and preparation for Crystallization




Make sure that the minimum layer thicknesses are maintained even after the minor adjustments.



Finish the margins with suitable grinding instruments.

Layering technique – Crystallization




Do **not** place the IPS e.max CAD restoration on metal firing pins and do not use a honey-comb firing tray for crystallization.




Fill the entire cavity with IPS Object Fix Putty or Flow and extend.

Place the firing tray in the furnace and start Crystallization

Layering technique – Preparation for veneering




Do **not** remove residue with Al_2O_3 or glass polishing beads.




Remove any residue with ultrasound in a water bath and/or with steam.



Layering technique – Preparation for veneering



Check marginal areas and slightly finish, if necessary.




Finished IPS e.max CAD MO framework.


Layering technique – Veneering with IPS e.max Ceram

Firing tray and pins

Use a honey-comb firing tray (Programat firing tray) and the corresponding support pins to fire the restorations (do not use an IPS e.max CAD Crystallization Tray or IPS e.max CAD Crystallization Pins). Round the top edges of the support pin to prevent the restoration from sticking to the pin. Another method of reducing this risk is to cover the pins with platinum foil or a small amount of IPS Object Fix Putty or Flow. Regularly clean the support pins. Do not use contaminated pins.



Layering technique – Wash firing – **Variant A with IPS e.max Ceram**




With ideal space, apply the wash layer with the required IPS e.max Ceram Deep Dentin, Dentin, Transpa Incisal and/or Impulse materials ...




... and fire using the stipulated firing parameters.

Layering technique – Wash firing – **Variant B with IPS Ivocolor**




Apply the wash using Glaze, Shades, and Essence ...




... and fire using the stipulated firing parameters.

Layering technique – 1st Dentin/Incisal firing



Contour the tooth shape with Dentin material.



Cut-back and build-up of the incisal area with incisal extension.



Design the incisal third using Impulse materials.

Layering technique – 1st Dentin/Incisal firing



Complete the layering procedure with Incisal and Transpa materials.



Subsequently, the restoration is fired using the firing parameters for the 1st Dentin/Incisal firing.

Layering technique – 2nd Dentin/Incisal firing (Corrective firing)



Compensate for the shrinkage using Dentin, Transpa and Incisal materials.



Subsequently, the restoration is fired using the firing parameters for the 2nd Dentin/Incisal firing.

Layering technique – Stain and Glaze firing








Completed, stained and glazed IPS e.max CAD MO restoration.

Cementation possibilities

	Adhesive cementation e.g. Variolink® Esthetic, Multilink® Automix	Self-adhesive cementation e.g. SpeedCEM® Plus	Conventional cementation e.g. Vivaglass® CEM, ZirCAD® Cement
Preparation requirements	Non-retentive preparation	Retentive preparation	Retentive preparation
Veneers	✓	–	–
Inlays, onlays (e.g. occlusal veneers, partial crowns)	✓	–	–
Minimally invasive crowns	✓	–	–
Crowns	✓	✓	✓
Three-unit bridges up to the second premolar as the terminal abutment	✓	✓	✓

Cementation possibilities

Firing										Polishing			
IPS e.max® CAD						IPS e.max® ZirCAD		IPS Empress® CAD		Tetric® CAD		Telio® CAD	
													
Lithium disilicate glass-ceramics (LS ₂)						Zirconium oxide ceramics (ZrO ₂)		Leucite glass-ceramics		Composite		PMMA	
Flexural strength	530 MPa ⁽¹⁾					MT Multi: 850 MPa ⁽²⁾ LT: 1200 MPa ⁽²⁾		185 MPa ⁽¹⁾		272 MPa ⁽²⁾		135 MPa ⁽²⁾	
Types of restorations	Veneers, inlays, onlays, minimally invasive crowns (min. 1.0 mm)		Crowns (min. 1.5 mm), 3-unit bridges up to the 2 nd premolar			Crowns, 3-unit bridges		Veneers, inlays, onlays, crowns		Occlusal veneers, veneers, inlays, onlays, crowns		Temporary crowns and bridges	Long-term crowns and bridges
Cementation method	adhesive		self-adhesive		conventional	adhesive	self-adhesive/conventional	adhesive		adhesive		temporary	adhesive
Blasting	–					Al ₂ O ₃ 25–70 µm, 1 bar or Al ₂ O ₃ 70–110 µm, 1.5 bar		–		Al ₂ O ₃ 25–70 µm, 1 bar or Al ₂ O ₃ 70–110 µm, 1.5 bar		Al ₂ O ₃ 25–70 µm, 1 bar or Al ₂ O ₃ 70–110 µm, 1.5 bar	
Etching	Option 1: Agitate Monobond Etch&Prime® for 20 s and allow it to react for another 40 s	Option 2: 20 s with IPS® Ceramic Etching Gel	Option 1: Agitate Monobond Etch&Prime® for 20 s and allow it to react for another 40 s	Option 2: 20 s with IPS® Ceramic Etching Gel	20 s with IPS® Ceramic Etching Gel	–		Option 1: Agitate Monobond Etch&Prime® for 20 s and allow it to react for another 40 s	Option 2: 60 s with IPS® Ceramic Etching Gel	–		–	
Conditioning	60 s with Monobond® Plus		60 s with Monobond® Plus		–	60 s with Monobond® Plus	–	60 s with Monobond® Plus		20 s with Adhese® Universal	30 s with Multilink® Primer A+B	–	2–3 min with SR® Connect
Cementation system	Variolink® Esthetic, Multilink® Automix ⁽³⁾		SpeedCEM® Plus		Vivaglass® CEM	Multilink® Automix	SpeedCEM® Plus, Vivaglass® CEM	Variolink® Esthetic, Multilink® Automix ⁽³⁾		Variolink® Esthetic	Multilink® Automix ⁽³⁾	Telio® Link	Variolink® Esthetic, Multilink® Automix

⁽¹⁾ Average biaxial flexural strength, Outcome after more than 10 years of ongoing quality testing, R&D Ivoclar, Schaan.

⁽²⁾ Typical mean value of biaxial flexural strength, R&D Ivoclar, Schaan.

⁽³⁾ Not recommended for veneers.



more information
www.cementation-navigation.com



Manufacturer
Ivoclar Vivadent AG
Benderstrasse 2
9494 Schaan/Liechtenstein
www.ivoclar.com

















Date information prepared: 2022-04-08, Rev. 2
702850/DE/EN

Shade Combination Tables

IPS e.max Crystall./Shades, Stains

To be used on "blue" and "tooth-coloured" IPS e.max CAD restorations
































A–D tooth shade	BL1	BL2	BL3	BL4	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
IPS e.max CAD Crystall./ Shade																				
		0					1				2				3				4	
IPS e.max CAD Crystall./ Shade Incisal																				
			I1						I2		I1						I2			
IPS e.max CAD Crystall./ Stains																				
			white			cream		sunset		copper		olive		khaki		mahogany				

IPS Ivoclar Shades, Essences

To be used on "tooth-coloured" IPS e.max CAD restorations



A–D tooth shade	BL1	BL2	BL3	BL4	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4	
IPS Ivocolor Shade																					
IPS Ivocolor Shade Incisal																					
IPS Ivocolor Essence																					
																					
					E 21 basic red			E 22 basic yellow			E 23 basic blue										

Crystallization and Firing Parameters

Crystallization MO, Impulse, LT, MT, HT

with or without the application of IPS e.max CAD Crystall./ materials

[illegible]

Crystallization LT, MT, HT

with or without the application of IPS e.max CAD Crystall./ materials

[illegible]

Crystallization and Firing Parameters

Speed crystallization (observe the block concept)

Max. 2 units *with or without* the application of **IPS e.max CAD Crystall./Glaze Spray** on an IPS Speed Tray

[illegible]

Corrective firing/stain firing/glaze firing

with IPS e.max CAD Crystall./ materials

[illegible]

Firing Parameters

Firing parameters for the staining technique
with IPS Ivocolor Shade, Essence, Glaze



	Stand-by temperature B [°C]	Closing time * S [min]	Heating rate t↗ [°C/min]	Firing temperature T [°C]	Holding time H [min]	Vacuum 1 V1 [°C]	Vacuum 2 V2 [°C]	Long-term cooling ** L [°C]	Cooling rate tl [°C/min]
Stain and Glaze firing	403	IRT/ 6:00	60	710	1:00	450	709	0	0

* IRT normal mode

** Note: If the layer thicknesses exceed 2 mm, long-term cooling to 500 °C is required.

Note: Due to their geometry, the restorations may feature varying layer thicknesses. When the objects cool after the firing cycle, the different cooling speeds in the areas with different thicknesses may result in a build-up of internal tension. In the worst case, these internal tensions may result in fracture of the restoration. By using slow cooling (long-term cooling L), these tensions can be minimized. For monolithic restorations (staining technique) featuring layer thicknesses of more than 2 mm, long-term cooling L must be used.

Firing Parameters

Firing parameters for the cut-back and layering technique
with IPS e.max Ceram/IPS Ivocolor Shade, Essence, Glaze



	Stand-by temperature B [°C]	Closing time * S [min]	Heating rate t ↗ [°C/min]	Firing temperature T1 [°C]	Holding time H1 [min]	Heating rate t ↗ [°C/min]	Firing temperature T2 [°C]	Holding time H2 [min]	Vacuum 1 11 12 [°C]	Vacuum 2 21 22 [°C]	Long-term cooling L [°C]	Cooling rate tl [°C/min]
Wash firing (foundation)	403	IRT/ 04:00	90	650	00:00	20	730	02:00	400/650	650/729	0	0
1 st Dentin and Incisal firing	403	IRT/ 04:00	90	650	00:00	20	730	02:00	400/650	650/729	0	0
2 nd Dentin and Incisal firing	403	IRT/ 04:00	90	650	00:00	20	730	02:00	400/650	650/729	0	0
Stain firing using IPS Ivocolor	403	IRT/ 06:00	60	710	01:00	–	–	–	450	709	0	0
Glaze firing using IPS Ivocolor	403	IRT/ 06:00	60	710	01:00	–	–	–	450	709	0	0
Add-On with Glaze firing	403	IRT/ 06:00	60	710	01:00	–	–	–	450	709	0	0
Add-On after Glaze firing	403	IRT/ 06:00	50	700	01:00	–	–	–	450	699	0	0

* IRT normal mode

IPS e.max CAD – Because it **WORKS**

The world's most used glass-ceramic.

Making
People
Smile

Die Inhalte dieser Präsentation sind urheberrechtlich geschützt.
Ohne vorherige ausdrückliche schriftliche Zustimmung der Ivoclar Vivadent AG dürfen diese Inhalte weder ganz noch auszugsweise vervielfältigt, verbreitet, öffentlich vorgeführt, öffentlich zugänglich gemacht, gesendet oder durch Bild- oder Tonträger öffentlich wiedergegeben werden.

The contents of this presentation are protected by copyright.
Reproduction, distribution, public performance, making available to the public, broadcasting, cable retransmission, or communication by video or audio recordings of these contents, whether in whole or in parts, is subject to prior express written permission by Ivoclar Vivadent AG.