**anaxBLEND**

**MANUAL**

**anaxDENT**

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**anaxBLEND** is a unique, new facing system comprising pastes of differing hardness and a variety of liquid materials (up to 75% porcelain fillers). This novel fusion allows faster and easier processes, no matter if you choose paste-like or flowing materials.

**anaxBLEND** can not only be used for layering and modelling, it can even be overpressed or injected into anaxform flasks. Rich, opaque dentin colors and different incisal and transparent materials ensure results with utmost esthetics.

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### **Stability and Esthetics**



Thanks to the concerted mixture of fillers, composites and colors, some of the advantages of anaxblend sound like high abrasion resistance, excellent elasticity, and brilliant esthetics. The high bending strength protects from fracture and flake-off. The excellent polishability and the thus resulting smooth surface provide the restorations made of anaxblend with a great natural gloss. anaxblend gains safety and the important decisive advantage from its low affinity to plaque due to the homogeneous structure. The anax bonding systems Metal Bonder, Zircon Bonder and Bonding Fluid offer an excellent compound.

### **Features**

The microfilled, low viscous, paste-like facing material anaxblend is suitable for the application of the dentin structure as well as for the conventional modeling of the complete facing whereas anaxblend Flow (highly viscous) was specially created for the brushing technique and for the injection technique in combination with the anaxform flasks. Dentin, incisal and transparent materials of the anaxblend Flow series can freely, i. e. like porcelain, be applied with a brush.

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


**Figure 01 - 02** Using Matrix Form 70 to make silicone rims from the wax-up (teeth with gingival parts and occlusal stops in the opposing arch). The rims are used to model the framework neatly and with reduced anatomy. (For your information: The facing material requires a minimum layer thickness of 1.2mm to meet the esthetical demands.) In case of metals, retention beads offer the according retention. **Figure 03 - 04** The finished metal framework after sandblasting with a 110-250  $\mu\text{m}$  abrasive, steam cleaning, and drying with oil-free compressed air. Please apply anaxblend Metal Bonder to achieve a good compound between the metal and the facing material. In case of zirconia frameworks, please use anaxblend Zircon Bonder. **Figure 05 - 06** Thinly apply the first layer of paste opaquer (max. 80  $\mu\text{m}$ ) with a disposable brush, then polymerize. Repeat this procedure 2 to 3 times until the desired opacity is achieved. In case of anaxgum Opaquers, this procedure must be repeated 3 to 4 times at least. With this material, only very thin layers should be applied since these opaquers are strongly pigmented and allow only a very low penetration depth of the polymerization light. After the polymerization, the surface should show an eggshell gloss. If a layer does not set properly, the material was applied too thickly. Please do not contaminate the opaquer surface by touching.

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**Figure 07 - 08** Opaque Dentin can be used to simulate the cervical area or to control the color intensity. Please polymerize in steps! Then apply Dentin Paste and let it fade to the incisal edge. Thanks to the high opacity and the intensive chroma of the dentin, a minimum layer thickness of 0.8 mm is sufficient to obtain a sufficient color. **Figure 09** Layer anaxblend incisal / transparent Flow in the incisal area as required. **Figure 10** Mamelons may be simulated by a mixed application of make-in MI-63 or anaxblend Colors brown with Transpa CL-O. Incisal lucency can be imitated with Enamel S57 Flow. **Figure 11** Step-by-step polymerization and control of the layer thickness by means of the rim. **Figure 12 - 13** Layer technique in the incisal area with different shades of the incisal pastes. Compared to the other porcelain materials, the incisal pastes distinguish by an increased abrasion resistance and hardness. Therefore, they are also suitable to have composite inlays and full crowns made by the dental technician. In this case, the rim can be used to verify the shape and the dimensions as well.

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**Figure 14 - 17** After removing the rim, the accurate course of the papillae becomes visible. Please do not unscrew the restoration for layering. Evenly spread anaxgum Gingiva dark pink as basic material from the course of the papillae up to the model. Supplement the restoration from interalveolar up to the mucolabial fold until reaching the final shape. Prepare the alveolar tubercles with anaxgum orange pink. Finish the margins of the papillae with anaxgum Gingiva light pink and brush it down in apical direction. Use anaxblend Dentin A2 Flow to build up the alveolar tubercle. Build up the course of the root of a canine with anaxblend Dentin B3 Flow to match the richer shade of canines. **Figure 18 - 22** Finished alveolar tubercle: Use anaxblend Dentin A2 Flow to produce the approximal connection of the ridge. The gingival surface gets a slight lucency by applying anaxblend Enamel S57 Flow below the course of the sulcus. Preparations to present the papillae and gingival stippling. Then polish it with pasta grigia II and a pointed polishing brush inserted in a handpiece. Add a high-gloss polish using a buff in a handpiece.

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
**Figure 01 - 03** Wax-up on the finished telescopic structure. **Figure 04 - 05** Fabrication of the silicone model by means of the anaxform glass flask by filling Matrix Form 60 into the bottom of the flask, injecting Matrix Flow 70 clear / 56 violet from basal into the telescopic crowns (dies), and embedding it into the still soft Matrix Form 60 putty.

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



**Figure 06 - 08** Trim the silicone model up to the cervical margin of the crown. Separate the silicone model (anaxdent Divorce) and let it dry.  
**Figure 09 - 12** Adapt injection canals (approx. 1.2mm) and adjust it in the center with the flask lid put on. Mix Matrix Cast Clear under vacuum and let it flow into the closed glass flask as a thin jet. Completely vulcanize it in a pressure pot (at 2 bar for 25 minutes) until clear.

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

**Figure 13 - 15** Divest the framework from the duplication, remove wax, clean, degrease, sandblast, steamclean, and then condition it (Metal Bonder / Zircon Bonder). **Figure 16** Thinly apply the first layer of paste opaquer (max. 80  $\mu\text{m}$ ) with a disposable brush, then polymerize. Repeat this procedure 2 to 3 times until the desired opacity is achieved. When using the pink anaxgum opaquer, this procedure is to be repeated 3 to 4 times. With this material, only very thin layers should be applied since this opaquer is strongly pigmented and allows only a very low penetration depth of the polymerization light. After the polymerization, the surface should show an eggshell gloss. If a layer does not set properly, the material was applied too thickly. Please do not contaminate the opaquer surface by touching. **Figure 17 - 19** The closed flask ready for the inject technique with anaxblend Flow.

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**Figure 20 - 21** Hold the flask in an angle of 45°, apply Dentin Flow 3g / 20g big block syringe, and inject into the last injection canal from bottom to top. Constantly press material out of the syringe until the hollow is completely filled. For polymerization times see the table on page 31. **Figure 22 - 24** Remove the injection canals and reduce the dentin body to create space for the incisal material. **Figure 25 - 26** After cutting back the dentin body, the surface must extensively be wetted with anaxblend bond LC to renew the dispersion layer (adhesion). Then polymerize for 60 seconds in a light-curing unit. Color the dentin body with the anaxblend Colors stains and add effect materials as required. Then polymerize.

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**Figure 27 - 29** Reposition the restoration in the flask, close the lid, add the according anaxblend flow enamel, and inject into the last injection canal from bottom to top until the hollow is completely filled. Divest after the polymerization. For polymerization times see the table on page 31. **Figure 30** Remove the injection canals and finish the restoration. **Figure 31** Wax modellation for the Acryline pink acrylic basis. **Figure 32** Condition Acryline pink with bond LC after being polymerized and processed to provide adhesion (dispersion layer) between the PMMA acrylic Acryline to the composite anaxgum. Individualize the gingival parts with anaxgum Pastes / Flow and finish. Use a brush or buff in combination with pasta grigia II to provide a high-gloss polish.




	max. Schicht / max. layer	Heraf flash / UniXs	Spectra LED	Triad 2000	Spectramat	Solidite V
Opaquere	max. 80µm	90 sec	30 sec	6 min	2 min	2 min
Dentin	2,00 mm	180 sec	3 min	6 min	5 min	5 min
Enamel	2,00 mm	180 sec	3 min	6 min	9 min	9 min
Transpa	5,00 mm	180 sec	3 min	6 min	5 min	5 min
anaxgum	2,00 mm	180 sec	3 min	6 min	5 min	5 min

	anaxBLEND Enamel Paste	anaxBLEND Paste	anaxBLEND Flow	anaxGUM Paste
Härte / hardness	574 MPa	405 MPa	365 MPa	405 MPa
E-Modul / modulus of elasticity	10.040 MPa	6.890 MPa	6.530 MPa	6.890 MPa
Biegefestigkeit / bending strength	143 MPa	94 MPa	140 MPa	94 MPa
Druckfestigkeit / compressive strength	455 MPa	365 MPa	360 MPa	365 MPa
Wasserlöslichkeit / water solubility	0,7 µg / mm <sup>3</sup>	2,6 µg / mm <sup>3</sup>	2,8 µg / mm <sup>3</sup>	2,6 µg / mm <sup>3</sup>
Wasseraufnahme / water absorption	15,2 µg / mm <sup>3</sup>	19,3 µg / mm <sup>3</sup>	21,1 µg / mm <sup>3</sup>	19,3 µg / mm <sup>3</sup>

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The perfect bonding compound between zirconia and composites / PMMA.

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Liquid bonding agent for the adhesion of metals to acrylics and composites.

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For the chemical compound of composite to composite / PMMA / Bis-Acryl.

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Light-curing modelling liquid for the application on anaxblend.

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The effective silicone separating agent.

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Gel to prevent the inhibitionlayer at the final curing.

## ZUBEHÖR ACCESSOIRES

**Küvetten / Flasks**

Glassküvetten Set small und big / Universalküvette  
Flask set out of glass, small and big / universal metal flask

**Silikone / Silicones\***

MATRIX Form 60 / Form 70  
MATRIX Flow 70 clear / 56 violett  
MATRIX Cast clear  
\*see Acrylic Brochure

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Dentin Flow 3g: A1 / A2 / A3 / A3,5 / B3 / C2 / D3

Enamel Flow 3g: S57 / S58 / S59 / S60 / Transpa NT / bond LC 10ml / model

LC 10ml / anaxbrush No. 4 / Cover Gel 3g / Pasta Grigia II / anaxrack

**anaxBLEND XL KIT Paste**

Opaquer 3g: A1 / A2 / A3 / A3,5 / B2 / B3 / C2 / D3

Dentin Paste 4g: A1 / A2 / A3 / A3,5 / B3 / C2 / D3

Enamel Paste 4g: S57 / S58 / S59 / S60 / Transpa NT / bond LC 10ml /

model LC 10ml / anaxbrush No. 4 / Cover Gel 3g / Pasta Grigia II / anaxrack

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

**Framework materials**

Please model the framework neatly and anatomically in the common technique. The facing material requires a minimum layer thickness of 1.0mm to meet the esthetical demands. The condition for an accurate connection between the framework and the facing material is a chamfer with a sharp margin.

**Precious and non-precious alloys and titanium**

If palladium-free casting alloys are used, the minimum layer thickness for telescopic cases of 0.4 mm may not be under-run (see indications of the alloy manufacturer). We recommend the use of an alloy with a modulus of elasticity of 85GPa at least. Sandblast the surface with 110 µm aluminum oxide beads (at max. 2bar) (in a jet angle of 45°, with a distance of 3cm). Sandblast non-precious alloys at 3bar. Then clean the restoration with a steamcleaner and dry it with oil-free compressed air.

**Chemical compound**

anaxdent Metal Bonder

anaxdent Metal Bonder serves as bonding agent between metals to PMMA and composites. Metal Bonder corresponds with the related monomers of the anaxblend opaquer to achieve an excellent adhesion of both components. Apply the Bonder with anax brush No. 2 once or twice and let it air-dry for 10 seconds. If required, add a second / third layer. Please cover the Metal Bonder bottle immediately after being used. Optimal bond strength is achieved by applying anaxblend Opaquer without contamination due to contact and without long delays followed by the application of the facing acrylic.

anaxdent Zircon Bonder

anaxdent Zircon Bonder serves as bonding agent between zirconia and composites / PMMA. Zircon Bonder corresponds with the related monomers of the anaxblend opaquer to achieve an excellent adhesion of both components. Apply the Bonder with anax brush No. 2 once or twice and let it air-dry for 10 seconds. If required, add a second / third layer. Please cover the Zircon Bonder bottle immediately after being used. Optimal bond strength is achieved by applying anaxblend Opaquer without contamination due to contact and without long delays followed by the application of the facing acrylic.

For care instructions, please ask for our separate information.

anaxdent GmbH  
Olgastrasse 120a  
D - 70180 Stuttgart  
Germany

+49 (0)711 62 00 92 0  
anaxdent.com

anaxdent North America  
100 W. Main Street  
Suite 300  
Ardmore, OK 73401  
USA

US Toll Free 1-855-anaxdent (262-9336)  
anaxdent.com

**anaxdent GmbH** D-Stuttgart +49(0)711 6200920 [anaxdent.com](http://anaxdent.com)  
**anaxdent North America** US-Ardmore +1 855 262 9336