
anaxDENT

MANUAL

anaxGUM

by Christian Coachman

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anaxGUM is a tissue-colored composite to be used for the esthetic reconstruction of gingival parts. It can either be applied by the dental technician at the dental lab or directly in the patient's mouth by the dentist. **anaxGUM** – for the simple way back to a natural appearance of the gums.

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


The natural appearance of the gums (red esthetics) should have the same importance than the creation of a perfect dental restoration (white esthetics). anaxgum was developed to match the different hues of the natural gingiva as accurately as possible. anaxgum is a well coordinated system of pastes (Gingiva Paste) and liquid modelling materials (Paint/Flow). Its manifold variation possibilities include a wide spectrum for the individual reproduction of the patient's gingiva. It can be used for screw retained implant substructures, in crown and bridge technique, for removable restorations as well as for the individualization of full dentures.

anaxgum was especially designed for the direct use in the mouth. Restorations from the replacement of missing papillas and closing black triangles up to the complete imitation of a gingivoplasty can be realized directly in the patient's mouth. The high density of the composite guarantees its easy polishability and, thus, ensures long-term hygiene of the restoration.

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Picture 01 Preoperative situation with extensive circumferential bone loss. **Picture 02** The diagnostic wax-up shows the amount of soft tissue missing and the ideal tooth shape. It was quite challenging to deal with the space because of the mesiodistal distance which was smaller for the left central incisor than for the right central incisor. This problem was solved by planning the future crown in the buccal position as it was in the preoperative situation. **Picture 03** Immediate implant placement, palatally positioned to realize a screw-retained restoration. This is mandatory for a prosthetic restoration of this type. A vertical augmentation was not attempted. **Picture 04** Customized tissue-colored ceramic shade tabs were produced for the shade selection of the gingiva color. **Picture 05** Metal-ceramic restoration after glaze bake and polishing. **Picture 06** 1st phase of the two-stage gingiva adaptation technique: At the first insertion appointment, the patient's soft tissue is not in an ideal condition due to the previous treatment such as numbing and reshaping. This complicated the adaptation of the shape, color, and texture extremely. In this case, the best solution was to insert the restoration with the finished ceramic part and to add the pink composite resin in a second appointment after the soft tissue has healed. **Pictures 07-12** 2nd phase: The restoration is prepared for the application of anaxgum: mechanical retention, sandblasting, acid etching, steam cleaning, silane, adhesive.

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

Picture 13 First layer: anaxgum Paint. **Picture 14** Light polymerization **Picture 15** The prepared restoration in situ - ready for direct build-up with anaxgum. **Picture 16** Schematic illustration showing the amount of pink ceramic to be placed and the space to be used for the pink composite resin (yellow line). **Pictures 17-18** anaxgum in different colors and translucencies. **Picture 19** The composite resin is applied by means of a spatula, starting with a darker color as background. **Picture 20** Adding a lighter color of the composite resin. **Picture 21** The margins are blended with a flat brush (anaxbrush No.2).

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

Picture 22 Creating grooves as well as the illusion of a gingival sulcus by means of a fine-tipped probe. **Picture 23** Light-cure after each layer. **Pictures 24-26** Customizing the color with light-curing stains (anaxgum Paint). **Picture 27** A brush is used to create the superficial texture. **Picture 28** The finished composite resin build-up. **Picture 29** The restoration is removed to finish the pink composite resin outside the mouth. Due to the interproximal extension of the artificial gingiva, the path of insertion should have two directions: first horizontal and then vertical. To accomplish this, the implant should have either an external connection or a very short internal connection. **Picture 30** The thin edges of the tissue-colored composite are to be trimmed, and the undercuts should be removed to improve durability and hygiene.

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





Picture 31 Ideal design of the pink composite resin after the application in the mouth. **Picture 32** After removing the restoration from the mouth, the thin extension on the edge of the pink composite resin must be trimmed. **Picture 33** Removal of the edge with a diamond bur in an angle of 45 degrees. **Picture 34** The edge is reduced but not removed completely. The undercut is eliminated by adding an additional layer of composite resin providing extra pressure on the soft tissue and smoothing the edge. **Picture 35** The undercuts and thin edges are just gone. **Picture 36** At this stage, the goal is to blend the transition between the artificial and the natural soft tissue. Since the removal of the pink composite resin edge creates a visible and unattractive interface, this must be modified. **Picture 37** The visible interface occurs mainly because of the different light reflection of the artificial part and of the natural gingiva. **Picture 38** A rounded diamond bur is used to adapt the profile of the artificial gingiva at the natural gingiva. **Picture 39** The unnatural angle between the artificial and the natural gingiva is removed. **Picture 40** The light reflection has now a similar direction on the natural and on the artificial soft tissue and, thus, provides the illusion of continuity. This mimetic effect improves even more in the salivating oral environment.

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Picture 41 Undercuts and concavities must be filled with pink composite to ensure a good oral hygiene in this place. Increased bonding is achieved by applying bond LC on the according areas. **Picture 42** Final emergence profile after reshaping and adding composite resin at chairside. **Picture 43** During the final polishing with Pasta Grigia II, care must be taken not to remove the gingiva-like texture. **Pictures 44-45** Hygiene procedures. It is necessary to guide the floss around the artificial gingiva on both sides and, thus, completely clean the interface between the natural and the artificial gingiva.

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
Picture 01 Application of the opaquer anaxgum PO I and anaxblend. **Picture 02** Apply a sufficient amount of dark pink as basic material. Interleaved supplementation with dark pink up to the reflections until the final shape is reached. **Picture 03** Add anaxgum orange to prepare the alveolar tubercles. **Picture 04** Refine the papilla margin with light pink and wipe it down towards the apex with a brush. **Picture 05** Build the alveolar tubercle using anaxblend Flow Dentin A2. The warmer hue of canines can be matched by building up the course of the root with anaxblend Flow B3. **Picture 06** The finished alveolar tubercle: Use anaxblend Flow Dentin A2 to create the interproximal connection of the ridge. **Picture 07** Apply anaxblend Flow S 57 below the sulcus line to slightly brighten up the gingival surface. **Picture 08** Finish to work out the papillae and the stippling. **Picture 09** Polished by means of a handpiece with Pasta Grigia II and a pointed polishing brush. **Picture 10** Buffing provides high gloss.

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With the help of the temperature control you improve the ability of modellation and the ease of flow.

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Fields of application:

- esthetic rehabilitation of gingival parts in case of dental restorations with severe bone loss
- facing of the gingival area of removable combined restorations (telescopic/conical crowns and attachmentretained dentures)
- superstructures for implant cases
- characterization of dentures

Processing time: 1-3 min depending on the lighting conditions

Polymerization: For the polymerization, a light curing unit with an emission spectrum between 350 and 500 nm is required. Please take into consideration that only properly functioning lamps provide the physical features required. Therefore, the light intensity has to be checked regularly according to the manufacturer's instructions.

Polymerization time per layer (max. 2 mm): Heraflash and UniXS: 3 min / Spectramat: 5 min / Spektra LED: 3 min / SolidiliteV: 5 min / GC Labo light LV-III: 5 min / Flashlight: 40 sec