

VK200

Steam generator



Item code: IV800105

Voltage	230	V
Frequency	50	Hz
Power input	1950	W
Maximum operating temperature	135	°C
Tank capacity	2.0	dm ³
Boiler capacity	2.2	dm ³
Water consumption	0.2	l/min
Maximum operating pressure	0.4 (4)	MPa (bar)
Safety valve set pressure	0.55 (5.5)	MPa (bar)
Noise level [EN-3746-LpA]	81.9	dB(A)
IP protection	X3	IP
Height	410	mm
Length	340	mm
Depth	330	mm
Weight	11	kg
Weight inclusive of packaging	13.5	kg

VK200

Steam generator

The latest of EFFEGI BREGA products. In line with the range, the appliance is strong, practical and long lasting. Connection to the water mains is not required. Quick and easy to install: just pour decalcified water into the tank and plug it in; the appliance will be ready in approximately 10 minutes, then all you will have to do is to top up the water tank whenever indicated by the yellow pilot light.

DECALCIFIED WATER TANK

- Manually-filled water tank, whereby the device can be moved around.
- Easy filling: 65 mm Ø inlet hole.
- The appliance can be topped up without switching it off.

GUN

SUPPLY: Steam - Mixture - Hot water

- Practical and thermally insulated.
- Ergonomic, easy to use.
- 3 control buttons.
- Soft touch buttons.



GUN HOSE

- Practical and strong.
- No remaining pressure.
- Extremely flexible.

3-BUTTON GUN

- Blue button: hot water.
- Red and blue button: mixture.
- Red button: steam.

PILOT LIGHTS

- The left pilot light indicates that the heating element is switched on.
- The right pilot light indicates that the water level in the boiler is low and the motor pump is running: boiler is being filled.
- The light goes on to warn that water in the tank is low. The appliance remains switched on and pressurized. Simply top up to continue with the work.

COMPONENTS

- AISI 316 L stainless steel boiler.
- Boiler insulation resistant to 500 °C.
- Varnished stainless steel casing.
- Varnished aluminium lid.
- Copper tubing.
- Brass fittings.
- Safety valve set at 5.5 bar.
- Stainless steel heating element.