**GEBRAUCHSANWEISUNG
NOVOplast M**

DE

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**INSTRUCTION
NOVOplast M**

EN

After receipt of the product, please ensure that you file the lot number and product name in your administrative system lay down. For each work that you produce, please make a note of all the materials used along with all the used lot numbers. Please follow the important guideline in the following instructions for your own protection as well as of your patients and the environment. Used by qualified personnel only.

Indication

Pattern resin for use as a modelling acrylic for multiple indications in dental engineering, especially modul casting. Novoplast m burns out without residue. For modelling with the crown and bridge casting technique, for fixation in the solder technique, for the manufacture of lingual and palatal clamps, links, metal plates, etc., for conus and telescope crowns. Short tests in the mouth possible (polymerisatet)

Material

Self cure pattern resin; Shades: Red, Clear

Properties

Very easy processing and handling, colour – true, best mechanical properties.

Processing**Preparation:**

Insulate the surface of the plaster with standart alginates. Vaseline it's suitable for insulating polyurethane modelling materials, or you can also use our isolation system from Duropont Flash (Isolation). As a rule, no further insulation is required on polished metal surface. If necessary, apply a water thin layer of vaseline or Duropont Flash isolation.

Application:

Can be applied directly to the working model. Work can be tried directly in the mouth. This material is only worked according to the brush methode. Place the NOVOPLAST powder in the one of the shallows of the mixing pad. Into the next recess, put few drops of NOVOPLAST liquid. Dip a clean brush first into the liquid, then into the powder. A pearl of resin has now formed at the end of the brush, saturated with monomer liquid. This mixture can then be applied onto the plaster model which has previously been coated with an isolating agent. The resin hardens after approx. 6 minutes and can, if necessary, be immediately corrected with a hard metal mill. However, when the resin has sufficiently hardened, special attention must be given to the edges and the inner surfaces. To resume working, the model crown must be placed on the stump again. To achieve an even thickness, sand off the work with a fine crosscard mill or a diamond grinder of a minimal thickness of 0,3 – 0,4 mm. It is recommended to control the even thickness of the model crown. Embedding, burning out and casting should be done by known methods. In small single crowns, burning out can be done as in wax models. Bigger crowns should be kept in an oven for one hour at a temperatur of 250°C/482°F. Only after then can the temperature be increased.

Important information:

The liquid contains methyl methacrylate (mma) and is irritant as well as easily flammable (Flash point +10°C, 48°F). The liquid is harmful to the eyes, the respiratory system and the skin. Possible skin reactions by direct contact. Keep in a safe, well ventilated place. Do not inhale vapors. Keep away from fire – do not smoke. This product must be disposed of safely, not through the drains. Protect from electrostatic sources. Wear a mask when grinding and use air suction. Always use spatula when handling. If a patient is known to be allergic to any of the ingredients in Novoplast M, the material must not be used.

Durability and storage

Store the material in a cool, dark, well-ventilated place. Storage temperature: 2–30°C / 36–86 °F. Do not use the material after the expiry date. Keep out of the reach of children.

Quantities

Powder: 80g, 100g, 200g, 500g, 900g
Liquid: 50ml, 100ml, 250ml, 500ml 1.000ml