I-BOND 02

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GB

Nickel Chrom based metal alloy for porcelain to metal restorations, contains no beryllium and precious alloys.

I-Bond 02 is a biocompatible nonprecious nickel chrome based alloy. **I-Bond 02** contains no beryllium and fulfills the recommendations of the standards EN ISO 22674 and EN ISO 9693-1. The low Vicker´s hardness of **180** (HV 10) leads to good and easy milling and grinding behaviour and gives perfect polished surfaces.

COMPOSITION (mass %)	TECHNICAL DATA (Orientation values)	
Ni 64.3	Type	3
Cr 24.2	Density	8,2 g/cm ³
Mo 10	Vickers hardness HV 10	180 (HV 10)
Si 1	Coefficient of thermal expansion	25 - 500°C 13.9 x 10-6K-1 20 - 600°C 14.0x 10-6K-1
Co, Fe, Ce, Nb <1	Melting interval	1250-1340°C
	0,2%- Elongation limit	330 MPa (N/mm ²)
	E-modul	ca. 205000 MPa (N/mm ²)
	Ductile yield A5	15%
	Casting temperature	ca. 1440°C

Recommendations for Use

Wax-up:

Do the wax-up as usual. The wall thickness in wax must not be thinner than **0.35 mm**. Lead wax sprues indirectly. For wax sprues use round wax wires with **Ø 2.5 mm** for single crowns and **Ø 2.5-3.0 mm** for bridges. For greater prosthetic constructions with more than 4 pontics use a distribution bar with **Ø 4.0- 4.5 mm**.

Investing and Casting:

For investing use phosphate bonded investments for crown and bridge work **Intervest K+B Speed** or **Interfine K+B Speed**. Preheat the investment to about **850 to 900 °C**. Hold the end temperature of investment for a minimum **30 minutes**. Refer to manufacturer´s instructions for use for the casting machines. For **I-Bond 02** use an individual ceramic crucible to prevent contamination with other alloys. Clean crucible after every use.

When melting by induction heating start casting as soon as the ingots have collapsed and oxide net cracks. For melting by flame heat rotate the reductive zone of the flame around ingots. Start casting as soon as the bath

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begins to vibrate. Allow the cylinder to cool down slowly to the room temperature and then deflask.

Ceramic:

Use ceramic, suitable for porcelain to metal restoration.

- Sand-blast the surface by use of a pencil-blaster. Blast all of the surface on which the porcelain is to be applied with **Al₂O₃ - Interlox 250 µm** and clean the surface (ultrasonically cleaning, under vapour or boiling water, diethyl acetate).
- If oxide firing is carried out to check the surface, it need to be fired on **960°C** at vacuum for **5-10 minutes**. After cooling the frame has to show slightly dark greenish colour. Oxide layer should be removed with **Al₂O₃ - Interlox 250 µm**. The frame is cleaned under vapour or boiling water
- The opaque is then applied on the surface by a first thin wash firing followed by a second covering opaque.
- Firing has to be done according to the ceramic manufacturer´s instructions.
- After each firing step (dentine firing, build-up and glazing) cool down in a normal manner.**
- Remove ceramic mechanically. Hydrofluoric acid (HF) corrodes the metal.

Finishing:

After firing of the metal-ceramic prosthetic work sand-blast the visible surface of the frame with **Al₂O₃ - Interlox 50 µm** and then polish it with suitable grinding and polishing instruments for nickel alloys. At the end use **Universal polishing paste for Co-Cr-Mo alloys REF 460** and polish up to high gloss. Clean polished surface with vapour or boiling water.

Soldering and Welding:

Soldering before firing of the frame can be carried out with suitable base-metal-solder and high temperature flux. For welding with laser use suitable base-metal welding wires.

Final Safety Notes!

Metal dusts are harmful. Use dust extractor and wear gloves. Consider allergic hypersensitivities to contents of the alloy.

Reciprocal Actions!

In case of occlusal or approximal contact of different alloys electrochemically based reactions may very rarely occur.

Warranty!

Whether given verbally, in writing or by practical instructions, our recommendations for use are based upon our own experience and trials and can only be considered as standard values. Our product are subject to further development. Therefore alterations in construction and composition are reserved.

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