# Bellavest<sup>®</sup> SH

Shock-heat or conventionally heatable, phosphate-bonded precision casting investment material for all crowns, bridges and pressable/press-to-metal ceramics Bellavest SH, phosphate-bonded dental casting investment material: Type 1 (for the production of inlays, crowns, bridges and other find restorations). Class 2 (recommended for burging aut during sheek besting)

**7** BEGO

en

bridges and other fixed restorations), Class 2 (recommended for burning out during shock heating) **Safety instructions** Please read and follow the instructions in the insert "Safety instructions and general instructions for BEGO investment materials"! This material contains quartz which causes lung damage when breathed in during prolonged or repeated exposure. We recommend sufficient ventilation or wearing a PF2 protective mask as suitable protection measures. Preparation Pressable/Press-to-metal ceramics Crown and bridge alloys • Fix and sprue the wax-up according to • Wax the sprued copings on the BEGO *base socket mould former* the instructions of the ceramics so that the distance to the mould edge and top surface is at manufacturer. least 5 mm (1/a''). Spray a thin coat of *Aurofilm* wetting agent and blow completely dry. Plastic copings (e.g. Pattern Resin or Palavit G) must be thinly coated with wax. • Use BEGO *fleecy inlay strips*: • **1 strip** for metal mould rings in sizes 1+3, 2 strips on top of each other for sizes 6 + 9 as well as for all non-precious alloys.

#### Investment



- 2 strips on top of each other for sizes 1+3,
  2 strips on top of each other for sizes 6 + 9 as well as for all non-precious alloys. Handling: The strips must be approx<sup>1</sup>/<sub>2</sub> cm longer than the circumference of the mould ring. Moisten strips slightly. Press strips in mould ring so that they overlap and are flush with the top edge of the mould ring. Slip over the wax-up and press the lower edge of the mould ring into the base socket mould former.
  Liquid: BegoSol<sup>®</sup> HE (Frost-sensitive! Storage and transport temperature: +5 °C to +35 °C/10 °F to 95 °F)
  Before mixing, rinse out the clean mixing bowl with water and wipe off. Mixing bowls that are not clean or are dry withdraw moisture from the investment material!
  Mix the liquid and powder using a spatula for 30 seconds (or mix mechanically using a spatula). Then mix for 60 seconds in a mixing unit under a vacuum at 350 rpm. Keep under vacuum for additional 30 seconds without stirring.(Mixing without mixing unit: 2 minutes on the vibrator.)
  Available working time: approx. 5 minutes (20 °C/70 °F, 50 % liquid).
- At higher room temperatures the working time will be reduced!
  Fill crowns carefully with a fine instrument. Fill the mould ring on the vibrator at the lowest vibration level.
- Do not vibrate any more after filling!
- If heating is to be carried out without a ring, remove the ring used for investment as soon as possible after complete setting of the investment material (at 20 °C/70 °F after approx. 15 minutes); metal mould rings cannot be removed.
   Pressing moulds must be left to cool for 25 30 mins after their being filled until the investment material has completely set. Any deformation, caused, for example, by moving the mould or premature removal from the mould, can cause micro cracks in the investment material, which can cause the mould to break during pressing.
- For shock heating, comply with the time window foreseen for insertion (20-30 minutes after mixing is initiated) and the insertion temperature (900 °C)!

Mixing ratio		100 g Bellavest® SH : 25 ml liquid			
Mould size	90 g bags	100 g bags	160 g bags		
1	1/22,5 ml	1/25 ml	_		
3	2/45 ml	2/50 ml	1/40 ml		
6	4/90 ml	4/100 ml	2/80 ml		
9	6/135 ml	5/125 ml	3/120 ml		



for pressable ceramics
 (lowering and staining technique)

(layering and staining techniques)		for all in
Inlay MO and OD	60-70%	
Inlay MOD	70-80%	The concent according to
Crowns, veneers and bridges	75-85%	dilute to les


- Press-to-metal ceramic (Metal): 70 % for all indications
- Press-on ceramics: (zirconium dioxide) for all indications: max. 40 %

The concentrations are standard values and can be adapted according to the working conditions and object size. Never dilute to less than 20 % !

%	20%	30 %	40 %	50 %	60%	70%	80 %	90 %
HE/H <sub>2</sub> 0	8/32 ml	12/28 ml	16/24 ml	20/20 ml	24/16 ml	28/12 ml	32/8 ml	36/4 ml
for aroun and bridge alloye								

• 1	for	crown	and	bridge	alloys
-----	-----	-------	-----	--------	--------

Modelling:	made of wax without pressure	made of wax with pressure (4 bar)	made of plastic without pressure (e.g. Pattern Resin)	made of plastic with pressure (4 bar)
Inlays and partial crowns	35 %	40 %	_	_
<b>Crowns, bridges and primary parts</b> in precious metal in precious metal-to-ceramic alloys	45 % 50 %	50 % 60 %		
Secondary parts in precious metal cone, ring telescope, full telescope, groove-shoulder attachment	_	_	45-75%	50-80%
Crowns and bridges in non-precious (Ni-Cr) metal-to-ceramic alloys	75 – 85 % 70 – 75 %	80 – 90 % 75 – 80 %		
Non-precious double crowns (external parts)	_	_	90-100%	_



# Bellavest<sup>®</sup> SH

Shock-heat or conventionally heatable, phosphate-bonded precision casting investment material for all crowns, bridges and pressable/press-to-metal ceramics Bellavest SH, phosphate-bonded dental casting investment material: Type 1 (for the production of inlays, crowns, bridges and other fixed restorations), Class 2 (recommended for burning out during shock heating)

#### Preheating



	Shock heating	Conventional heating	
Setting time after investment	20–30 minutes	at least 30 minutes	
Insertion temperature	900 °C / 1,650 °F	Room temperature (or 250 °C / 500 °F) *	
Holding level	_	250 °C / 500 °F (with 5 °C (9 °F) / min )**0	
Final temperature		250 °C / 500 °F (with 7 °C (12 °F)/min)**	
Precious metal Precious metal-to-ceramic alloys Non-precious metal Pressable ceramics	700 °C / 1,290 °F 850 °C / 1,560 °F 900 °C / 1,650 °F up to 900 °C / 1,6	700 °C / 1,290 °F 850 °C / 1,560 °F 900 °C / 1,650 °F 550 °F (Follow manufacturer's instructions!)	
Hold times for holding level and final temperature	30-60 minutes (depending on size and number of moulds)		

\*/\*\* Only for furnaces with conventional control / with computer control.

**Press-to-metal ceramics:** for pressing to non-precious metal alloys, preheat Bellavest<sup>®</sup> SH moulds **rapidly only**; for precious-metal alloys, preheat **normally or rapidly**!

### Shock heating

**Only for mould sizes 1 to 6** • Roughen mould bottom slightly after setting • Place moulds upright in the furnace (funnel former facing down) and without direct contact to the heating source (use spacer or ceramic plate) • **Always comply with setting time and insertion temperature**!



**Risk of injury in connection with shock heating!** Place all moulds in the furnace within 10 seconds and then keep the furnace door closed for 15 minutes!

Inserting further moulds leads to temperature decrease and thus to considerable extension of the preheating process. After casting/pressing allow the moulds to cool down until warm to the touch, in a protected and designated location,

do not quench in water! Investment materials contain quartz. Do not inhale dust! Danger of lung damage (silicosis, lung cancer). To avoid dust during deflasking, place the moulds in water, once they have cooled down completely

#### After casting



#### Data



### Availability and recommendations




Whether given verbally, in writing or by practical instructions, our recommendations for use are based upon our own experience and trials and can only be considered as standard values. Our products are subject to constant further development. Therefore alterations in construction and composition are reserved.

Info: Phone +49 421 2028-380 www.bego.com

	BegoS	ol® HE	
	50%	80%	
Available working time at 20 °C / 70 °F	approx. 5 min	approx. 4.5 min	
Total expansion	1.7%	2.2%	
Minimum shelf life	2 years		
Characteristic mater accordance with DIM			
Beginning of setting (Vicat time)	approx. 10.5 min	approx. 10 min	
Compressive strength (after 2 hours)	4.2 MPa	5.1 MPa	
Linear thermal expansion	0.8%	0.9%	

after casting, until they are thoroughly moistened.

Expansion curve Bellavest® SH



This product was made according to the specifications of DIN EN ISO 15912 and meets its requirements in all respects.

		Carton		Carton		
	90 g bag			12.96 kg (144 bags)	- 54257	
	00 g bag 60 g bag	5.0 kg (50 bags) 4.8 kg (30 bags)	- 70060 - 54247	12.8 kg ( 80 bags)	- 54252	
BegoSol® HE		1000 ml (1 bottle)	- 51095	5000 ml (1 canister)	- 51096	
Aurofilm(100 ml)BEGO base socket formerSize 3(4 pieces)Size 6(4 pieces)Size 9(4 pieces)	52019 52627 52628 52629	BEGO <i>metal mould ri</i> . Size 3 (4 pieces) Size 6 (4 pieces) Size 9 (4 pieces)	ngs 52422 52423 52424	45 mm (3 x 30 m)	52409 52408	
Manufacturer	REF	Article number	Use b		Warning	
Date of manufacture	LOT	Batch number		ve the ctions for use		

