# Variolink® Esthetic DC

## EN Instructions for Use

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## English

### Description

Variolink® Esthetic is a colour-stable, adhesive luting system for the permanent cementation of ceramic and composite resin restorations. Variolink Esthetic is offered in a purely light-curing (Variolink Esthetic LC) and a dual-curing version (Variolink Esthetic DC).

The special filler composition gives Variolink Esthetic a very high radiopacity.

### Shades

Variolink Esthetic is available in the following five shade gradations:

Shade		Effect
Light +	Translucency approx. 5% (opaque)	Lighter / more whitish
Light	Translucency approx. 10%	
Neutral	Translucency approx. 17%	None
Warm	Translucency approx. 12%	
Warm +	Translucency approx. 8.5% (opaque)	Darker / more yellowish

The gradation of Variolink Esthetic shades is based on the effect that a luting composite exerts on the brightness value of the final restoration. Variolink Esthetic Neutral does not affect the brightness value. At the same time, it demonstrates the highest translucency and is thus neutral in shade. Light and Light + make the restorations lighter, while Warm and Warm + create a darker overall shade impression. In order to evaluate the overall effect of the restoration in conjunction with the various Variolink Esthetic shades prior to permanent cementation, we recommend using the Variolink Esthetic try-in pastes.

### Working time

The working and curing times depend on the ambient temperature.

The following times apply as soon as Variolink Esthetic DC has been extruded from the automix syringe:

When used purely self-curing	Room temperature 23°C ± 1°C/ 73°F ± 2°F	Intraoral 37°C ± 1°C/ 99°F ± 2°F
Working time	approx. 5 min	approx. 2 min
Setting time (incl. working time)	approx. 8 min	approx. 4 min

### Mixing ratio

Variolink Esthetic DC is always dispensed from the automix syringe in the optimum ratio.

## Composition

The monomer matrix of Variolink Esthetic is composed of urethane dimethacrylate and further methacrylate monomers. The inorganic fillers are ytterbium trifluoride and spheroid mixed oxide. Initiators, stabilizers and pigments are additional ingredients.

The particle size is 0.04–0.2 µm. The mean particle size is 0.1 µm.

The total volume of inorganic fillers is approx. 38%.

## Indication

- Adhesive luting of glass-ceramic, lithium disilicate glass-ceramic and composite restorations (inlays, onlays, partial crowns, crowns, bridges).
- Restorations made of opaque ceramics, e.g. oxide ceramics, can only be permanently cemented if an adhesive is additionally used that is separately light-cured.
- The following overview provides recommendations for the selection of the suitable Variolink Esthetic shade:

Shade	Effect	Indications	
		Inlays / Onlays	Crowns / Bridges
Light +	considerably lightens	--*	✓
Light	slightly lightens	✓	✓
Neutral	no shade effect	✓	✓
Warm	slightly darkens/ more yellowish	✓	✓
Warm +	considerably darkens/ more yellowish	--*	✓

\* The use of the shades Light+ and Warm+ may lead to visible restoration margins when cementing inlays/onlays.

## Contraindication

The application of Variolink Esthetic is contraindicated

- if a dry working field cannot be established or the stipulated working procedures cannot be applied;
- if a patient is known to be allergic to any of the ingredients of Variolink Esthetic.

## Side effects

Systemic side effects are not known to date. In individual cases, allergic reactions to individual components have been reported.

## Interactions

Phenolic substances (e.g. eugenol, wintergreen oil) inhibit polymerization. Consequently, the application of products containing these components, e.g. mouth rinses and temporary cements, must be avoided.

Disinfectants with an oxidative effect (e.g. hydrogen peroxide) may interact with the initiator system, which in turn may impair the curing process.

Therefore, do not disinfect the preparation and the syringe using oxidative agents. Disinfection may be carried out, for instance, by wiping with medical alcohol.

Alkaline jet mediums applied on dentin (e.g. Airflow) may compromise the effect of self-etching adhesives.

Hemostatic substances may inhibit polymerization and/or lead to discoloration. Therefore, the directions for use of these substances must be observed.

## Application

For more detailed information also refer to the Instructions for Use of the products used in conjunction with Variolink Esthetic.

### 1 Removal of temporary restoration and cleaning of the cavity

Remove possible residue of the temporary luting composite from the cavity or the preparation with a polishing brush and an oil- and fluoride-free cleaning paste (e.g. Proxyl<sup>®</sup> fluoride-free prophylactic paste). Rinse with water spray.

Subsequently, dry with oil- and moisture-free air. Avoid overdrying.

#### Note:

Cleaning with alcohol can lead to dehydration of dentin.

### 2 Try-in of the restoration and isolation

Seat the restoration using the desired Variolink Esthetic try-in paste and check the shade, fit and occlusion of the restoration.

Care should be taken when checking the occlusion of fragile and brittle ceramic objects before they are permanently cemented, as there is a risk of fracture.

If necessary, make adjustments with fine diamonds at medium speed and with slight pressure and adequate water cooling. Polish ground surfaces.

When an adhesive luting protocol with composites is used, safe isolation of the operating field – preferably with a rubber dam, e.g. OptraDam<sup>®</sup> Plus, or alternatively with cotton rolls and a saliva ejector – is required. If the dental hard tissue is contaminated with saliva or blood during try-in, it needs to be cleaned again as described in section 1.

### 3 Pre-treatment of the restoration

3.1 If a restoration which has been pre-treated in the dental laboratory is contaminated with saliva or blood during try-in, the contaminated bonding surface of the restoration needs to be cleaned as follows:

- Following try-in, thoroughly rinse the restoration with water spray and dry with oil-free air.
- Shake Ivoclean before use and completely coat the bonding surface of the restoration with the product using a microbrush or bush.
- Leave Ivoclean to react for 20 s, then thoroughly rinse with water spray and dry with oil-free air.
- Next, prime the bonding surface of the restoration with a suitable bonding agent (e.g. Monobond<sup>®</sup> Plus).  
Make sure to observe the Instructions for Use of the bonding agent used.

3.2 Restorations that have not been pre-treated by the dental laboratory should be rinsed with water spray and dried after the try-in procedure.

3.2.1 Subsequently, condition the bonding surface of the restoration as follows:

- a) Glass-ceramic restorations (e.g. IPS Empress®)
    - Etch with 5% hydrofluoric acid (e.g. IPS® Ceramic Etching Gel) for 60 s or according to the instructions of the manufacturer of the restorative material.
    - Thoroughly rinse with water spray and dry with oil-free air.
  - b) Lithium disilicate glass-ceramic restorations (e.g. IPS e.max® Press, IPS e.max CAD)
    - Etch with 5% hydrofluoric acid (e.g. IPS Ceramic Etching Gel) for 20 s or according to the instructions of the manufacturer of the restorative material.
    - Thoroughly rinse with water spray and dry with oil-free air.
  - c) Zirconium oxide (e.g. IPS e.max ZirCAD, Zenostar®) or aluminium oxide ceramics
    - Sandblast the internal surface of the restoration (sandblasting parameters according to the instructions of the manufacturer of the restorative material).
    - If necessary, clean the restoration in an ultrasonic unit for about one minute.
    - Thoroughly rinse with water spray and dry with oil-free air.
    - **IMPORTANT!** In order to create a strong bond, do not clean the zirconium oxide surfaces with phosphoric acid.
  - d) Composite or fibre-reinforced composite restorations
    - Blast the internal surface of the restoration (according to the instructions of the manufacturer of the restorative material).
    - If necessary, clean the restoration in an ultrasonic unit for about one minute.
    - Thoroughly rinse with water spray and dry with oil-free air.
- 3.2.2 Apply Monobond Plus with a brush or a microbrush to the pre-treated surfaces, let it react for 60 s and then disperse it with a strong stream of air.
- 3.3. Alternatively, restorations made of glass-ceramics and lithium disilicate glass-ceramics (e.g. IPS Empress, IPS e.max PRESS, IPS e.max CAD) may be conditioned with Monobond Etch & Prime® as follows, independent of their pre-treatment:
- Following try-in, thoroughly rinse the restoration with water spray and dry with oil- and moisture-free air.
  - Apply Monobond Etch & Prime on the bonding surface using a microbrush and agitate it into the surface for 20 seconds. Allow to react for another 40 seconds.
  - Thoroughly rinse off Monobond Etch & Prime with water spray until the green colour has disappeared. If any residue remains in the microporosities, the restoration may be cleaned with water in an ultrasonic unit for up to 5 minutes.
  - Dry the restoration with a strong stream of oil- and moisture-free compressed air for approx. 10 seconds.

#### 4 Pre-treatment of the preparation and application of the adhesive (e.g. Adhese® Universal)

Please observe the Instructions for Use of the adhesive used.

##### 4.1 Conditioning with phosphoric acid gel (optional)

The bond to enamel can be improved by selectively etching the enamel or by applying the "etch & rinse" technique. Unprepared enamel surfaces must be conditioned with phosphoric acid (e.g. Total Etch). Please observe the Instructions for Use for the phosphoric acid gel.

###### a) Selective enamel etching

Apply phosphoric acid gel (e.g. Total Etch) onto the enamel and allow it to react for 15–30 s. Then rinse thoroughly with a vigorous stream of water for at least 5 s and dry with compressed air until the etched enamel surfaces appear chalky white.

###### b) Etch & rinse technique

Apply phosphoric acid gel (e.g. Total Etch) onto the prepared enamel first, and then to the dentin. The etchant should be left to react on the enamel for 15–30 s and on the dentin for 10–15 s. Then rinse thoroughly with a vigorous stream of water for at least 5 s and dry with compressed air until the etched enamel surfaces appear chalky white.

##### 4.2 Application of Adhese Universal

- Starting with the enamel, thoroughly coat the tooth surfaces to be treated with Adhese Universal.
- The adhesive must be scrubbed into the tooth surface for at least 20 s. This time must not be shortened. Applying the adhesive on the tooth surface without scrubbing is inadequate.
- Disperse Adhese Universal with oil- and moisture-free compressed air until a glossy, immobile film layer results.  
**Important note:** Avoid pooling as this may compromise the accuracy of fit of the definitive restoration.
- Light-cure Adhese Universal for 10 s using a light intensity of  $\geq 500 \text{ mW/cm}^2$  (e.g. Bluephase® Style).

#### 5 Seating of the restoration

- For each application place a new mixing tip on the syringe. Dispense Variolink Esthetic DC from the automix syringe and apply the desired quantity directly to the restoration.  
As the luting material will cure in the used mixing tip, it serves as a seal for the remaining contents of the syringe until needed again (replace with a new tip before the next use).
- Variolink Esthetic DC should be processed quickly after extrusion from the automix syringe and the restoration seated into place. Apply Variolink Esthetic DC directly to the internal surface of the restoration.  
**Note:** Direct application of Variolink Esthetic DC to a preparation that has been pre-treated with adhesive may - depending on the adhesive used - lead to an acceleration of the curing process and affect the accuracy of fit of the restoration.
- Seat the restoration and **fix/hold it in place during excess removal.**
- Remove all excess luting material.
  - a) Wiping technique  
Remove excess material immediately with a brush, dental floss or scaler. Make sure to remove excess material in time, especially in areas that are difficult to reach (proximal areas, gingival margins, pontics).
  - b) Quarter technique (max. 6 bridge abutments) - Light-curing of excess and subsequent removal  
Light-cure excess material with the polymerization light (e.g. Bluephase Style) **for 2 s** per quarter surface (mesio-oral, disto-oral, mesio-buccal, disto-buccal) at a distance of max. 10 mm.  
Thereafter, excess cement is easy to remove with a scaler. Make sure to remove excess material in time, especially in areas that are difficult to reach (proximal areas, gingival margins, pontics).

## 6 Polymerization

- As with all composite systems, Variolink Esthetic is subject to oxygen inhibition. This means that the surface layer (approx. 50 µm) does not polymerize during curing, as it comes in contact with atmospheric oxygen. To prevent this, cover the restoration margins with glycerine gel/air block (e.g. Liquid Strip) immediately after excess removal.
- Polymerize Variolink Esthetic in segments, starting with the proximal margins:

Light intensity	Exposure time per mm ceramic and segment	
500 - 1,000 mW/cm <sup>2</sup>	20 s	
≥ 1,000 mW/cm <sup>2</sup>	10 s	e.g. Bluephase Style

- Rinse off glycerine gel/airblock (e.g. Liquid Strip).

## 7 Finishing of the completed restoration

- Check occlusion and functional movements and make adjustments if necessary.
- Finish the cement lines with finishing diamonds if necessary.
- Smooth out the cement lines using finishing and polishing strips and polish them with suitable polishing instruments (e.g. OptraPol®).
- If necessary, finish the restoration with suitable polishers (ceramic: e.g. OptraFine®; composite resin: e.g. OptraPol).

## Warning

- Prevent any contact of uncured Variolink Esthetic with the skin / mucous membrane and eyes.
- Uncured Variolink Esthetic may have a slight irritating effect and cause sensitization to methacrylates.
- Customary medical gloves do not provide protection against the sensitizing effect of methacrylates.

## Shelf life and storage

- Store Variolink Esthetic at 2–28 °C/36–82 °F.
- Do not use Variolink Esthetic after the expiry date.
- Do not disinfect syringes with oxidizing disinfectants.
- In order to ensure a tight seal of Variolink Esthetic DC automix syringe, the mixing tip is left on the syringe after use.
- Expiry date: see note on syringes and packages.

## Keep out of the reach of children!

### For use in dentistry only.

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

**Ivoclar Vivadent Pty. Ltd.**

1 – 5 Overseas Drive, P.O. Box 367, Noble Park, Vic. 3174, Australia

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[www.ivoclarvivadent.us](http://www.ivoclarvivadent.us)

**Ivoclar Vivadent Shanghai**

Trading Co., Ltd., 2/F Building 1, 881 Wuding Road, Jing An District

200040 Shanghai, China

Tel. +86 21 6032 1657, Fax +86 21 6176 0968

[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

**Ivoclar Vivadent Marketing Ltd.**

Calle 134 No. 7-B-83, Of. 520, Bogotá, Colombia

Tel. +57 1 627 3399, Fax +57 1 633 1663

[www.ivoclarvivadent.co](http://www.ivoclarvivadent.co)

**Ivoclar Vivadent SAS**

B.P. 118, F-74410 Saint-Jorioz, France

Tel. +33 4 50 88 64 00, Fax +33 4 50 68 91 52

[www.ivoclarvivadent.fr](http://www.ivoclarvivadent.fr)

**Ivoclar Vivadent GmbH**

Dr. Adolf-Schneider-Str. 2, D-73479 Ellwangen, Jagst, Germany

Tel. +49 7961 889 0, Fax +49 7961 6326

[www.ivoclarvivadent.de](http://www.ivoclarvivadent.de)

**Wieland Dental + Technik**

GmbH & Co. KG, Schwenninger Strasse 13, D-75179 Pforzheim, Germany

Tel. +49 7231 3705 0, Fax +49 7231 3579 59

[www.wieland-dental.com](http://www.wieland-dental.com)

**Ivoclar Vivadent Marketing (India) Pvt. Ltd.**

503/504 Raheja Plaza, 15 B Shah Industrial Estate

Veera Desai Road, A, ndheri (West), Mumbai, 400 053, India

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[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

**Ivoclar Vivadent Ltd.**

12 Omega St, Rosedale, PO Box 303011 North Harbour

Auckland 0751, New Zealand

Tel. +64 9 914 9999, Fax +64 9 914 9990

[www.ivoclarvivadent.co.nz](http://www.ivoclarvivadent.co.nz)

**Ivoclar Vivadent Polska Sp. z o.o.**

Al. Jana Pawla II 78, 00-175 Warszawa, Poland

Tel. +48 22 635 5496, Fax +48 22 635 5469

[www.ivoclarvivadent.pl](http://www.ivoclarvivadent.pl)

**Ivoclar Vivadent Marketing Ltd.**

Prospekt Andropova 18 korp. 6/, office 10-06, 115432 Moscow, Russia

Tel. +7 499 418 0300, Fax +7 499 418 0310

[www.ivoclarvivadent.ru](http://www.ivoclarvivadent.ru)

**Ivoclar Vivadent Marketing Ltd.**

Qlaya Main St., Siricon Building No.14, 2nd Floor, Office No. 204

P.O. Box 300146, Riyadh 11372, Saudi Arabia

Tel. +966 11 293 8345, Fax +966 11 293 8344

[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

**Ivoclar Vivadent Pte. Ltd.**

171 Chin Swee Road, #02-01 San Centre, Singapore 169877

Tel. +65 6535 6775, Fax +65 6535 4991

[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

**Ivoclar Vivadent S.L.U.**

C/ Ribera del Loira nº 46, 5ª planta, 28042 Madrid, Spain

Tel. + 34 913 757 820, Fax + 34 913 757 838

[www.ivoclarvivadent.es](http://www.ivoclarvivadent.es)

**Ivoclar Vivadent AB**

Dalvägen 14, S-169 56 Solna, Sweden

Tel. +46 8 514 939 30, Fax +46 8 514 939 40

[www.ivoclarvivadent.se](http://www.ivoclarvivadent.se)

**Ivoclar Vivadent Liaison Office**

: Tesvikiye Mahallesi, Sakayik Sokak, Nisantas' Plaza No:38/2

Kat:5 Daire:24, 34021 Sisli – Istanbul, Turkey

Tel. +90 212 343 0802, Fax +90 212 343 0842

**[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)**

**Ivoclar Vivadent Limited**

Ground Floor Compass Building, Feldspar Close, Warrens Business Park,  
Enderby, Leicester LE19 4SE, United Kingdom

Tel. +44 116 284 7880, Fax +44 116 284 7881

**[www.ivoclarvivadent.co.uk](http://www.ivoclarvivadent.co.uk)**

**Ivoclar Vivadent, Inc.**

175 Pineview Drive, Amherst, N.Y. 14228, USA

Tel. +1 800 533 6825, Fax +1 716 691 2285

**[www.ivoclarvivadent.us](http://www.ivoclarvivadent.us)**