


Duceram® Kiss

Directions for Use
 Gebrauchsanweisung
 Mode d'emploi
 Istruzioni per l'uso
 Instrucciones de uso

THE DENTAL
 SOLUTIONS
 COMPANY™

 **Dentsply
 Sirona**

Instructions for use	1	en
Verarbeitungshinweise	6	de
Indications de mise en oeuvre	11	fr
Note per la lavorazione	16	it
Consejos para su utilización	21	es
تعليمات بخصوص التحضير	26	ar
Упътване за употреба	31	bg
Pokyny pro zpracov	36	cs
Bearbejdselsindikationer	41	da
Υποδείξεις επεξεργασίας	46	el
Napomene za obradu	51	hr
Feldolgozsi tmutat	56	hu
加工について	61	ja
가공처리 주지사항	66	ko
Apdorojimo nuorodos	71	lt
Nordes sakar ar prstrdi	76	lv
Instructies voor de bewerking	81	nl
Informacje dotyczce obrbki	86	pl
Instrues de utilizao	91	pt
Instruiuni de utilizare	96	ro
Рекомендации по работе с материалом	101	ru
Beredningsanvisningar	106	sv
İřlem talimatları	111	tr
加工说明	116	zh

112	General instructions Duceram Kiss
3	Shade conversion chart/Instructions for use
4	Instructions for use
5	General recommendations for firing

Release date: 2017-09

Duceram® Kiss

Indications for use

Duceram Kiss is indicated for veneering metal framework and copings for the preparation of crowns and bridges.

Product Information

- Duceram Kiss is a high-fusing ceramic material for veneering crowns and bridges with frameworks made of dental alloys with a CTE range of 13.8–15.4 $\mu\text{m}/\text{m} \cdot \text{K}$ (25–600°C).

Contraindications

- Suitable only for the indications listed above
- Duceram Kiss is contraindicated in bruxism or other types of parafunction.
- In addition, Duceram Kiss is contraindicated in situations where the interocclusal dimension is insufficient.

Warnings

Undesired side effects of these medical products are extremely rarely to be expected with proper processing and application. Immunoreactions (e.g. allergies) and/or local unpleasant perceptions (e.g. taste irritations or oral mucosa irritations), however, cannot completely be excluded on principle. If you know any undesirable side effects – also in doubtful cases –, then please communicate them to us.

In case of hypersensitivity of the patients against Duceram Kiss veneering porcelain or one of the ingredients, this medical product must not be used or only used under strict supervision by the treating physician/dentist.

Known cross-reactions or interactions of the medical product with other medical products or materials existing in the mouth must be taken into account by the physician/dentist when using the medical product.

Please pass on all the above-mentioned information to the treating physician/dentist, when processing this medical product for an individually manufactured object. Notify the dentist or physician in charge of all factors described above if you use this medical product for a custom construction.

- Do not inhale abrasive dusts
- Paste fluid are harmful if swallowed

Precautions

When working with these materials, make sure to comply with the Instructions for Use and the pertinent Material Safety Data Sheets (MSDS).

- For professional use only

Adverse reactions

We are not aware of any risks and/or side effects related to Duceram Kiss veneering porcelain.

Technical data

- CTE dentine: 13,0 $\mu\text{m}/\text{m}\cdot\text{K}$ (25–600 °C)
- Dental ceramics, type 1, class 1 according to DIN EN ISO 6872
- Metal-ceramic bond, flexural strength and chemical solubility according to DIN EN ISO 9693 / 6872
- Use only with alloys with a solidus temperature of at least 1030 °C

Alloy selection

- Duceram Kiss is compatible with high-gold and reduced-gold alloys as well as non-precious alloys. Ask the alloy manufacturer about the alloy's composition and coefficient of thermal expansion. If the cooling times (below) are observed, bonding alloys with a CTE of 13,8–15,4 $\mu\text{m}/\text{m} \cdot \text{K}$ (25–600 °C) can be recommended.

Contraindicated	< 13,8
No tempering/long-term cooling	13,8 to 14,5
3 min. tempering/long-term cooling	14,6 to 15,4
Contraindicated	> 15,4

Duceram Kiss does not discolour on alloys containing silver. Nevertheless, it is recommended to clean the ceramic ovens and all carriers regularly.

Market launch: March 2004

Transport and storage conditions

- Protect the liquids from temperature below freezing.
- Store powders and pastes protected from moisture.
 -  Keep dry
 -  Keep away from sunlight

Please pay attention to the following symbols on the product labels:

- REF Catalogue number
- LOT Batch code
-  Use by
-  Pay attention to the instructions for use
-  Do not reuse
-  Date of manufacture

Combinable Liquids

- Bonder/Powder Opaques:
 - Ducera® Liquid B
 - Ducera® Liquid OCL universal
- Paste Opaques:
 - Fluid Paste opaque
- Shoulder porcelain:
 - Ducera® Liquid Quick
- Dentines/Incisals, etc.:
 - Ducera® Liquid SD
 - Ducera® Liquid Form
 - Ducera® Liquid Blend
- Stain/Glaze:
 - Ducera® Liquid Stain improved
- Isolating Fluid:
 - Ducera® Sep Isolating Fluid

Ceramic furnace

For best results, make sure that the required firing temperatures and times are observed strictly. If necessary, adjust the firing parameters of the ceramic furnace as required.

Names/Explanations of important materials

Power Chroma (PC 1–6)

The Power Chromas are highly chromatic, fluorescent and intensive powders for individual shade adjustments. The power chromas can be used to increase the chroma in the cervical, palatal, and occlusal areas. The powders can be used pure or in a 1:1-mixture. By mixing the Power Chromas with Stand by you can create an excellent mamelon dentine. The shade conversion chart is meant as a guidance.

Stand by

Strongly opalescent, almost transparent multi functional powder. It is possible to use Stand by pure or to mix it with all other powders of the Kiss concept. Hence the Stand by powder has a key function.

Opal Effekt Sunrise/Opal Effekt Sunset

Opalescent effect powder for yellow as well orange/reddish incisal areas. Very well suited for increasing the chroma at the 2nd or 3rd dentine firings. Sunset is mainly used for A-shades, whereas Sunrise is used mainly for B-shades. Both powders can be diluted by using Stand by.

Opal Effekt Sky/Opal Effekt Ocean

Opalescent effect powders for discreet and for strong, deep blue incisal areas – can be diluted by using Stand by.

Opal Effekt Fog

Opalescent effect powder for grayish incisal areas – can be diluted by using Stand by.

White Surface

Whitish opalescent effect powder to enhance occlusal cusps in the posterior region as well as for palatal/lingual ridges in the anterior area – can be diluted by using Stand by.

Final Kiss

Low-fusing, transparent correction powder (Final Kiss) – firing temperature 880 °C.

Table 1: Colour conversion chart, Duceram® Kiss

Shade	A1	A2	A3	A3,5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4	
Standard layering																	
Opaque	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	
Dentine	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	
Incisal	1	2	3	3	5	1	1	4	6	1	5	5	6	2	4	4	
Individual layering																	
Opaque	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	
Opaque Orange	For characterizing occlusal, cervical and palatal regions.																
Opaque Bleach	For extremely bright/bleached teeth. Normally used in connection with bleach dentines only.																
Opaque Gum	For gingival aspects.																
Shoulder SM/F SM	1	2	2+3	2+4	3+4	1	1+3	3	3+5	1	1+4	2+4	4	1+4	2+4	3+4	
Dentine	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	
Power Chroma 1	1+2	2	2+5	3+5	4+6	1	1+3	2+3	3+6	1+6	2+6	3+6	5+6	1+6	2+6	3+6	
Power Chroma 2																	
Power Chroma 3																	
Power Chroma 4																	
Power Chroma 5																	
Power Chroma 6																	
Flu Inside 1	x	x	Mix	Mix		x	x	Mix		x	Mix			x	Mix		
Flu Inside 2					x				x			x	x				x
Opal incisal 1	x	x	Mix	Mix		x	x	Mix		Mix	Mix			Mix	Mix		
Opal incisal 2					x				x			x	x				x

The Kiss colour wheel shows you the right powders for the different shades at a glance.

1 2 3 4 Framework preparation

The framework should be designed in a small anatomic shape for supporting the porcelain and to ensure an even layering thickness of the veneering material.

Use only cross-cut tungsten-carbide cutters to finish the framework. The use of chip-removing tools is recommended. Using stones invariably leads to microscopic overlaps, which in turn will lead to bubble formation in the ceramic veneering material. This is particularly true of soft alloys (copper-free and palladium-free biological alloys).

After finishing with tungsten-carbide cutters, unless otherwise recommended by the alloy manufacturer, sandblast the metal frameworks with alumina (100–150 µm) at a pressure of 2 bars (non-precious alloys: up to 4 bars and 250 µm alumina).

To avoid tensions in the Duceram Kiss veneering porcelain, it is necessary to round off any sharp corners and edges of the frame.

Then the frameworks are cleaned by steam or in a clean ultrasonic bath (see fig. 1–4).

Soldering and laser-welding

For soldering and laser-welding metal frameworks, follow the comprehensive instructions and notes of the alloy manufacturer.

5 6 Applying the opaque

The Duceram Kiss powder opaques cover the framework in a layer of uniform thickness in the basic colour selected. Two opaque firings are recommended both for the paste and for the powder opaque. When using alloys with a low melting range (copper-free and palladium-free alloys such as BiOcclus Kiss) the first firing should be the neutral paste instead of the first opaque. When the neutral paste is used, the temperature for the first opaque firing is reduced to 900 °C (if the firing time is increased by 1 min). This serves to avoid distortion of the metal framework. The neutral paste is not suitable for veneering non-precious frameworks.

The Duceram Kiss paste opaque and the neutral paste can be applied easily with a standard paste opaque brush, without diluting.

If the paste changes its consistency over time and becomes less easy to apply, the original consistency may be restored by adding a small quantity of paste opaque fluid.

Alternatively, the framework may be coated with a very thin paste opaque fluid layer.

Caution: Use only very small quantities of paste opaque fluid. Using too much paste opaque fluid may result in cracks and bubbles during opaque firing.

The powder opaque is mixed with the Ducera OCL universal liquid (see fig. 5–6).

For processing non-precious alloys the application of the NE-Bonder is recommended (powder/paste). The NE-Bonder is mixed with Ducera® Liquid B or Ducera® Liquid OCL universal and applied on the framework to cover it with a uniform thin layer. Then the NE-Bonder is fired at 980 °C, which ensures that it is sintered to the framework in a homogeneous fashion and the subsequent vitrification is performed in a uniform manner. This in turn ensures that the non-precious alloy bonds securely to the ceramic veneer. Please also observe the processing parameters for firing, which are specially adjusted to non-precious alloys (see dentsplysirona.com).

7 8 9 Aesthetic Line Basic

In the standard layering technique (opaque/dentine/incisal), highly aesthetic veneers can be created within a minimum of time. Start by building up the dentine core and then use the cut-back technique to obtain room for the incisal layer. Rebuild the cut-back portion of the crown with incisal material. This is followed by the first dentine firing (see table 2, 3). Complete to the original contour using the appropriate dentines and incisors. This is followed by the second dentine firing (see table 2, 3) and the finishing of the restoration. Finally, this is followed by the glaze firing (see table 2, 3), which may be per-

formed with or without glaze or Duceram Kiss stains for special tooth characteristics (see fig. 7–9).

10 11 12 Ceramic shoulder

- When creating a ceramic shoulder, use a carbide bur to reduce the crown margin so it ends 0,5–0,8 mm above the lowest part of the chamfer or ledge. Also, ensure that the crown margin tapers off softly to avoid tension in the ceramic material.
- Sandblast the framework inside (especially the margin) and outside carefully, as described initially, and afterwards clean the framework (steam-cleaner).
- Mark the preparation border with a graphite-free pencil and seal die as usual, e.g. with superglue.
- Apply Ducera-Sep ceramic insulation generously. Expose the insulation to the air and apply a second layer.
- Mix the shoulder margin material (SM 1–5) with the SD Quick modelling liquid. Afterwards complete the margin area by using the shoulder margin porcelain. Apply the margin porcelain in the cervical quarter. Let the material dry, maybe by using a heat source, and lift the coping carefully from the die. Bake the shoulder margin after the recommendation.
- The gap created by the shrinkage of the ceramic can be complemented by a second firing or after the glaze by using the final shoulder margin materials (F-SM 1–5) (see fig. 10–12).

13 14 15 16 17 18 Aesthetic Line Individual

With the individual layering technique, the power chroma and opal effect powders allow you to create natural-looking and highly sophisticated restorations. With the Power Chromas, you can first line the framework with a highly chromatic, fluorescent basic layer for individual tooth shade characteristics. This

is followed by building up dentine layer, as usual. After the cutback, complete the incisal area with the corresponding opal effect enamel and use e.g. Opal Effect Sky or Ocean (in the incisal areas) or OE Sunrise or Sunset (for the cervical and body areas) as supplements (see fig. 13–18).

Table 2: General firing recommendations – Duceram® Kiss

		Pre-heating °C	Drying time min	Heating Rate °C/min	Final temperature °C	Holding time min	Vacuum hPa	Tempering
Oxide firing		For details, please refer to the processing instructions for the respective alloys.						
Biological alloy program	Neutral paste	575	7:00	55	900	3:00	50	–
	Paste opaque	575	7:00	55	900	3:00	50	–
	Powder opaque	575	5:00	55	900	3:00	50	–
Conventional alloys	Paste op. 1+ 2	575	7:00	55	930	2:00	50	–
	Powder op. 1+ 2	575	5:00	55	930	2:00	50	–
No long-term cooling. kühlung Example: Degudent Kiss	Margin 1	575	7:00	55	920	1:00	50	–
	Margin 2	575	7:00	55	920	1:00	50	–
	Dentine 1	575	6:00	55	910	1:00	50	–
	Dentine 2	575	4:00	55	900	1:00	50	–
	Glaze	575	3:00	55	890	1:00	–	–
	Correction	575	4:00	55	880	1:00	50	–
Long-term cooling from a CTE 14,6 µm/m·K	Final Shoulder	450	4:00	55	660	1:00	50	–
	Dentine 1	575	6:00	55	910	1:00	50	3 min/850 °C
	Dentine 2	575	4:00	55	900	1:00	50	3 min/850 °C
	Glaze	575	3:00	55	890	1:00	–	3 min/850 °C

Table 3: Firing recommendations for non-precious alloys – Duceram® Kiss

Firing	Pre-heating °C	Drying time min	Heating rate °C/min	Final temperature °C	Holding time min	Vacuum hPa	Tempering	Long-time cooling
Oxide firing	For details, please refer to the processing instructions for the respective non-precious alloys.							
Bonder (powder/paste)	575	7:00	55	980	2:00	50	–	–
Paste opaque	575	7:00	55	950	2:00	50	–	–
Powder opaque	575	5:00	55	950	2:00	50	–	–
Margin 1+ 2	575	7:00	55	930	1:00	50	–	–
Dentine 1	575	6:00	55	920	1:00	50	3 min/850 °C	up to 600 °C
Dentine 2	575	4:00	55	910	1:00	50	3 min/850 °C	up to 600 °C
Glaze	575	3:00	55	890	1:00	–	3 min/850 °C	up to 600 °C
Correction	575	4:00	55	880	1:00	50	3 min/850 °C	up to 600 °C
(Final Kiss)								
Final Shoulder	450	4:00	55	660	1:00	50	–	–

The tempering phase must be omitted for all non-precious alloys with a CTE value of 14.2 µm/m·K (25–600 °C) or less. Stress relief cooling down to base temperature should be performed independent of the CTE value. Please also follow the alloy manufacturers' recommendations. To compensate for the low thermal conductivity of non-precious alloys, it is recommended to increase the firing temperature or to extend the firing time for all bridges spanning more than five units.

The values indicated here are recommended values and serve only as a clue. Deviations of the burning results are possible. The burning results depend on the respective output per oven and are due to the manufacturer and to age. Therefore the recommended values must be individually adapted at each burning. We recommend a test burning to control the oven. All indications have been carefully elaborated and tested by us, but they are passed on without any guarantee.

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DeguDent GmbH
Dentsply International Inc.