ISO-K

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Separating liquid for all types of acrylic**Introduction**

Dear customer,

Thank you for purchasing ISO-K. CANDULOR ISO-K is an aqueous alginate-based separating liquid used for isolating acrylics versus plaster or plaster versus plaster.

When used correctly, ISO-K will make work much easier and will leave the surface of the polymerised acrylic smooth and plaster-free.

Instructions for Use

Precondition for a good result is complete wax extraction (if necessary use wax remover) and a well watered plaster model/matrix/counterpart. To remove any remaining traces of wax and grease, it is advisable to boil out the model using water only, i.e. without washing-up liquid.

Apply ISO-K to the warm or already cooled model with a soft brush. Attention: the model should not be too hot when applying ISO-K! Apply the separating liquid evenly without forming "puddles". After complete drying of the 1st separating layer, a 2nd layer should always be applied. It is extremely important to ensure that the first applied layer has already dried before applying the second layer. Do not use compressed air to remove excess isolation!

Processing tip:

The plaster model/matrix/counterpart should be watered for at least 15 minutes in warm water (40-45°C).

After pouring ISO-K into a small receptacle, it is applied with a soft-haired brush which has been slightly dampened beforehand with water.

After use, rinse the brush thoroughly under running water.

Note

Any formed sediment does not affect the separating effect. Streaks and sediment can be removed simply by shaking the separating liquid prior to use.

Storage

Permitted ambient temperature 2 - 28°C / 36 - 82°F

Shelf life

3 years, subject to correct storage

"Please read the label!"

Ordering details

559537 – ISO-K blue 500ml

559538 – ISO-K blue 1000ml

559539 – ISO-K blue 5000ml

559540 – ISO-K pink 500ml

559541 – ISO-K pink 1000ml

559542 – ISO-K pink 5000ml

Warning

This material has been developed solely for use in dentistry. Processing is according to the Instructions for Use. The manufacturer cannot be held liable for any damages arising as a result of failure to observe the Instructions for Use or the stated area of application. The same applies in the event that the product is mixed or processed with other manufacturers' products. The user shall be solely responsible for testing the material with respect to its suitability prior to use for any purpose other than those explicitly stated in the Instructions for Use.

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CANDULOR AG