WiroFine



Shock-heat or conventionally heatable, phosphate-bonded precision investment material for partial denture frames.

WiroFine, phosphate-bonded dental casting investment material:

Type 2 (for the production of complete or partial dentures or other removable restorations),

Class 2 (recommended for burning out during conventional or rapid heating up)



Safety instructions

Please read and follow the instructions in the insert

"Safety instructions and general instructions for BEGO investment materials"! This material contains quartz and cristobalite which causes lung damage when breathed in during prolonged or repeated exposure. We recommend sufficient ventilation or wearing a FFP2 protective mask as suitable protection measures.

WiroFine can be heated up rapidly ("shock heat") or conventionally.

The possible insertion temperatures are room temperature approx. 21 °C / 70 °F (conventional) or 700 °C to 1000 °C / 1290 °F to 1830 °F (shock heat).

General instructions



• Liquid for shock and conventional preheating:

BegoSol® K (Frost-sensitive! Storage and transport temperature: +5°C to +35°C/41°F to 95°F).

Liquid only for conventional preheating:

BegoSol® (storage and transport temperature: $-10 \,^{\circ}\text{C}/14 \,^{\circ}\text{F}$ to $+35 \,^{\circ}\text{C}/95 \,^{\circ}\text{F}$).

- Before mixing, rinse out the clean mixing bowl with water and wipe off.
 Mixing bowls that are not clean or are dry withdraw moisture from the investment material!
- Ideal processing at 21 23 °C, higher processing temperatures reduces the working time.
 Processing width 21 °C/70 °F: approx. 3.5 min
 Processing width 23 °C/73 °F: approx. 3.0 min
- First put in liquid and add powder, mix thoroughly with a spatula for at least 15 seconds.
 Then mix in a mixing unit for 60 seconds under vacuum conditions, as far as possible.
 (Processing without mixing unit: mix for 2 min on the vibrator.)
- Suitable for one-piece casting technique.
- Duplication can be carried out in gel moulds or in silicone moulds.
 When working with a pressure compaction unit, silicone moulds and the duplicate model must be made under the same conditions (2 4 bar)! Recommendation: Allow to set under pressure for 10 minutes. Duplicate in gel moulds only without pressure!
- Fill duplication mould on the vibrator and then remove it immediately from the vibrator.

Duplication

Investment

g	Mixin	g ratio 100	100 g Powder : 20 ml Mixing liquid		

for 2 duplicate models	WiroFine	Liquid	Aqua dest.	Mixing liquid Total Concentratio	
Liquid: BegoSol® K*	1 x 400 g	56 ml	24 ml	80 ml	70 %
Liquid: BegoSol® **	1 x 400 g	40 ml	40 ml	80 ml	50 %

^{*} for shock and conventional preheating

Mixing

• In case that a higher expansion is recommended, the model can be produced with a higher liquid-concentration of 80 %. The liquid-concentration for mould/ring remains unchanged at 70 %!

	Gel (Castogel®, Wirodouble	e®, WiroGel® M)	Silicone (Wirosil®, Wirosil® plus)		
Removal from mould	after 45 – 60 mir	nutes	after 30 – 60 minutes		
Surface treatment	urface treatment Durol Durofluid		k	Durol E	
Drying	30 minutes	10 minute	c	45 minutes	

Surface treatment	Durol	Durofluid*	Durol E
Drying	30 minutes	10 minutes	45 minutes
	(250°C/500°F)	(70-100°C/180-210°F)	(150 °C / 300 °F)
Dipping/spraying	briefly 3 times	spray with thin film	long 1 time
	(approx. 2 seconds)	5 minutes	(approx. 4 seconds)
Subsequent drying	5 minutes (250 °C / 500 °F)	(70-100°C/180-210°F)	1 minute (150 °C /300 °F)

^{*} do not use for duplication with gel

- Please note: Only those models fabricated in gel moulds need immersion hardening with Durol or Durol E.
 Models duplicated with silicone only have to be dried for 10 mins. before applying Durofluid.
- Before investing, prepare the wax-up with Wiropaint plus fine investment material or Aurofilm wetting agent (please follow the processing instructions for the products).
- Fill mould ring on the vibrator and then take away from vibrator immediately. Recommendation: Allow mould ring to set under pressure for 10 minutes.
- Remove the mould ring 10 minutes after investment!

Mixing		Mixing ratio 100 g Powder : 20 ml Mixing liquid					
for 1 mould		WiroFine	Liquid	Aqua dest.	Mixing liquid Total Concentration		
Liquid: BegoSol® K*		1 x 400 g	56 ml	24 ml	80 ml	70 %	
Liquid: BegoSol® **		1 x 400 g	40 ml	40 ml	80 ml	50 %	

^{*} for shock and conventional preheating

^{**} only for conventional preheating

^{**} only for conventional preheating

WiroFine



Shock-heat or conventionally heatable, phosphate-bonded precision investment material for partial denture frames.

WiroFine, phosphate-bonded dental casting investment material:

Type 2 (for the production of complete or partial dentures or other removable restorations),

Class 2 (recommended for burning out during conventional or rapid heating up)



Preheating



	Conventional preheating	Shock heat	
Setting after investment	at least 30 min	20 min	
Insertion temperature *	room temperature (21 °C)	700-1000°C / 1290-1830°F	
Holding levels	250 °C / 500 °F (5 °C/min / 9 °F/min) ** 570 °C / 1060 °F (7 °C/min / 12° F/min) **	- -	
Final temperature	950°C-1050°C/1740°F-1920°F		
Holding times	30—60 min (depending on size a	and number of moulds)	

^{**} Shock heat: After insertion you can heat up to the final temperature immediately.





Risk of injury in connection with shock heating. Caution: Danger of darting flame!

Place all moulds in the furnace within $10 \ \text{seconds} - \text{then} \ \text{keep}$ the furnace door closed for $15 \ \text{minutes}$!





After casting, allow the moulds to cool down until warm to the touch in a protected and designated location; do not quench in water!

Investment materials contains quartz. Do not inhale dust! Danger of lungs harms (silicosis, lung cancer).

To avoid dust during deflasking, place the moulds in water after they have cooled down completely after casting until they are thoroughly moistened.

Data



Processing time at 21 °C/70 °F approx. 3.5 min Characteristic material values in accordance with DIN EN ISO 15912


(70 % BegoSol® K)

Beginning of setting (Vicat time) approx. 6.0 min

Compressive strength

(after 2 hours) approx. 11 MPa Linear thermal expansion approx. 0.8 % Flowability approx. 140 mm

This product was manufactured according to the specifications of DIN EN ISO 15912 and meets its requirements.



Availability and recommendations



WiroFine	1 carton $18 \text{ kg} = 45 400 \text{ g bags}$	- REF 54345
	1 carton $6 \text{ kg} = 15 400 \text{ g bags}$	- REF 54344
	1 carton $6 \text{ kg} = 30 \text{ 200 g bags}$	- REF 54348
BegoSol® K	1 bottle $= 1000 \text{ ml}$	- REF 51120
	1 canister = 5000 ml	- REF 51121
BegoSol®	1 bottle $= 1000 \text{ ml}$	- REF 51090
	1 canister = 5000 ml	- REF 51091

	Castogel®	52052 (6 kg)	Wirosil®	52001 (2 kg)	Durol E	52148 (1000 ml)
	Wirodouble®	52050 (6 kg)	Wirosil® plus	54854 (2 kg)	Durol	52111 (1000 ml)
	WiroGel® M	54351 (6 kg)	Wirosil® duplicating	50070 /	Durofluid	52008 (100 ml)
ı			flask system	52072 (small)	Wiropaint plus	51100 (200 ml)
ı				52083 (large)	Aurofilm	52019 (100 ml)

Tel. +49 421 2028-380 www.bego.com

Whether given verbally, in writing or by practical instructions, our recommendations for use are based upon our own experience and trials and can only be considered as standard values.

Our products are subject to a constant further development. Therefore alterations in construction and composition are reserved.



Manufacturer



Article number



Use by



Worning



Date of manufacture



Batch number



Observe the instructions for use

^{**} Heating rate only applies to furnaces with computer control.