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Instructions for Use

1. Measuring of the Tooth Bow

- After determination of the adequate bow size, the corresponding nets can easily be fixed to the lingual surfaces (PET bow) respectively be adapted (stainless steel bow). After that the bow is corrected to the required size.

2. Fastening of the Bow

- The teeth to be splinted are cleaned and etched with a phosphoric acid agent (e.g. Micro-Etch) at the potential contact surfaces of the bow.
- Now apply a thin layer of e.g. thin layer of composite to the complete etched surface.

PET Bows:

- Apply a little of the thin layer of composite mix in the centre of the lingual surface of the first tooth to be splinted at the beginning of the bow.
- Now press the bow to the tooth with the help of an instrument, e.g. a spatula. Make sure that the material is extruded through the net.
- Allow the composite to cure.
- Proceed the same way with each tooth to be splinted.

Stainless Steel Bows:

- Apply some thin layer of composite mix in the centre of the lingual surface of all teeth to be splinted.
- After that the additionally adapted bow or the wings of the provisional Maryland bridge, which has been made in the dental laboratory before (two visits), are completely pressed to the prepared lingual surfaces with an instrument. Make sure that the material is extruded through the perforation.
- In case of a single visit making of a provisional Maryland bridge, e.g. with the help of the existing tooth crown, it is recommended to attach the crown interdentially as "intermediate bridge element" first, followed by attaching the bow as "bridge structure" step by step respectively completely to the lingual surfaces.
- Finally the fixed bow is covered with bonding for optimisation of the surface.

3. Fields of Application

Splinting and stabilisation of anterior teeth in case of:

- re-implantation and transplantation of teeth
- traumatically or periodontally considerably loose teeth
- surgical treatments
- diastema
- orthodontic treatments
- making of fixed immediate or interim bridge provisions in case of single tooth diastema in the anterior or pre-posterior teeth area (e.g. prior to planned implant provisions).

For the two visit treatment, an additional fixing element can be tightly connected with the net bow thus allowing a more stable fixing of individual prosthetic teeth.

Recommendation: intraoral use not longer than 30 days

The following criteria help you choosing the correct bow material:

PET bows (thickness 0.25 mm):

- Allow a physiological self-flexibility of the teeth in case of non-required interdental blocking of the respective anterior tooth area.
- So a possible disuse-atrophy of the alveolar bone as well as of the tooth pulp because of lacking functional stimuli can be avoided.

Stainless Steel bows (thickness 0.20 mm):

- For required immobilisation of the teeth, e.g. for treatment of root or alveolar process fractures.
- As bridge structure for the making of an immediate and/or intermediate bridge provision, as provisional Maryland bridge or in case of single tooth loss.

Warning

Perfect Splint C is a disposable product. It was not designed to be re-used. Sterilizing it could have a negative impact on the product features. Re-use of the article causes a high risk of infection.

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HAGER & WERKEN GmbH & Co. KG