TriVest

Premium Universal Speed Einbettmasse

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Preparing

Pressable

Fix and sprue the wax-up according to the instructions of the ceramics manufacturer and keep a distance of min. 3 mm / 1/8 " to ring base and top.

Crown and bridge alloys

Wax the sprued copings on ring base former and keep a distance of min. 3 mm / 1/8 " to ring base and top.

Spray a thin coat of **Hydent** Surfactant and slightly blow dry.

Mixing/Pouring

- TriVest Liquid press / cast is Frost-sensitive! Temperatures below -5 °C / 23 °F may cause irreversible cristalization!
- Before mixing, rinse out the clean mixing bowl with water and wipe off.
-> Mixing bowls that are not clean or are kept dry, will withdraw moisture from the investment material!
- Mix the liquid and powder using a spatula for 30 seconds (or use automatic pre-mixing). Then mix for 120 seconds at 350 rpm using a vacuum mixer.
- Available working time: approx. 5 - 6 minutes (20 °C / 68 °F, 50 % Liquid concentration).
-> Higher room temperatures will reduce the working time!
- Fill crowns carefully with a fine instrument. Fill the mould ring on the vibrator at the lowest vibration level.
-> Do not vibrate any further after filling!
- If heating is to be carried out without using a metal ring, remove the rubber sleeve used for investing as soon as possible after complete setting of the investment material (at 20 °C / 68 °F after approx. 15 minutes); metal mould rings cannot be removed. Pressable moulds must be left to cool for 25 minutes after being filled, until the investment material has completely set.
Any deformation caused, for example, by moving the mould or early removal from the mould, may cause micro cracks, which may cause the mould to burst during pressing.
- For speed heating, it is mandatory to stay within the time window of inserting to furnace (25 – 30 minutes after mixing process was initiated) and the insertion temperature of (850 °C)!
-> If you missed this time window, traditional heating process is strongly recommended!
- Preheat the ring for at least 60 minutes.

Mixing ratio

Only use distilled water for diluting!

Liquid-concentration TriVest Liquid press

Pressables

Inlays / partial crowns	45 - 50 %
Crowns, veneers and bridges	50 - 60 %

100 g TriVest : 22 ml Liquid

1	1 / 22 ml
3	2 / 44 ml
6	4 / 88 ml
9	5 / 110 ml

Liquid-concentration TriVest Liquid cast

Crown and bridge alloys

precious alloys	%
Inlays / partial crowns	45 - 55%
Crowns, bridges and primary telescopic crowns	50 - 60 %
Secondary telescopic crowns, CS attachments	55 - 65 %

non-precious alloys	%
Crowns, bridges and primary telescopic crowns	70 - 80 %
Secondary telescopic crowns	75 - 85 %

%	45 %	50 %	55%	60 %	65 %	70 %	75 %	80%	85 %
TriVest Lq. / H ₂ O	10 / 12 ml	11 / 11 ml	12 / 10 ml	13 / 9 ml	14 / 8 ml	15 / 7 ml	16 / 6 ml	17 / 5 ml	18 / 5 ml

TriVest

Preparing

Refractory dies / models

Prepare duplicate - silicon mould as usual.

Spray a thin coat of **Hydent** Surfactant and slightly blow dry.

Mixing/Pouring

- TriVest allows various porcelains to be built directly onto the refractory model; has outstanding dimensional stability and provides sufficiently high strength to withstand repeated firings.
- Thanks to its unique granular formula, TriVest provides a smooth model surface.
- TriVest is designed for use with most metal bonding porcelains which have average coefficients of thermal expansion of $10.0 - 13.5 \times 10^{-6}$ at $100-500^\circ\text{C}$ ($212-932^\circ\text{F}$). It will not cause the porcelain to crack or separate from the model during firing.
- After hand-mixing thoroughly, vacuum mix uniformly for 120 seconds.
-> The smallest possible mixing container (5-6 cm internal diameter) should be used. The use of a larger container will not allow for a uniform consistency, resulting in improper physical properties such as thermal and setting expansion.
- After mixing the refractory, immediately pour the mixture into the impression while using a gypsum vibrator (maximum pouring time after mix: about 2 minutes).
-> Some wetting agents may roughen the surface of the refractory model.
- Remove the model from the impression 1.5 hours after pouring.

Degassing

- Gradually heat the model in the ring furnace from room temperature up to 700°C ($1,292^\circ\text{F}$) within 1 hour. Maintain temperature for 10 minutes.
-> Too rapid rise in temperature may fracture the model due to too rapid evaporation of moisture inside the refractory model.
- Next, place the model in the preheated (700°C / $1,292^\circ\text{F}$) porcelain furnace and raise to $1,050^\circ\text{C}$ ($1,922^\circ\text{F}$) at a rate of 55°C (131°F) per minute under atmospheric pressure. Maintain this temperature for 5 minutes. Remove the model and bench-cool to room temperature.
- It is not necessary to bench-cool in the porcelain furnace.
-> Do not degass the model at a temperature higher than $1,050^\circ\text{C}$ ($1,922^\circ\text{F}$). Otherwise it may affect its dimensional precision. Degassing the model can also be performed in a porcelain furnace, but phosphate-bonded refractory materials give off small amounts of ammonia gas, which may affect the porcelain furnace. Therefore a ring furnace is recommended to heat the model up to 700°C ($1,292^\circ\text{F}$).
- Soak the die before each operation.
- Application and firing of the ceramic are made as per manufacturers manual of the ceramic used.
- After the last firing remove TriVest with a precision blaster using $50\text{ }\mu\text{m}$ glassbeads at a pressure of max. 2 bars.

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Concentration du TriVest Liquid cast

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Safety instructions

- Investment materials contain quartz. Do not inhale dust! Risk of damage to lungs (silicosis, lung cancer).
Recommendation: Use gas mask type FFP 2 – EN 149:2001. Open bag along the perforation or cut open with scissors and avoid generation of dust when filling into mixing bowl. Rinse out empty bag with water before crumpling up. Following eye contact: Remove contact lenses, hold eyelids open and rinse thoroughly with clean, running water for at least 15 minutes. Consult an eye specialist.
- Remove dust at workplace only when moist.
- To avoid dust during deflasking, place the mould in water after it has completely cooled down after casting until it is thoroughly moistened.
- Use suction extraction system with fine dust filter when blasting.
- Remove dried-up mixing liquid only when moist (fine dust). Seal bottle securely after use.
- The setting reaction of the investment-material generates heat. To avoid injuries, please use gloves or suitable muffle tongs.
- Combustion gases resulting from the preheating of casting moulds must be discharged via a suitable smoke outlet.
- Place all moulds in the furnace within 10 seconds – then keep the furnace door closed for 15 minutes!
- Never touch hot moulds with your hands. Use suitable mould tongs to avoid burns.
- After casting, allow the moulds to cool down until warm to the touch in a protected location!

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Contains: Quarz (SiO₂), Cristobalit
H372: Causes damage to lungs through prolonged or repeated exposure: Inhalation
P260: Do not breathe dust/fumes/gas/mist/vapours/spray.
P285: In case of inadequate ventilation wear respiratory protection.
P501: Dispose of contents/container in accordance with local/ regional/national/international regulation.

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General instructions

- The ideal processing temperature for TriVest is 20 °C .
To keep it constant at a higher ambient temperature, use a temperature control cabinet if necessary. At a high ambient temperature also place the mixing bowl and mixing bowl in the temperature control cabinet.
- Observe shelf life dates of investment material and mixing liquid! Do not use without checking after the specified date.
- Do not use crystallized mixing liquid any longer. The higher the concentration of mixing liquid, the greater the expansion.
- Data on the concentration of the mixing liquid are based on our own experience and tests and are to be regarded as a reference value only.

Application-related recommendations made by us, whether given verbally, in writing or by way of practical instructions, are based on our own experience and tests and may therefore be regarded only as general guidelines.

Our products are subject to continuous development. We therefore reserve the right to change the design and composition of our products. Work results may be influenced by such parameters as temperature, mixing intensity, investment system, etc.

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