# MIRAFIT CORE®

.....  
Gebrauchsanweisung  
Instruction for Use  
Mode d'emploi  
Istruzioni per l'uso  
Instrucciones de Uso  
.....

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### **Important information**

If you use a temporary resin material, you have to apply a separating medium (Vaseline) onto the core to avoid an adhering of the temporary material.

### **Core build-up**

1. Isolate the tooth with a rubber dam.
2. Apply 37% etching gel onto the exposed dentin and enamel surfaces. Leave for 20 seconds.
3. Rinse the tooth thoroughly with water to remove the etching agent. Dab or blow dry lightly in order to remove excess moisture. Do not dry too long, the surface should remain slightly moist.
4. Apply two consecutive layers of primer/bonding agent (e.g. Mirafit Bond) onto the dentin and enamel surfaces within approx. 10 seconds. Air-dry gently to remove the solvent. Do not dry too long, the surface should shine slightly.
5. Light-cure according to the manufacturer's instructions (approx. 10 sec. for Mirafit Bond).

### **Chemical curing**

6. When using a matrix band: Place a band around the prepared tooth. Apply Mirafit Core material. It has a preparation time of 2 – 2.5 minutes and cures within 3.5 – 4 minutes.
7. When using a prefabricated core build-up cap (e.g. Mirafit BC), fill Mirafit Core material into the cap. Attach the filled cap to the prepared tooth and leave for 4 minutes until the material is cured.
8. Prepare the core as per the usual procedure.

### **Light curing**

For light curing of A2, A3 and Opaque White, cure tooth surfaces for 40 seconds. Curing depth for light curing is 3 – 4 mm. Possibly uncured material cures chemically within 4 minutes.

### **Cementing of the root post using Mirafit Core**

1. Apply 37% etching gel (e.g. Micro-Etch) into the prepared root canal; leave for 20 seconds.
2. Rinse the root canal thoroughly with water and air-dry gently. Absorb excess water with paper points. Do not dry too long.
3. Dispense one drop of bonding (e.g. Mirafit Bond) into the root canal. Leave for 20 seconds.
4. Remove excess with a paper point. Do not dry too long.
5. Apply 2 layers of Primer / bonding agent (e.g. Mirafit Bond) onto the post. Light-cure according to the manufacturer's instructions (approx. 10 seconds for Mirafit Bond).
6. Insert Mirafit Core material into the root canal, place the post, remove excess Mirafit Core material.

### **Chemical curing**

7. Let the Mirafit Core material cure for 4 minutes.

### **Light curing**

Cure occlusal surfaces for 40 seconds. Curing depth for light curing is 3 – 4 mm. Uncured material cures chemically within 4 minutes.

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HÄGER & WERKEN GmbH & Co. KG