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### Allgemeine Grundlagen

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## Application and processing

### General principles

All instruments must be cleaned, disinfected and sterilised prior to every use, and in the case of instruments supplied non-sterile, also before first use. Effective cleaning and disinfection is an essential requirement for effective sterilisation. Special notes on cleaning/sterilisation are given in the instructions for use. In addition, the operating instructions for your practice equipment must be complied with.

As part of your responsibility for the sterility of the instruments for application purposes, please observe as a matter of principle, that only sufficiently equipment and product-specific validated processes are employed for cleaning/disinfection and sterilisation, that the equipment used (disinfector, steriliser) has been maintained and checked regularly, and that the validated parameters are complied with for every cycle.

**Please also observe** all valid legal regulations as well as the regulations on hygiene for the practice or the hospital. This applies in particular to the different requirements regarding the inactivation of prions. For your personal safety, always wear gloves when handling contaminated instruments.

### Cleaning and disinfection

If possible, an automated process (disinfector) should be employed for cleaning and disinfection of the instruments. A manual process – even when using an ultrasonic bath – should only be employed in case of non-availability of an automated process due to the considerably lower efficacy and reproducibility. Pretreatment is to be carried out in both cases.

### Pretreatment

Coarse contamination must be removed from the products directly after use (within a maximum of 2 hours). Never allow residues to dry in! Immediately after use on the patient, place the instruments in the instrument rack filled with a suitable cleaning/disinfectant agent for interim storage and pre-disinfection/cleaning (storage period 2 hours maximum). A clean instrument rack with new foam disc must be used for every patient.

Then clean the instruments from contamination under running water or in a disinfection solution; the disinfectant should be aldehyde-free (otherwise fixation of blood contamination), have tested efficacy (e.g. VAH or FDA approval or CE certificate), be suitable for the disinfection of instruments and be compatible with the instruments (see section "Material resistance").

Only use a clean soft brush or a clean soft cloth, which is only used for this purpose, for the manual removal of contamination, but never metal brushes or steel wool.

**Please note** that the disinfectant used for pretreatment only serves for personal protection and cannot replace the subsequent disinfection stage to be performed – after cleaning.

### Automatic cleaning/disinfection

#### -Thermal disinfection

When selecting a disinfectant, the following should be observed,

- that the disinfectant offers tested efficacy as a matter of principle (e.g. VAH or FDA approval or CE certification according to DIN EN ISO 15883),
- that, if possible, a tested programme for thermal disinfection is used (min. 10 min. at 93°C or A0-value <3000) (chemical disinfection bears the risk of disinfectant residues on the instruments),

- that the programme used is suitable for the instruments and includes sufficient rinsing cycles,
- that only sterile water or low-germ as well as low-endotoxin water is used for post-rinsing (e.g. High Purified Water HPW)
- that the disinfectant is maintained and checked regularly.

When selecting the cleaning agent to be used, the following should be observed,

- that this is suitable in principle for cleaning the instruments,
- that – in as far as thermal disinfection is not employed – an additional suitable disinfectant with proven efficacy is used (e.g. VAH or FDA approval or CE certification) and that this is compatible with the cleaning agent
- that the chemicals used are compatible with the instruments (see section "Material resistance"). The concentrations specified by the manufacturer of the cleaning agent, and, if applicable, of the disinfectant, must be complied with.

### Process

1. Sort the pre-cleaned instruments on a wash tray for root canal instruments. The cleaning of loose instruments is not permissible.
2. Place the wash tray in the disinfectant.
3. Start the programme.
4. At the end of the programme, remove the wash tray from the disinfectant.
5. Check and package the instruments as immediately as possible after removal (see section Control, maintenance and packaging, if applicable, additional post-drying at a clean location). Instruments and products which cannot be cleaned in the wash tray must – if possible – be disassembled. Please also note that the instruments/products may not touch each other.

### Manual cleaning and disinfection

When selecting the cleaning agent and disinfectant to be used, the following should be observed,

- that this is suitable in principle for cleaning or disinfection of the instruments,
- that the cleaning agent – if it can be used – is suitable for ultrasonic cleaning (no foam formation),
- that a disinfectant with proven efficacy is used (e.g. VAH or FDA approval or CE certification) and that this is compatible with the cleaning agent used
- that the chemicals used are compatible with the instruments (see section "Material resistance"). Combined cleaning/disinfection agents should only be used in case of low pre-contamination (no visible soiling) of the instruments. The concentrations and exposure times specified by the manufacturer of the cleaning and disinfectant agents must be complied with. Only use freshly prepared solutions, only sterile water or low-germ as well as low-endotoxin water (e.g. Aqua purificata (PW)) and only filtered air for drying.

### Process

#### 1. Cleaning

- a. Sort the pre-cleaned instruments on the wash tray for root canal instruments. The cleaning of loose instruments is not permissible.
- b. Place the instruments or the wash tray respectively, horizontally into the cleaning bath for the specified exposure period, such, that the instruments are adequately covered (if applicable, with ultrasonic support or careful brushing with a soft brush).
- c. Then remove the instruments from the cleaning bath and rinse these thoroughly for at least 1 minute with water.

### 2. Disinfection

- a. Place the cleaned and checked instruments in the wash tray into the disinfection bath for the specified exposure period, such, that the instruments are adequately covered.
- b. Then remove the instruments from the disinfection bath and rinse these thoroughly for at least 5 minute with water.
- c. Check, dry and package the instruments as immediately as possible after removal (see section Control, maintenance and packaging). Instruments and products which cannot be cleaned in the wash tray must – if possible – be disassembled. Please also note that the instruments/products may not touch each other.

### Checking

Check all instruments after cleaning or cleaning/disinfection. Instruments with the following deficiencies are to be sorted out immediately:

- plastic deformation
- instrument bent
- untwisted threads
- damaged cutting surfaces
- blunt blades
- missing size mark
- corrosion

Information on limitations to the frequency of use is given in the section "Re-use". Instruments which are still contaminated must be cleaned and disinfected again.

### Maintenance

Reassemble the disassembled instruments again. Instrument oils may not be used.

### Packaging

Please package the instruments into the endo sterilisation trays and then into disposable sterilisation pouches (disposable packaging) which meet the following requirements:

- compliance with DIN EN 11607/ANSI AAMI ISO 11607
- suitable for steam sterilisation (temperature resistant up to at least 138°C sufficient vapour permeability)

### Sterilisation

Use only the sterilisation methods listed below; other sterilisation methods are not permitted.

#### Steam sterilisation

- fractional vacuum/pre-vacuum method (at least 3 vacuum cycles) or gravity displacement method (product must be sufficiently dry)
- steam steriliser according to DIN EN 13060 or DIN EN 285
- validated according to DIN EN ISO 17665 (valid installation and operation qualification (IQ and QR) and product-specific performance qualification (PQ))

- maximum sterilisation temperature of 138°C; plus tolerance according to DIN EN ISO 17665
- sterilisation time (exposure time at sterilisation temperature) at least 20 minutes (at 121°C or 5 minutes at 132°C/134°C).

The rapid sterilisation method or the sterilisation method of unpacked instruments is not permitted. In addition, do not use any hot air sterilisation, no radiation sterilisation, no formaldehyde or ethylene oxide sterilisation and no plasma sterilisation.

### Storage

After sterilisation, the instruments must be stored in the sterilisation package and kept dry and dust-free.

### Material resistance

When selecting the cleaning agents and disinfectants agents:

- please ensure that the disinfectant does not contain phenol strong acids (ph<6) or strong bases (ph>8) and
- that no anticorrosion solutions are used on steel instruments, i.e. solutions must not contain aldehyde or di- and triethanolamine.
- Hydrogen peroxide solutions (H2O2) will damage hand instruments as well as plastic stands.
- all chemicals used must be suitable for the respective instrument.

Never clean the instruments and sterilisation trays with metal brushes or steel wool.

All instruments and sterilisation trays may not be subjected to temperatures higher than 138°C!

### Boxes:

Please ensure that the instrument box does not touch the inner walls or the bottom of the autoclave. Temperatures there can be exceeded the pre-set temperature, which can lead to material deformation. Only use the designated storage and cleaning inserts to store the instrument boxes. Do not use the boxes as containers for inserting instruments in cleaning/disinfection solutions. This could make the plastic brittle.

### Re-use

Instruments can be reused several times- with proper care an if they are not damaged and contaminated. Each re-use or application of non-validated methods is the sole responsibility of the user.

**All liability is disclaimed for failure to follow these instructions or the use of non-validated methods for the re-use of instruments.**

**Always ensure that that sterile packaging/ wrapping is undamaged.**

Instruments/product	Material	Special instructions on cleaning/sterilisation	Re-use	Possible damage in care instructions for use are not followed
M+W root canal instruments	Stainless steel and temperature-resistant plastic		Clean and undamaged instruments can be used up to 8 to 10 times.  Re-use of undamaged pluggers is not limited.	Cracks on plastic handle, corrosion on working part and/or shaft

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Müller & Weygandt GmbH • Reichardsweide 40 • D-63654 Büdingen • Tel.: +49 (0)6042 – 88 00 88