

a bit better



denture teeth
since 1947



enta®

The power of innovation & profitability



As one of the few full suppliers of artificial teeth made of acrylic and ceramic, **enta** looks back on more than 60 years of experience in the development, manufacturing and processing of denture teeth.

The objective: To link up outstanding quality and the largest possible shape and color selection with a first rate cost-benefit ratio.

To achieve this, **enta** has always been performing scientific studies in collaboration with leading scientists and universities – for example with Prof. Dr. W. Kalk of Groningen University or Dr. N. Postema of Nijmegen University.

As a result of that, **enta** was the first to find out, among other things, that full ceramic dentures have decisive health benefits for patients. In addition, **enta** worked out the leading processing and arrangement concept named **"Lingualized Occlusion Concept"** or invented the utilization of only one tooth line for all

An overview of the **enta** brand family:

	
enta[®] cryl Denture teeth made of sustainable acrylic	enta[®] ceram Artificial teeth made of ceramic
enta[®] cryl basic	enta[®] ceram classic
enta[®] cryl classic	enta[®] ceram solution
enta[®] cryl optiform	enta[®] ceram optiform
enta[®] cryl polychrome	enta[®] ceram polychrome

Highest quality in acrylic & ceramic

enta fulfills all individual requirements from one single source: With a complete portfolio of colors and shapes as well as a free selection of the prosthetic teeth material (acrylic, ceramic or a combination of both).

An overview of the advantages of the materials:

	enta cryl	enta ceram
Advantage	acrylic dentures	Artificial teeth made of ceramic
Highest aesthetics – especially with high-quality work (implant- supporting denture, for example)	-	✓
Relatively easy processing	✓	-
Customization possible (for example, firing, stains)	✓	✓
Highest biocompatibility	-	✓
High availability	✓	✓
No change of the joint because of the wear of the teeth	-	✓
None or only very slight wear	-	✓
None or only very slight discolorations	-	✓
Not much advice needed	✓	-
Little pre-planning necessary	✓	-
Long-term stabilization of the tooth position in the oral situation	-	✓

A direct comparison of the concepts:

Concept/Approach	Anatomical concept according to Gysi	Non-anatomical concept according to Sears	„Lingualized Occlusion Concept“ according to Payne
Arrangement	<ul style="list-style-type: none"> • Orientation to shape and size of natural dental elements: • Anatomical molars with cusp inclinations of 33 degrees 	<p>Monoplane occlusion, i.e. totally flat molars without cusps</p>	<ul style="list-style-type: none"> • For the upper denture: Molars with 33 degree inclination and relatively large palatal cusps • For the lower denture: Molars with less steep inclinations (Optiform, f. ex.)
Aesthetics	<p>Very aesthetic</p>	<p>Aesthetic cutbacks in favor of other characteristics</p>	<p>Highly aesthetic, as the natural shape of the elements is maintained</p>
Wearing comfort and chewing characteristics	<ul style="list-style-type: none"> • Easy grinding of the food between the steep cusps of the artificial elements • Sits comfortably because the patient finds a point of support in the (central) occlusion. 	<ul style="list-style-type: none"> • The dental prosthesis can absorb horizontal forces better owing to the missing cusp inclinations. • Thus, the concept is especially suitable for patients with severely resorbed jaws and offers a more stable dental prosthesis 	<ul style="list-style-type: none"> • The steeply inclined cusps allow good food penetration. • The concept of the basin-shaped fossae in the lower elements prevents large horizontal movements on the lower denture when there are lateral movements. • This is especially suitable when horizontal forces are unwelcome, f. ex., with severely resorbed jaws, flabby ridges, parafunctions or implants.
Ease of processing	<p>Relative easy arrangement and adjustment</p>	<ul style="list-style-type: none"> • Easy adjustment of the occlusion with changes in the horizontal and/or perpendicular connection as a result of a resorption of the processus alveolaris • Plus easy processing with „Angle's“ Verbindungen (Class II and III) and when the connection to the jaw makes a cross bite 	<ul style="list-style-type: none"> • Because there is only one contact plane between the lower and upper elements, they are relatively easy to arrange, check and repair. • Furthermore, grinding processes are simplified and a cross bite arrangement prevented



“Lingualized Occlusion Concept” & Continuing education courses

The “**Lingualized Occlusion Concept**” as developed by **enta** in collaboration with leading scientists and universities describes a specific contouring and arrangement of the artificial elements.

It combines the advantages of the so far highly mutual controversial concepts (anatomic and monoplane concept).

The “**Lingualized Occlusion Concept**” achieves optimal results with regard to aesthetics, wearing comfort, chewing characteristics as well as stability, bio- compatibility, easy arrangement and high ease of processing.

enta also offers sound continuing education courses for interested dentists and dental technicians. More information and current dates are found in.



▲ Contact points between upper and lower elements in the „**Lingualized Occlusion Concept**” concept:

To achieve this, the central fossae of the lower elements are widened through selective grinding; during arrangement, the upper elements are edged somewhat towards the oral cavity



Characteristics & market positioning

There are only very few suppliers in the ceramic denture sector. Over 60 years of experience and passionate R&D makes enta a unique specialist in this field.

enta's product range is geared especially to large laboratories that combine the highest demands on quality, ease of processing and personal service with an eye for profitability.

Your enta contact

If you still have any questions or would like to arrange for a personal consulting session? We will look forward to your call or e-mail.



Enta b.v.
Brouwerijbaan 10
4615 AA Bergen op Zoom
The Netherlands

Phone: +31 (0)164.23.74.54
Fax: +31 (0)164.23.53.12
E-mail: info@enta.nl
Internet: www.enta.nl