# PERMAPLAST LH FLOW

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## GB INSTRUCTIONS FLOWABLE COMPOSITE

M+W Permaplast LH Flow is a light curing, flowable, radiopaque, low viscosity composite. Complies with the guidelines described in EN 24049.

### Composition

Monomer matrix: diurethane dimethacrylate, butanediol dimethacrylate, isopropylidene-bis[2(3-hydroxy-3(2)-(4-phenoxy)propyl)bismethacrylate

Total filler: 55 % by weight glass filler (mean particle size: 0.7 µm), pyrogenic silicic acid (mean particle size: 0.012 µm)

### Indications

- Fissure sealing
- Extended fissure sealing on molars and premolars
- Fillings in Black's class V cavities (cervical caries, eroded areas in roots, wedge-shaped defects)
- Minimally invasive fillings in Black's class I and II cavities in areas not exposed to severe occlusal loads
- Minimally invasive fillings in Black's class II cavities
- Restoring defects in enamel
- Blocking out undercuts
- Minimal adjustments to the contours and shade of the enamel

### Contraindications

Unpolymerized composite may cause skin allergies. The user must take adequate precautions (e.g. gloves). Should the operator or patient be known to be allergic to one of the constituents listed under „Composition“, do not use this material.

### Side-effects

To prevent possible reactions of the pulp in cavities where the dentine is exposed, the pulp must be protected adequately (e.g. calcium hydroxide preparation).

### Interaction with other substances

As phenolic substances (such as eugenol) inhibit polymerization, do not use cavity liners containing such substances.

### Application

#### Preparatory measures

Before commencing the treatment, clean the tooth with non-fluoride polishing paste. Use a Vita® shade guide to select the shade while the tooth is still moist.

#### 1. Cavity preparation

Prepare the cavity minimally invasively as generally required for adhesive techniques. All enamel margins in the anterior region must be bevelled. Do not bevel the margins in the posterior region and avoid slice preparations. Spray the cavity with water to clean it, remove all debris and dry it. The cavity must be isolated. It is advisable to place a rubber dam.

#### 2. Pulp protection / Cavity liner

If an enamel-dentine adhesive is used, no cavity liner is required. In very deep cavities those areas in close proximity to the pulp must be coated with a calcium hydroxide material.

#### 3. Approximal contact areas

For cavities with approximal sections, place a transparent matrix and fix it in place.

### 4. Etching

Apply Etching gel absolutely precisely, beginning with the enamel margins. Allow it to react for at least 15 seconds before spraying thoroughly with water and drying the etched surfaces with oil-free compressed air. Avoid dehydrating the dentine. After drying, the etched surfaces of the enamel appear chalky/white and must not be contaminated prior to applying the bonding agent. Should they be contaminated with saliva, rinse and dry them again – re-etch if necessary.

### 5. Applying the bonding agent

Refer to the Bond LC instructions for details.

### 6. Applying M+W Permaplast LH Flow

Place thin layers (max. 2 mm) of Flowable composite directly in the cavity using the curved application tips supplied with the material – for hygienic reasons, they are for single-use only. When applying the material with the tip, ensure that no air bubbles become entrapped. Ensure that the prepared tooth surfaces are wetted thoroughly. Light cure each layer with a commercially available polymerization unit (e.g. M+W Superlite) for 40 seconds, holding the light guide as close as possible to the surface and must not be touched or removed. If forms the bond with the layer of composite applied later.

### 7. Trimming

M+W Permaplast LH Flow can be trimmed and polished immediately after curing using finishing diamonds, flexible disks, silicone polishers and polishing brushes.

### Please note

- When placing time consuming restorations to prevent the composite curing prematurely the dental light should be moved away from the site temporarily or the composite covered with foil impervious to light.
- For hygienic reasons, the curved application tips supplied with the material must only be used once!
- Use a light curing unit with an emission spectrum of 350 – 500 nm for curing this material. As the required physical properties can only be achieved if the lamp is functioning correctly, its luminous intensity must be checked regularly as described by the manufacturer.

### Storage

Do not store above 25°C. Avoid direct sunlight. Close the syringes tightly immediately after use. The material should be at room temperature before use. Retract the plunger of the syringe slightly to prevent the apertures becoming blocked.

Do not use after the expiry date (refer to label on syringe). For use by dentists only. Keep out of reach of children. This product was developed specifically for the described range of applications. It must be used as described in the instructions. The manufacturer is not liable for damage caused by handling or processing the material incorrectly.

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Müller & Weygandt GmbH • Reichardsweide 40 • D-63654 Büdingen • Tel.: +49 (0)6042 – 88 00 88