Brugsanvisning
Käyttöohjeet
Bruksanvisning
Productinformatie
Οδηγίες Χρήσεως
Kullanma Talimatı
Инструкция
Instrukcja stosowania

Instructions for Use
Gebrauchsinformation
Mode d'emploi
Istruzioni d'uso
Instrucciones de uso
Instruções de Uso
Bruksanvisning

- Rodstift af glasfiberforsterket komposit
- Lasikuituvahvisteinen yhdistelmämuovinen juurikanavanasta
- Rotstift av glassfiberforsterket komposit
- Wortelstift van glasvezelversterkt composit
- Ενδονοριακοί αδενες σύνθετης ρητίνης ενισχυμένης με ίνες γυαλιού
- Cam lifilye güçlendirilmiş kompozit kök kanal postu
- Корневой штифт из стекловолокна, усиленный композитом
- Wkłady korzeniowe z materiału kompozytowego wzmocnionego włóknami szklanymi


- Glass-fibre reinforced composite root canal post
- Wurzelstift aus glasfaserverstärktem Composite
- Tenon radiculaire en composite renforcé aux fibres de verre
- Perno radicolare in composito rinforzato con fibre di vetro
- Espigas radiculares de composite reforzado com fibra de vidro
- Poste intra-radicalar de composito reforçado com fibra de vidro
- Rotkanalsstift av glassfiberförstärkt komposit

Instructions for Use


Description

FRC Postec® Plus is a radiopaque root canal post made of glass-fibre-reinforced composite with a conicity of 5°/18°. The adhesive technique used for cementation provides a durable bond between FRC Postec Plus, the core build-up and tooth. The elastic behaviour is comparable to that of dentin. The high translucency allows the fabrication of highly aesthetic, metal-free restorations.

Size 0

Size 1

Size 3

Reamer

High-quality rotary instruments substantially contribute to the success of the treatment. The dimensions of the reamers correspond to the sizes of the FRC Postec Plus posts and thus ensure accuracy of fit in the root canal. For easier identification, the instruments are colour-coded:

- blue → Size 0
- white → Size 1
- red → Size 3

Composition

FRC Postec Plus is a light-conducting, radiopaque root canal post made of glass fibres. The polymer matrix is composed of aromatic and aliphatic dimethacrylates. It also contains ytterbium trifluoride.

Indication

- **Anterior region**
 Size 0: Teeth # 32 - 42
 Size 1: Teeth # 12, 22, 32 - 42
 Size 3: Teeth # 11, 13, 21, 23, 33, 43
- **Posterior region**
 Depending on the clinical situation, both Size 1 and Size 3 can be used in the posterior region.

Contraindication

- Patients suffering from bruxism or suspected bruxism
- Deep overbite
- Less than 2-3 mm supragingival tooth structure
- Circular, isogingival damage
- If patients are known to be allergic to any of the ingredients, FRC Postec Plus should not be used.

Side effects and other interactions

None known to date.

Clinical instructions

For more detailed information, please refer to the separate Instructions for Use of the mentioned products.

1. Preparing the tooth

- Successfully completed endodontic treatment
- Select the FRC Postec Plus suitable for the tooth to be restored
- Prepare the coronal residual tooth structure to accommodate the crown (minimum dentin wall thickness ≥ 1.5 mm) and determine the reconstruction method (direct or indirect technique by means of a post and core build-up) depending on the residual dentin structure.
- Determine the required preparation depth and mark the distance on the corresponding instruments by using silicone stoppers.
 - **The remaining root canal filling should demonstrate an apical length no shorter than 4 mm.**
 - **The length of the post within the root canal must at least correspond to the coronal length of the prosthetic restoration.**
- Remove the root canal filling (e.g. Gates-Glidden or Peeso reamer)
- 1,000 - 5,000 rpm (instrument speed, green-angled handpiece)
- For posts of Size 0, use instruments with a maximum diameter of 0.5 mm.
- For posts of Size 1, use instruments with a maximum diameter of 0.7 mm.
- For posts of Size 3, use instruments with a maximum diameter of 0.9 mm.
- Prepare the canal down to the required depth using the reamer.
- 1,000 - 5,000 rpm (instrument speed, green-angled handpiece)
- Rinse and clean the reamer during preparation.
- Final preparation of the tooth to accommodate the post and core build-up. If necessary, apply antirotation preparation (e.g. retention box, canal inlay)
- Rinse the root canal (e.g. H2O2, 3 %)

2. Further procedure for the direct method

2.1 Try-in and conditioning of FRC Postec Plus

- Check proper fit of the post.
- Determine the coronal length and shorten the post using rotary diamond grinders (extraorally), if necessary. Do not use scissors or cutting pliers!
- After try-in, clean the post with phosphoric acid etching gel (e.g. Total Etch). Allow the etching gel to react for 60 seconds, then thoroughly rinse with water and dry.
- Silanate the post (e.g. Monobond-S). Allow a reaction time of 60 seconds, then carefully dry with an air syringe. Do not touch the surface with your fingers after that.

2.2 Adhesive cementation with a dual- or self-curing luting composite

- Use for example Multilink Automix, MultiCore Flow or Variolink II
- For conditioning the dentin and cementing the post, please refer to the respective Instructions for Use of the cementation system in use. (If a dual-curing luting composite is used, polymerize from the occlusal surface for at least 60 seconds holding the curing light close to the post).
- **Important: When luting composite is used in combination with an adhesive system, the luting composite should not be applied to the root canal which has previously been coated with adhesive. The luting composite and adhesive excess may intermix, which can lead to premature curing. As a result, it may not be possible to move the post into the desired position.**

2.3 Core build-up

- Directly sculpt the build-up using a highly filled composite such as MultiCore Flow, MultiCore HB or Tetric Ceram.
- For conditioning the dentin and applying the composite, please refer to the Instructions for Use of the corresponding composite.

2.4 Preparing the tooth core according to the planned crown

- Observe circular embrasure of the core.
- The preparation margin must be located approx. 2-3 mm in the natural tooth structure.

3. Impression taking and provisional restoration

- Only use eugenol-free cements to place the temporary, since eugenol may inhibit the polymerization of the luting composite.
- If resin-based temporary cements are used, isolate the preparation with Vaseline to prevent the cement from chemically bonding with the composite core build-up.
- During long temporization periods, the provisional should be re-cemented several times to ensure an effective seal at all times.

4. Next appointment

- Incorporate the permanent restoration (e.g. IPS Empress System) according to the corresponding Instructions for Use.

5. Removing the post in case of a revision

- If the removal of the FRC Postec Plus root canal post is indicated, this can be done using rotary instruments.
- Due to the similar moving behaviour of the rotary cutting instruments within the glass-fibre-reinforced root canal post and the root dentin, this procedure has to be carried out with caution to control the risk of root perforation. Check the individual steps by means of X-rays.
- Working steps:
 - Make an X-ray in the right angle technique to determine the length and size of the post.
 - Remove the build-up component and determine the position of the root canal post
 - Drill a preliminary hole in the centre of the post using an extra-long round bur with a diameter of 0.6 mm.
 - After that, remove the root canal post using the FRC Postec Plus reamer.

Important Note

- If the stipulated clinical indications cannot be observed, the FRC Postec Plus should not be used. Any use for indications other than those stipulated may jeopardize the success of the subsequent prosthetic restoration.
- Do not grind FRC Postec Plus and do not apply mechanical retentions → predetermined breaking points
- Cementation should be exclusively carried out according to the principles of adhesive luting. We recommend Multilink Automix, MultiCore Flow or Variolink II from Ivoclar Vivadent.
- Adhesive cementation requires keeping the bonding surfaces dry. If possible, place a rubber dam. In all other cases, provide relative isolation (cotton rolls, saliva ejector).

Warnings

Grinding FRC Postec Plus generates glass-fibre dust, which may cause itching of the skin. Therefore, the use of gloves is recommended. Avoid the inhalation of grinding dust.

Care instructions

- **Reamer:**
 - The product is supplied in non-sterile form. Please disinfect or sterilize the instruments before each use.
- **FRC Postec Plus:**
 - The root canal post can be disinfected with medical alcohol.
 - The root canal post can be sterilized in the autoclave for 5 minutes at 134 °C/273.2 °F and 2.1 bar.

Store out of children's reach.
For use in dentistry only.

Date information prepared: 05/2010, Rev. 3

Manufacturer

Ivoclar Vivadent AG, FL-9494 Schaan/Liechtenstein

The product has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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Ivoclar Vivadent UK Limited
Ivoclar Vivadent, Inc.