Cercon® ceram Kiss

Directions for Use
 Gebrauchsanweisung
 Mode d'emploi
 Istruzioni per l'uso
 Instrucciones de uso

THE DENTAL
 SOLUTIONS
 COMPANY™

Instructions for use	1	en
Verarbeitungshinweise	6	de
Indications de mise en oeuvre	11	fr
Note per la lavorazione	16	it
Consejos para su utilización	21	es
تعليمات بخصوص التحضير	26	ar
Упътване за употреба	31	bg
Pokyny pro zpracová	36	cs
Bearbejdselsindikationer	41	da
Υποδείξεις επεξεργασίας	46	el
Napomene za obradu	51	hr
Feldolgozási útmutató	56	hu
加工について	61	ja
가공처리 주지사항	66	ko
Apdorojimo nuorodos	71	lt
Norādes sakarā ar pārstrādi	76	lv
Instructies voor de bewerking	81	nl
Informacje dotyczące obróbki	86	pl
Instruções de utilização	91	pt
Instrucțiuni de utilizare	96	ro
Рекомендации по работе с материалом	101	ru
Beredningsanvisningar	106	sv
İşlem talimatları	111	tr
Технологія використання	116	uk
加工说明	121	zh

01

02

03

04

05

06

07

08

09

112	General instructions Cercon ceram Kiss
3	Shade conversion chart/Instructions for use
4	Instructions for use
5	Instructions for use/General recommendations for firing

Date: 2017-09

Cercon® ceram Kiss

Indications for use

Cercon ceram Kiss is designed for use on zirconia (zirconium oxide) in single tooth or bridge type restorations. Applications include both anterior and posterior locations.

Product Information

- Cercon ceram Kiss is a veneering porcelain that was especially and exclusively developed for veneering crown and bridge zirconium oxide frameworks (Y-TZP) with a CTE of 10.1–10.5 $\mu\text{m}/\text{m}\cdot\text{K}$ (25–500 °C), preferably made out of Cercon base.
- Cercon ceram Kiss is also suitable for veneering zirconium oxide frames that were primary pressed over with Cercon ceram press.

Contraindications

- Cercon ceram Kiss is not suitable for veneering titanium or alloy frameworks with a similar CTE.
- Cercon ceram Kiss is not compatible with other porcelains.
- Cercon ceram Kiss is contraindicated in the case of bruxism, or of other parafunctions.
- Furthermore Cercon ceram Kiss is also contraindicated by insufficient occlusal distance.

Warnings

Undesired side effects of these medical products are extremely rarely to be expected with proper processing and application. Immunoreactions (e.g. allergies) and/or local unpleasant perceptions (e.g. taste irritations or oral mucosa irritations), however, cannot completely be excluded on principle. If you know any undesirable side effects – also in doubtful cases –, then please communicate them to us.

In case of hypersensitivity of the patients against Cercon ceram Kiss veneering porcelain or one of the ingredients, this medical product must not be used or only used under strict supervision by the treating physician/dentist.

Known cross-reactions or interactions of the medical product with other medical products or materials existing in the mouth must be taken into account by the physician/dentist when using the medical product.

Please pass on all the above-mentioned information to the treating physician/dentist, when processing this medical product for an individually manufactured object.

- Do not inhale abrasive dusts.
- Paste opaque liquid: Harmful if swallowed

Precautions

When working with these materials, make sure to comply with the Instructions for use and the pertinent Material Safety Data Sheets (MSDS).

- For professional use only

Adverse reactions

We are not aware of any risks and/or side effects related to Cercon ceram Kiss veneering porcelain.

Technical data



- CTE, Dentine (25–500 °C): 9.2 $\mu\text{m}/\text{m}\cdot\text{K}$
- Dental ceramics, type 1, class 1 according to DIN EN ISO 6872

Introduction to the market: February 2005

Transport and storage

- Liquids: Store containers tightly closed at temperatures above 10 °C.
- Protect others and pastes from moisture.

Symbols on product labels

REF	Product code
LOT	Batch or lot number
	Use before
	Consult instructions for use
	Keep dry
	Do not reuse
	Date of manufacture

Combinable liquids

- Paste liner:
Paste opaque liquid
- Shoulders:
Ducera® Liquid Quick
- Dentine/incisal etc.:
Ducera® Liquid SD
Ducera® Liquid Form
- Stains/glazes:
Dentsply Universal Liquid
- Isolation:
Ducera® Sep isolating fluid
Ducera® High Temp isolating fluid

Ceramic furnace

For best results, make sure that the required firing temperatures and times are observed strictly. If necessary, adjust the firing parameters of the ceramic furnace as required.

Names/Explanations of important materials

Power Chroma (PC 1–6)

The Power Chromas are highly chromatic, fluorescent and intensive powders for individual shade adjustments. The power chromas can be used to increase the chroma in the cervical, palatal, and occlusal areas. The powders can be used pure or in a 1:1-mixture. By mixing the Power Chromas with Stand by you can create an excellent mamelon dentine. The shade conversion chart is meant as a guidance.

Stand by

Strongly opalescent, almost transparent multi functional powder. It is possible to use Stand by pure or to mix it with all other powders of the Kiss concept. Hence the Stand by powder has a key function.

Opal-Effekt Sunrise/Opal-Effekt Sunset

Opalescent effect powder for yellow as well orange/reddish incisal areas. Very well suited for increasing the chroma at the 2nd or 3rd dentine firings. Sunset is mainly used for A-shades, whereas Sunrise is used mainly for B-shades. Both powders can be diluted by using Stand by.

Opal-Effekt Sky/Opal-Effekt Ocean

Opalescent effect powders for discreet and for strong, deep blue incisal areas – can be diluted by using Stand by.

Opal-Effekt Fog

Opalescent effect powder for grayish incisal areas – can be diluted by using Stand by.

White Surface

Whitish opalescent effect powder to enhance occlusal cusps in the posterior region as well as for palatal/lingual ridges in the anterior area – can be diluted by using Stand by.

Final Kiss

Low-fusing, transparent correction powder (Final Kiss) – firing temperature 680 °C.

Table 1: Colour conversion chart, Cercon® ceram Kiss

Shade	A1	A2	A3	A3,5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
Standard layering																
Opaque	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x
High Chroma Dentin	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x
Incisal	1	2	3	3	5	1	1	4	6	1	5	5	6	2	4	4
Individual layering																
Opaque	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x
Opaque Orange	For characterizing occlusal, cervical and palatal regions.															
Opaque Bleach	For extremely bright/bleached teeth. Normally used in connection with bleach dentines only.															
Opaque Gum	For gingival aspects.															
Shoulder SM/F SM	1	2	2+3	2+4	3+4	1	1+3	3	3+5	1	1+4	2+4	4	1+4	2+4	3+4
Dentine	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x
Power Chroma 1																
Power Chroma 2																
Power Chroma 3	1+2	2	2+5	3+5	4+6	1	1+3	2+3	3+6	1+6	2+6	3+6	5+6	1+6	2+6	3+6
Power Chroma 4																
Power Chroma 5																
Power Chroma 6																
Flu Inside 1	x	x	Mix	Mix		x	x	Mix		x	Mix			x	Mix	
Flu Inside 2					x				x			x	x			x
Opal incisal 1	x	x	Mix	Mix		x	x	Mix		Mix	Mix			Mix	Mix	
Opal incisal 2					x				x			x	x			x

The Kiss colour wheel shows you the right powders for the different shades at a glance.

1 Preparation of the framework

The zirconium-oxide frameworks should be milled in a small anatomic shape for supporting the porcelain and to ensure an even layering thickness of the veneering material. Fine corrections, e.g. adjustment in the marginal area, are carried out by using fine diamond burs.

To avoid tensions in the Cercon ceram Kiss veneering porcelain, it is necessary to round off any sharp corners and edges of the frame.

To avoid micro-cracks, it is necessary to use a water-cooled turbine for griding the framework. To avoid overheating of the framework material do not use a high grinding pressure.

Frames that have cracks and/or holes are not suitable for further processing. The zirconium oxide frames are sandblasted by using with aluminium oxide (grain size 110–125 µm, pressure 3–3.5 bar). Afterwards the frameworks have to be cleaned by using a steam-cleaner or in a clean ultrasonic bath (see fig. 1).

2 Liner

The orange and gum paste liners can be used for the individual characterization of the framework and for gingival sculpting.

With the dentin bleach, the bleach liner enables light colors to be reproduced.

For using opaque frameworks, single packages of liners are also available in all Vita colors (see Fig. 2).

3 4 5 6 Aesthetic Line Basic

In the standard layering technique you can proceed the same way like the build-up diagram that is already known from Duceragold Kiss and Duceram Kiss.

With two more powders, Dentine and Incisal, it is possible to produce a high aesthetic veneer within the shortest time. Please follow figs. 3–6.

First, build up the dentine core, afterwards use the Cut-back-Technique to obtain room for the incisal material. Then complete the restoration with the corresponding incisal powder and proceed with the first dentine firing (see chart 2).

Adjust the shrinkage of the first bake by using the respective dentine and incisal materials. The second dentine firing follows (see chart 2) and thereafter finish the restoration as usual. Finally, continue with the glaze firing (see chart 2). It is possible to do the glaze with or without glaze powder as well as some characteristics by using the corresponding stains (see figs. 3–6).

7 8 9 The porcelain shoulder

Minor corrections may be made using Final

Shoulder (F SM 1–5) following the glaze firing

- The final shoulder materials can be used for repairs with a low firing temperature after assigning the color. The shoulder materials for fabricating a shoulder are also available individually (Fig. 7–9).

10 11 12 13 14 15 Individual aesthetic line

With the individual layering technique, by using the power-chroma and the opal-effect powders, you have the possibility to create natural and highly sophisticated restorations. In accordance with an individual custom shade you can dispense the liner, as described above, and build up directly on the framework.

But keep in mind, if you use Cercon base without the liner or only one liner firing that you increase the value and on the other hand by using Cercon base colored without liner you increase the chroma of dark shades. Beside that consider the change of the basic shade and the transparency by working without liner.

When working without any liner, apply a high fluorescent and high chromatic material, first, like the power chroma. Afterwards proceed with the build-up of the dentine.

After the Cut-back, complete the incisal area with the corresponding opal effect enamel as well as with e.g. in the incisal areas OE Sky or Ocean and for the cervical and body area by using OE Sunrise or Sunset, the opal effect powders (see figs. 10–15).

General information on High Chroma Dentin

Cercon ceram Kiss High Chroma Dentins expand the product range of Cercon ceram Kiss veneer ceramics to include the option of veneering the highly translucent Cercon ht framework material in a simple dentin-incisal stratification in the colors A1-D4.

The general processing instructions for Cercon ceram Kiss apply.

More information is found in the product description / instructions for use for Cercon® ceram Kiss.

Smooth Chroma

With Smooth Chroma, there is now a simple instrument for reducing the color-defining framework intensity if required. This can be necessary for pontics or for strongly discolored stumps. The firing process is identical to dentins.

Preparing for cementation

- Before the final cementation, sandblast the interior carefully with Alu-oxide 50 –110 µm at approx. 2.5 bar. In case of adhesive luting you have to etch the Cercon ceram Kiss, respectively the layered porcelain margin area.

Information on working after glaze firing

- If reworking of Cercon ceram Kiss veneering ceramic should be necessary after glaze firing (grinding in the occlusion/contact point), the whole restoration should subsequently undergo glaze firing once more. This procedure can sinter-seal any possible flaws on the surface. In the case of intra-oral re-working of a restoration cemented in position, it is vital for the corresponding surfaces to be polished off with a suitable ceramic rubber burnisher.

Table 2: General recommendations for firing – Cercon® ceram Kiss

Firing	Pre-heating °C	Drying time min	Heating rate °C/min	Final temp. °C	Holding time min	Vacuum hPa	Long-term cooling min
Paste liner 1*	575	8:00	55	970	1:00	50	–
Paste liner 2*	575	8:00	55	960	1:00	50	–
Shoulder 1*	450	6:00	55	850	1:00	50	–
Shoulder 2*	450	6:00	55	850	2:00	50	–
Dentine 1	450	5:00	55	830	1:30	50	–
Dentine 2	450	5:00	55	820	1:30	50	–
Glaze	450	3:00	55	800	1:00	–	6:00
Correction (Final-Kiss)	450	5:00	55	680	1:00	50	6:00
Final Shoulder (F-SM)	450	5:00	55	680	1:00	50	6:00

* Optional – only if needed

Slow cooling is mandatory after the last firing; this includes correction firings of restorations after try-in.

The firing temperature must be adapted to the number of units fired in the same cycle. 5 to 9 units require an increase by 5 °C to 10 °C; 10 or more units require an increase by 10 °C to 20 °C.

The values indicated here are recommended values and serve only as a clue. Deviations of the burning results are possible. The burning results depend on the respective output per oven and are due to the manufacturer and to age. Therefore the recommended values must be individually adapted at each burning. We recommend a test burning to control the oven. All indications have been carefully elaborated and tested by us, but they are passed on without any guarantee.

For up-to-date firing recommendations please visit www.dentsplysirona.com

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DeguDent GmbH
Dentsply International Inc.