

### Data Sheet

- 01 Product Name**  
TriLor® Arch Metal Free Dental System
- 02 Material Description**  
Fiber Reinforced Composite (FRC)
- 03 Manufacturer**  
Bioloren S.r.l.  
Via Alessandro Volta 59  
21047 Saronno (VA), Italy  
Phone: +39 029 67 6703261  
Email: info@bioloren.com  
Web: www.bioloren.com
- 04 Exclusive Agency**  
Harvest Dental Products, LLC  
905 Columbia Street  
Brea, California 92821  
Phone: +1 714 674 7400  
Email: hello@harvestdental.com  
Web: www.harvestdental.com
- 05 Product Description**  
TriLor- Arch is a new techno-polymer for fabricating implant supported bars and structures for definitive restorations on implants. Made of a thermoset resin with integrated multi-directional Fiber-Reinforced-Composite (FRC).
- 06 Product Character**  
Biocompatible  
Strong  
Lightweight
- 07 Material Properties**
- |                             |     |        |
|-----------------------------|-----|--------|
| Tensile strength            | 380 | MPa    |
| Flexural strength           | 540 | MPa    |
| Tensile elongation          | 2   | %      |
| Flexural modulus            | 26  | GPa    |
| Tensile modulus             | 26  | GPa    |
| Compressive strength        | 530 | MPa    |
| Charpy impact strength      | 300 | KJ/cm2 |
| Rockwell hardness (scale R) | 111 | HRR    |
| Barcol hardness             | 70  |        |
| Shore D. hardness           | 90  |        |
| Density                     | 1,8 | g/cm3  |
- 08 Biocompatibility Testing**
- |                        |                               |
|------------------------|-------------------------------|
| Genotoxicity           | ISO 10993-3                   |
| Acute toxicity         | ISO 10993-11:2006             |
| Hypersensitivity       | ISO 10993-10:2010             |
| Animal skin irritation | ISO 10993-10:2010             |
| Water absorption       | ISO 10477-2009                |
| Cytotoxicity           | ISO 10993-5;2009,10093-5:2000 |
- 09 Mechanical Testing**
- |                    |                |
|--------------------|----------------|
| Flexural hardness  | ISO 14125:2000 |
| Fracture toughness | ISO 14125:2008 |
- 10 Contraindications**  
Insufficient oral hygiene.  
High temperature process: techniques or protocols subject To temperatures above 150° C.

### Instructions for Use

- 11 Dental Indications**  
TriLor- Arch allows the analog creation of connection structures between implants that can be embedded in removable prosthetic devices, as well as frameworks for transferring the position of the implants into the mouth, or simple reinforcement bars for total definitive prosthesis.
- Full Arch
  - Partial or total mobile prostheses with attachments
  - Reinforced bar on implants for immediate or deferred load
- 12 Trilor® Arch produced in the following thickness:**  
NOTE: This arch cannot be modified in thickness, but in width only.)
- Arch Thickness (mm)**
- H 3,5 mm**  
Suggested for the solidarization of transfers and for immediate structures
- H 5,5 mm**  
Standard height, used for the most part for prosthetic frameworks.
- H 7,5 mm**  
Suggested in those cases where it is necessary to model the framework due to angled implants or implants positioned on different heights.
- 13 Framework Design**  
Since the most important use of TriLor- Arch is the creation of a connection framework among the implants, it is required that the working model have the analogs (indirect system of implants) in detected position with abutment inserted and screwed.
- It is important to choose the right thickness of TriLor- Arch according to the characteristics of the case and according to the indications.
- It is suggested to detect the distances among the abutments using a baseplate wax strip in order to fix the exact positions of the abutments by marking a hole through the wax.
- Record the detected position, overlaying the wax strip on TriLor- Arch and write with a pencil directly on the material. Position the wax strip over TriLor- Arch, mark the holes on underlying structure with a pencil.
- Pierce TriLor-Arch to create the housing of implant abutments according to the needs of shape and length decided for the case using tungsten ball drills. The more accurate the passage hole, the more stable and resistant the structure will be.

# TriLor® Arch Metal Free Dental System

## Fiber Reinforced Composite



### Important Warning

Before cutting of the abutment, verify that the height of abutments or of the cannula is the same or greater than the thickness of the Trilor Arch.

Draw with a pencil the shape of the bar on the arch.

Model with carbide tungsten drills using a micro-motor at max speed of 15000 rpm. To speed up this roughing phase you can also use the trimmer, finishing the hand-piece with carbide tungsten drills.

The Trilor-Arch can also be used partially (cut in sections) for the creation of structures in separated occlusal sectors (posterior or anterior).

### 14 Suggested Thicknesses

The minimum thickness of the free part between the two abutments shall be 7 mm<sup>2</sup> (3,5 x 2mm) (A).

The possible cantilever (C) shall be at maximum 10mm.

The minimum thickness near the holes for cementation with cannulas shall not go under 0,8mm (B).



Photo, courtesy of mdt. E. Riccomini

### 15 Surface Preparation

1. Sandblast the surface of Trilor® Arch using disposable aluminum dioxide at 110 micron at 2 bars of pressure
2. Clean with light air blasts (dry, no oil)
3. Clean with ethyl alcohol
4. Treat with silane and let evaporate for a few minutes
5. To link with metal components (abutments) it is suggested to consolidate Trilor® Arch doing a pre-treatment of surface of metals (sandblasting). Sandblast the abutments and then cement the framework with resin cement, preferably opaque to obtain a better mimetic effect.
6. The possible polishing of exposed areas of Trilor® Arch is carried out with silicon gums like the ones used for composites. Then proceed with the final finishing with the diamond paste.
7. Do not exceed 150° C
8. In case of total removable prostheses, having finished the modeling, bonding and superficial finishing, proceed directly tot the assembly of the aesthetic parts (resin teeth, ceramic teeth and so on) that can be linked to the Trilor® Arch structure using the pressure method or pouring with masks.

### 16 Batch/ Lot Number Traceability

Batch numbers of TriLor® discs are located on box label. Document Lot No. to ensure traceability.

### 17 Preventative Measure

Wear a protective dust mask, gloves, protective eye glasses, combined with a good suction system. After any treatment or processing, the prosthetic framework has to be cleaned and disinfected according to the national guidelines before being put "in situ". The Trilor Arch can be sterilized in autoclave at 121° C for 15 minutes or 131°C for 5 minutes.

### 18 Side Effects

No known side effects if used as intended.

### 19 Storage

Avoid exposure to bright light. Keep away from heat sources

### 20 Disposal

Can be disposed of with normal household garbage.

### 21 Ordering information

- 50000-35 TriLor® Arch 3,5mm 3/pk
- 50000-55 TriLor® Arch 5,5mm 3/pk
- 50000-75 TriLor® Arch 7,5mm 3/pk



### 22. Contact Harvest Dental

If you would like to place an order, or need technical Support, please contact us using the info below.

Harvest Dental Products  
800.706.7599  
hello@harvestdental.com  
[www.harvestdental.com](http://www.harvestdental.com)

**Thank you for reading. We appreciate your time.**