Harvard ZirconCore

Manufacturer

Harvard Dental International GmbH
 Margaretenstr. 2 – 4
 15366 Hoppegarten, Germany

Tel: +49 (0) 30/99 28 978-0
 Fax: +49 (0) 30/99 28 978-19
 info@harvard-dental-international.de
 www.harvard-dental.de

CE 0482
 Made in Germany

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Instructions for Use EN

Core build-up and post cementing composite

Harvard ZirconCore is a dual cure high radiopaque flowable microhybrid composite material with nano particles (nano-zirconium dioxide, nano-calciumfluoride) for core-build-up and cementation of posts. Due to its excellent mechanical properties final crown preparation can be carried out more precise. Harvard ZirconCore is also suitable for cementing of crowns and bridges, inlays and onlays if there are no high esthetic requirements.

Delivered in 1:1 Minimix-syringes Harvard ZirconCore can be easily dispensed and applied directly. Harvard ZirconCore exhibits a short setting time without high heat generation. The dual cure properties enables the dentist also to carry out cementations and core build-ups in cases where a light cure cannot be guaranteed to be sufficient.

Indications

Material for core build-ups and cementing of posts

Contraindications

The placement of Harvard ZirconCore is contraindicated if a dry working area or the recommended application technique are not possible. Also do not use Harvard ZirconCore if the patient is known to be allergic to any of the ingredients.

Side effects

In singular cases, Harvard ZirconCore may cause a sensitizing reaction in patients with a hypersensitivity to any of the ingredients. In these cases, the material should not be used. Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection areas close to the pulp should be covered with a thin layer of calciumhydroxide material.

Incompatibility with other materials

Do not use in combination with substances containing eugenol because eugenol inhibits the polymerization of the composite. Neither store the composite material in proximity of eugenol containing products, nor let the composite allow coming into contact with materials containing eugenol.

Application

1. Preparing the Minimix-Syringe

Remove the cap of the Minimix-syringe and throw it away (**do not use it again!**). It is replaced by a special 1:1 mixing cannula. Turn the cannula 90° until it locks in position. The Material is now ready for application.

2. Application

2.1. Isolation

Use of a rubber dam to isolate the tooth is strongly recommended.

2.2. Root canal preparation

Refer to directions of the selected post manufacturer.

Remove all existing old restorations and decay from the tooth. Prepare and clean the root canal with e.g. sodium hypochlorite solution, rinse and remove excess solution from the canal with a soft paper tip.

Conditioning of the root canal should be done with the self etching bonding agent Harvard Bond SE Dual:

One drop of Harvard Bond SE Dual PART A and one drop of Harvard Bond SE Dual PART B were combined in a mixing pallet and mixed for 5-10 seconds.

Note:
Do not interchange lids of the bottles, because this can lead to a cross-contamination of the liquids.

Apply the homogeneous mixture generously with a brush onto the slightly wet root canal walls for 15 seconds with agitation. Repeat procedure 1 – 2 times. All the dentine surfaces must kept wet with the primer over the specified time

Remove excess Harvard Bond SE Dual with dry paper points. Dry the root canal cautiously with oil free air for 15 seconds to remove all volatile components and to disperse the adhesive to an even layer. Light cure all areas that are available for a dental curing unit for 20 seconds. For a maximal adhesion it is strongly recommended to use the light cure mode.

If light cure is absolutely impossible, the adhesive will also cure in the autocure mode. After application of the adhesive air thin to remove all volatile components. There must remain a sticky layer.

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