Harvard Ionoresin Cem LC

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Harvard Ionoresin Cem LC

Resin modified glass ionomer cement for luting and lining, light-curing

Harvard Ionoresin Cem LC is a light-curing, resin modified glass ionomer cement for luting and lining.

As a glass ionomer cement it is simple to use. The curing is initiated by common visible light. Harvard Ionoresin Cem LC exhibits strong chemical adhesion bond to teeth and good biocompatibility. This glass ionomer cement is characterised by good mechanical properties.

Indications

Cementing of crowns, bridges, inlays and onlays (all types: metal, resin, ceramic fused to metal and ceramic)

Contraindications / side-effects

- In singular cases, the material may cause a sensitizing reaction in patients with a hypersensitivity to any of the ingredients. In these cases, the material should not be used further on.
- Direct and indirect pulp capping

Interaction with other materials

Avoid direct contact with products containing eugenol since eugenol impairs the setting of Harvard Ionoresin Cem LC.

Do not mix the powder or liquid of Harvard Ionoresin Cem LC with any other glass ionomer product.

Application

1. Preparation

Prepare the tooth to be restored in the usual manner.

For pulp protection areas close to the pulp should be covered with a thin layer of calcium hydroxide material (e.g. Harvard CalciumHydroxide) or MTA cement (e.g. Harvard MTA Universal). Refer to the instructions of the product used.

2. Dosing of powder and liquid

The powder to liquid ratio is 1.8 / 1.0. This can be obtained by mixing one level (orange) scoop of powder and 2 drops of liquid.

Note:

The spoon dosage is only an approximate reproduction of the nominal mixing ratio.

For accurate dispensing of Harvard Ionoresin Cem LC powder shake the bottle to loosen the powder. Overfill the spoon with the powder and level the powder using the scraper at the top of the bottle. Avoid compressing powder into the spoon. Place the powder on the mixing pad.

If any bubbles are present in the Harvard Ionoresin Cem LC liquid, lightly tap the bottle with the fingers holding it. Discard drops that contain bubbles and are obviously not full-sized.

For dispensing of Harvard Ionoresin Cem LC liquid turn the bottle vertically with the tip about 5 cm above the mixing pad. Steady your hand and squeeze the bottle gently to dispense one drop.

After use, tightly close both liquid and powder bottles to prevent exposure to moisture.

3. Mixing of powder and liquid

Use a small spatula to rapidly mix all the cement powder into the liquid. Usually one scoop powder/ 2 drops liquid should provide for sufficient amount of mixed cement. The mixed cement should have a homogenous consistency. Total mixing time is **30 seconds**.

4. Cementing

Avoid water and saliva contamination during application and setting of the cement. For a dry working area adequate isolation of the tooth is required.

Prepare the restoration according to the manufacturer instructions.

Mix the required amount of Harvard Ionoresin Cem LC. Apply the mixed cement into the lumens of the crown or bridge (approx. half) resp. to the bonding surface of inlays or onlays in a thin layer. Insert the restoration immediately.

Working time from start of mix at 23°C (74°F): **1:30 minuten**

Net setting time without any light is approx. 5:30 minutes (at 37°C / 99 °F).

Remove excess cement at the first setting stage.

Maintain isolation until the set of the cement is verified (ca. 5:30 minutes). Setting can be optimized with **20 seconds** light cure with a suitable dental light cure unit (wavelength range 400–500 nm, light intensity min. 1000 mW/cm²).

Note: Higher temperatures will shorten the working time, lower temperatures will prolong the working time.

An overextended working time will cause the loss of adhesion to the enamel and the dentine.

Notes/Warnings

- The material can also be used as liner under composite restorations.
- Unpolymerized material may have an irritant effect and can lead to sensitization against methacrylates.
- Avoid contact with skin, mucous membrane and eyes
- If the material comes into contact with skin, remove the material with absorbent cotton soaked in alcohol and rinse with water. If the material comes into contact with eyes, immediately rinse with copious amounts of water and seek medical advice if required.
- Commercial medical gloves do not protect against the sensitizing effect of methacrylates.

Storage information

Storage temperature 4-25°C (39-77°F). Protect from direct sunlight. Do not use after expiry date.

Warranty

Harvard Dental International GmbH warrants this product will be free from defects in material and manufacture. Harvard Dental International GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and Harvard Dental International GmbH's sole obligation shall be repair or replacement of the Harvard Dental International GmbH product.

Limitation of Liability

Except where prohibited by law, Harvard Dental International GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

For dental use only!

Store product out of reach of children!

Harvard Ionoresin Cem LC conforms to the DIN EN ISO 9917-2.

Keep the instructions for use for the duration of the application.

Order-No.	Article
7041216	15 g powder, shade Universal / 10 ml liquid dosage spoon, mixing pad