

Processing instructions

1. Intended purpose

The light-curing denture lacquer is used to chemically smooth the surface of a denture base after finishing.

2. Description of product and users

2.1 Brief description of product

Light-curing single-component material for sealing surfaces of denture bases, and primer for preparation of surfaces.

2.2 Users

For use by laboratory technicians in a dental laboratory.

3 Composition

New Outline Skin Primer: Methyl methacrylate

New Outline Skin glaze: Methyl methacrylate (MMA), urethane acrylates, initiators, stabilisers

4 Indications

- for processing surfaces of denture acrylics
- for processing surfaces of tray materials
- for processing surfaces of temporary crowns and bridge acrylics
- for processing surfaces of reworked and worn dentures

5 Contraindications

If a patient has allergies or hypersensitivities to a component of this product, it should not be used or used only under the strict supervision of the attending physician/dentist.

6 Warnings

New Outline Skin Primer

Danger. Contains methyl methacrylate. May cause an allergic skin reaction.

New Outline Skin Glaze

Danger. Liquid contains methyl methacrylate, 2-propenoic acid, reaction product with dipentaerythritol, diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide. May cause an allergic skin reaction. Harmful to aquatic life with long lasting effects.

7 Safety instructions

New Outline Skin Primer

Avoid breathing vapours / spray.

8 Interactions with other medicinal products

Phenolic substances such as Eugenol inhibit polymerisation. Therefore, do not use any material containing these substances.

9 Application / Preparation

Roughen the area of the surface to be lacquered prior to use or blast and clean with 50 µm aluminium oxide, remove adhesive dust. Pretreat the roughened surface with a cleaner according to the manufacturer's instructions. Apply a thin brush coat of

New Outline Skin Primer and leave to act for 4 min. Apply a thin coat of New Outline Skin Glaze and cure.

10 Notes on processing

A light polymerisation unit with an emission spectrum of at least 310-500 nm should be used. The physical characteristics required can only be achieved if the correct lamps are used. Regular checks of the light intensity in accordance with the manufacturer's specifications are therefore required. Do not use intermediate polymerisation units.

11 Troubleshooting / FAQ List

Fault	Cause	Corrective action
does not solidify	inadequate polymerisation	observe the polymerisation times
		- check lamp and replace if necessary - do not use intermediate polymerisation lamps
surface greasy	inadequate polymerisation	observe the polymerisation times
		check device / service device regularly
lacquer flakes off	inadequate adhesion bond	surface must be free of adhesive dust
		Observe the processing times

12 Polymerisation times

Light-curing unit	Time
HiLite / UniXS	90 sec.
Labolight LV-II / III	5 min
Solidilite EX	5 min

Processing instructions

13 Storage

Storage temperature 10-25°C / 50-77°F

Close the bottle tightly immediately after use.

14 Shelf life

The maximum shelf life is printed on the label of each bottle.

Do not use after the expiry date.

15 Side effects

With proper preparation and use of this medical device, adverse effects are extremely rare. Immune reactions (such as allergies) or local discomfort, however, cannot be ruled out completely.

Any serious adverse events associated with the use of this product should be reported to the manufacturer stated below and the competent authority.

16 Disposal

Any residues or packaging materials should be disposed of in accordance with local and/or statutory regulations.