TriVest CAM Refractory disc | Digitized

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Why TriVest CAM will blow your mind?

- FACT 1** While the market for esthetic veneers continues to grow, the standard of analog processing is still not economical at all!
- FACT 2** Digital processes leave little room for creativity and esthetics in the fabrication of veneers!
- FACT 3** The digital process chain enables significantly higher profit margins for dental technology!
- FACT 4** It was about time that HPdent created the digital “Refractory Die” process!

PROCESS WITHIN THE MACHINE

- ▶ The digital die model was designed, 3D printing has started
→ Die/model needs to be designed without undercuts
- ▶ Nest the dies in your CAM software also
- ▶ Please create at least 3 connectors per die
- ▶ Milling with zirconia strategy, if your software allows O-scaling
- ▶ Alternatively, use a composite strategy
→ Attention: always use dedicated ZrO_2 cutters
- ▶ We mill “dry” – but “wet” is also an option
- ▶ Once milled – cut out the dies, grind the connectors and place the dies in the model
- ▶ Possibly rework the anti-twist rim slightly

YOUR JOB

- ▶ Don't you dare doing a degas firing, we already did that for you!
- ▶ Start the wash firing immediately with your favorite ceramic
- ▶ Do the layering and fire your awesome veneer as usual
- ▶ Now sandblast the die with glass beads
- ▶ Place your veneer on the printed die
- ▶ Don't forget to write the invoice

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