SPLINT

UCAN
PRINT

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UCAN PRINT SPLINT

1. Intended use / Indication

The UCAN PRINT SPLINT resin is a light-curable polymerizable resin intended to be used in conjunction with extra-oral curing light equipment for the fabrication, by additive manufacturing, of dental objects such as mouthguards, nightguards, splints, repositioners and retainers.

2. Physical properties*

Shore hardness D (ISO 48-4)**: 80

Ultimate flexural strength MPa (ISO 20795-2)***: ≥ 50

Flexural modulus MPa (ISO 20795-2)***: ≥ 1500

Water absorption $\mu\text{g}/\text{mm}^3$ (ISO 20795-2): ≤ 32

Solubility $\mu\text{g}/\text{mm}^3$ (ISO 20795-2): ≤ 5

3. Ordering information

1 kg: $\lambda \leq 405 \text{ nm}$

transparent

REF: 768486

4. Contraindication

UCAN PRINT SPLINT is contraindicated ...

1. ... if a patient is known to be allergic to any of the ingredients.

2. ... for bruxism splints and splints with wall thickness $\leq 1.5 \text{ mm}$.

3. ... for denture bases.

4. ... for every application that is not part of the indication (see above).

5. Patient Target Group

Persons being treated in the context of a dental procedure.

6. Intended Users

Dentist, dental technicians

* These data come from measurements of a representative sample that was determined as part of our quality assurance.

** According to internal design and requirements specifications
on the basis of ...

*** not applicable

7. Requirements

Software - information available from:

exocad GmbH, Julius-Reiber-Strasse 37, 64293 Darmstadt, Germany

Hardware (3D Printing) - information available from:

ASIGA, 2/19-21 Bourke Road, Alexandria, NSW 2015, Australia

Hardware (Post Curing) - information available from:

pro3dure medical (CD-1, CD-2), Am Burgberg 13, 58642 Iserlohn, Germany

Additional specifications on www.candulor.com.

8. Material

UCAN PRINT SPLINT consists of functional (meth)acrylic resins, initiators, dyes and stabilisers.

9. Geometric presetting

Minimum wall thickness: 1.5 mm

10. Material parameters

Radiation penetration depth controllable by exposure time

50 µm

100 µm

11. Manufacturing process (fig. 1-10)

1. Prepare data (CAD & build preparation).
2. Choose process parameters (build style, etc.).
3. Transfer prepared data to 3D printing device.
4. Prepare 3D printing – shake bottle.
5. Fill resin tank of 3D printing device.
6. Build the parts.
7. Clean parts (with IPA $\geq 97\%$ or equivalent cleaning agent) approx. 4 min. in an ultrasonic bath or equivalent device – precleaning recommended.
8. Dry parts, until there are no residues of IPA or equivalent cleaning agent.
9. Post curing (4 min.): inert atmosphere recommended (use adequate light curing device).
10. Finish parts.

12. Finishing processes

Polishing

13. Notice

Please follow the instructions provided by the software manufacturer in terms of parameter settings and design recommendations. Please follow the instructions provided by the hardware manufacturer in terms of parameter settings/printing & post curing recommendations. To avoid detrimental effects on material quality do not expose the liquid material to irradiation under any circumstances. Deviations from the described manufacturing processes or storing conditions may lead to different mechanical and optical properties of the material. According to the EU Medical Devices Regulation, users/patients are obliged to report serious events with a medical device to the manufacturer and to the competent authority of the country in which they occurred.

Caution: The lot number and the best before date are indicated on each packaging. In case of claims please always indicate the lot number of the product. Do not use the product after expiry of the best before date. Obtain special instructions before use. Do not handle until all safety precautions have been read and understood. Avoid breathing dust/fume/gas/mist/vapours/spray. Contaminated work clothing should not be allowed out of the workplace. Avoid release to the environment. Wear protective gloves/protective clothing/eye protection/face protection.

IF ON SKIN: Wash with plenty of soap and water. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. IF exposed or concerned: Get medical advice/attention. If skin irritation occurs: Get medical advice/attention. If skin irritation or rash occurs: Get medical advice/attention. If eye irritation persists: Get medical advice/attention. Take off contaminated clothing and wash it before reuse. Collect spillage. Store locked up. Dispose of contents/container in accordance with local/regional/national/international regulations.

14. Hazard statements

Causes skin irritation. May cause an allergic skin reaction. Causes serious eye irritation. Suspected of damaging the unborn child. Toxic to aquatic life with long lasting effects.

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pro3dure medical GmbH
Candulor AG
Candulor Dental GmbH