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TA-Clean - high-performance cleaning material for zirconium oxide sintering furnaces

Cleanliness. Efficiency. Perfection.

Product description:

TA-Clean is a specially developed cleaning material for zirconium oxide sintering furnaces. It effectively removes impurities that arise during the sintering process. Regular use extends the service life of your furnace and ensures the long-term quality of your sintered zirconium oxide products.

Contamination of the furnace by foreign materials can lead to unpredictable changes in translucency and color, often whitish opaque green-yellow discoloration. Especially after previous sintering of different zirconium oxide generations or frequent use of coloring liquids. Regular cleaning is recommended for reliably unaltered sintering results.

Application:

1. add approx. 10g cleaning beads to the aluminium oxide crucible and place in the sintering furnace.
2. cleaning firing:
Heating: 55°C/min to 1000°C - 20min holding time.
Cooling: 40°C/min to 300°C.
3. allow the furnace to cool down and carefully remove the crucible. Discoloration of the beads indicates successful cleaning. Only use cleaning beads once.

For optimum cleaning, use regularly after every 10th sintering cycle.