IPS Classic[®]

[en] Instructions for Use

- Metal-ceramic
- Ceramic mixing liquid
- Stains and Glazes

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- Patients with permanent teeth

Intended users

- Dentists (clinical procedure)
- Dental laboratory technicians (fabrication of restorations in the dental laboratory)

Special training

No special training required.

Use

For dental use only.

Description

IPS Classic is a ceramic veneering material containing leucite. It is suitable for the fabrication of metal-ceramic restorations at firing temperatures higher than 900°C (1652°F). The product can be used for veneering alloys in the CTE range of $13.8 - 15.2 \times 10^{-6} / \text{K}$ ($25 - 500^{\circ}\text{C}$).

Product name	Product description
IPS Classic V Opaquer IPS Classic V I-Opaquer IPS Gingiva Opaquer IPS Classic V Powder Opaquer	Masking of metal frameworks and adjustment of the base shade
IPS Opaque Dentin V IPS Classic V Dentin	Veneering of metal frameworks
IPS Classic V Incisal IPS Classic Opal Incisal	Veneering of the occlusal and incisal areas of metal frameworks
IPS Classic V Transparent	Characterization of ceramic veneers on metal frameworks
IPS Impulse Transparent IPS Impulse Molar Incisal IPS Impulse Occlusal Dentin IPS Effect	Characterization of ceramic veneers on metal frameworks
IPS Gingiva	Veneering of metal frameworks in the gingival area
IPS Gingiva Modifier	Characterization of IPS Gingiva materials
IPS Classic Add-On Powder	Adjustment of ceramic restorations on metal frameworks
IPS Shade V	Characterization of IPS Classic restorations in the anterior and posterior region
IPS Stains P	Characterization of IPS Classic restorations in the anterior and posterior region
IPS Classic Glazing Paste	Glazing of IPS Classic restorations in the anterior and posterior region
IPS Classic Glazing/Staining Liquid	Mixing of IPS Classic Glazing Paste for ceramic restorations in the anterior and posterior region
IPS Classic Powder Opaquer Liquid	For mixing IPS Classic V Powder Opaquer
IPS Classic Build-Up Liquid (N, L)	For mixing IPS Classic V, IPS Classic, IPS Impulse and IPS Gingiva layering materials

Indications

- Missing tooth structure in anterior and posterior teeth
- Partial edentulism in the anterior and posterior region

Types of restorations:

- Crowns
- 3-14-unit bridges

Contraindications

The use of this product is contraindicated if the patient is known to be allergic to any of its ingredients.

Limitations of use

- Untreated bruxism (a splint is indicated after incorporation)
- Mixing and processing together with other dental ceramics
- Failure to observe the necessary minimum connector and framework thicknesses
- Exceeding or falling short of the stipulated veneering layer thicknesses
- Failure to observe the layer thickness ratio between the framework and layering ceramics
- Veneering of dental alloys not within the stipulated CTE range
- Veneering of titanium and zirconium oxide frameworks
- The final restoration must not be reused.

Side effects

There are no known side effects to date.

Interactions

There are no known interactions to date.

Clinical benefit

- Reconstruction of chewing function
- Restoration of esthetics

Composition

System components powders and pastes		Ingredients
Opaquer	IPS Classic V Opaquer IPS Classic V Intensive Opaquer IPS Gingiva Opaquer	<ul style="list-style-type: none"> - Leucite glass-ceramic - Glass - Zirconium oxide - Inorganic pigments - Glycols (butandiol, glycerine) - Hydrophilic, pyrogenic silicic acid
	IPS Classic V Powder Opaquer	<ul style="list-style-type: none"> - Leucite glass-ceramic - Zirconium oxide - Inorganic pigments
Gingiva	IPS Gingiva IPS Gingiva Modifier	<ul style="list-style-type: none"> - Leucite glass-ceramic - Opalescent glass - Inorganic pigments
Dentin	IPS Opaque Dentin V IPS Classic V Dentin	<ul style="list-style-type: none"> - Leucite glass-ceramic - Opalescent glass - Inorganic pigments
Incisal	IPS Classic V Incisal IPS Classic Opal Incisal IPS Classic V Transparent	<ul style="list-style-type: none"> - Leucite glass-ceramic - Opalescent glass - Inorganic pigments
Characterization	IPS Impulse Transparent IPS Impulse Occlusal Dentin IPS Impulse Molar Incisal	<ul style="list-style-type: none"> - Leucite glass-ceramic - Opalescent glass - Inorganic pigments
	IPS Effect	<ul style="list-style-type: none"> - Leucite glass-ceramic - Glass - Inorganic pigments
Stains and glazes	IPS Shade V IPS Stains P	<ul style="list-style-type: none"> - Leucite glass-ceramic - Glass - Inorganic pigments - Glycol (propylene glycol) - Hydrophilic, pyrogenic silicic acid
	IPS Classic Glazing Paste	<ul style="list-style-type: none"> - Leucite glass-ceramic - Glass - Glycols (butandiol, glycerine) - Polyvinylpyrrolidone - Highly dispersed silicic acid
Corrections	IPS Classic Add-On Powder	<ul style="list-style-type: none"> - Leucite glass-ceramic - Glass

System components ceramic mixing liquids	Ingredients
IPS Classic Powder Opaquer Liquid	<ul style="list-style-type: none"> - Water - Glycol (propandiol) - Acetic acid - Metallic salt
IPS Classic Build-Up Liquid N (Normal)	<ul style="list-style-type: none"> - Water - Metallic salt
IPS Classic Build-Up Liquid L (Long)	<ul style="list-style-type: none"> - Water - Glycol (butandiol) - Metallic salt
IPS Classic Glazing/Staining Liquid	<ul style="list-style-type: none"> - Glycol (ethylene glycol)

2 Application

Notice

Observe the following specifications:

- Wall thickness for metal frameworks: Crowns min. 0.3 mm; abutment crowns min. 0.5 mm.
- Ceramic layering thickness: min. 0.8 mm; max. 1.7 mm.

Notice

Only use the system components mentioned in chapter Intended use. Proceed as follows:

- Mix powder materials with mixing liquid before use.
- Extrude and mix paste materials before use. Use mixing liquid for diluting.

Notice

Conduct all firing cycles according to the instructions of the manufacturer. Firing parameters can be found in chapter Technical specifications.

Fabricating the framework

Fabricate and oxidize the framework according to the instructions of the alloy manufacturer.

Applying and firing the opaquer

1. Evenly apply paste or powder opaquer on the framework.
2. Conduct opaquer firing.
3. Repeat the steps 1 to 2.

Layering and firing dentin and incisal material

1. If space is limited: Apply opaque dentin material on the opaquerized framework.
2. Completely layer the restoration with dentin material.
3. Cut back the restoration.
4. Add incisal material to the incisal area.
5. Optional: Customize restorations using characterization material.
6. Conduct dentin and incisal firing.
7. Repeat the steps 4 to 6.

Layering and firing gingiva material

1. Layer the gingiva area.
2. Conduct gingiva firing.
3. Repeat the steps 1 to 2.

Characterizing, glazing and firing

1. Create an anatomical shape and surface texture in the tooth- and gingiva-coloured areas.
2. Customize the restoration using stains.
3. Conduct stain firing.
4. Apply glaze material using a brush.
5. Conduct glaze firing.

Applying and firing corrections

Notice! Add-on material is neutral in colour. Colour the add-on material by mixing it with either dentin or incisal material in the ratio of 1:1.

1. Apply add-on material.
2. Conduct corrective firing.

Seating the restoration

Notice! Seating is exclusively done by the dentist.

Seat the restoration using conventional or adhesive methods.

3 Technical specifications

Properties of IPS Classic according to ISO 6872:2015 + Amd 1:2018 and ISO 9693:2019. IPS Classic has been classified as a Type I, Class 1 dental ceramic.

Physico-chemical properties of the components of IPS Classic according to ISO 6872:2015 + Amd 1:2018 and ISO 9693:2019.

Product components/properties	Radioactivity	Biaxial flexural strength σ	Chem. solubility	Adhesive bond tb				
	[Bq/g U ²³⁸]	[MPa]	[$\mu\text{g}/\text{cm}^2$]	[MPa]				
	ISO 6872:2015 + Amd 1:2018			ISO 9693:2019				
IPS Classic V Opaquer IPS Classic V Intensive Opaquer IPS Gingiva Opaquer	≤ 1.0	≥ 50	< 100	> 25				
IPS InLine Powder Opaquer								
IPS Gingiva IPS Gingiva Modifier IPS Opaque Dentin V IPS Classic V Dentin IPS Classic V Incisal IPS Classic Opal Incisal IPS Classic V Transparent IPS Impulse Transparent IPS Impulse Occlusal Dentin IPS Impulse Molar Incisal								
IPS Effect								
IPS Classic Shade V IPS Classic Stains P								
IPS Classic Glazing Paste								
IPS Classic Add-On Powder								
								n.a.

Thermal properties of the components of IPS Classic according to ISO 6872:2015 + Amd 1:2018 and ISO 9693:2019.

Product components/properties	Coefficient of linear thermal expansion (CTE)				Glass transition temperature τ_g	Thermocycling test	
	2 Firing cycles [10 ⁻⁶ /K]	4 Firing cycles [10 ⁻⁶ /K]	Average value [10 ⁻⁶ /K]	Measuring range [°C]			[°C]
						ISO 9693:2019	
IPS Classic V Opaquer IPS Classic V Intensive Opaquer IPS Gingiva Opaquer	13.6 ± 0.5	13.7 ± 0.5	13.7 ± 0.5	25–500	620 ± 20	applicable	
IPS InLine Powder Opaquer	13.2 ± 0.5	13.4 ± 0.5	13.3 ± 0.5		605 ± 20		
IPS Gingiva IPS Gingiva Modifier IPS Opaque Dentin V IPS Classic V Dentin IPS Classic V Incisal IPS Classic Opal Incisal IPS Classic V Transparent IPS Impulse Transparent IPS Impulse Occlusal Dentin IPS Impulse Molar Incisal	12.6 ± 0.5	13.4 ± 0.5	13.0 ± 0.5		585 ± 20	applicable ¹ / n.a. ¹	
IPS Effect	12.6 ± 0.5	13.5 ± 0.5	13.1 ± 0.5		595 ± 20	n.a.	
IPS Classic Shade V IPS Classic Stains P	8.3 ± 0.5	8.2 ± 0.5	8.3 ± 0.5		575 ± 20		
IPS Classic Glazing Paste	9.4 ± 0.5	9.4 ± 0.5	9.4 ± 0.5		580 ± 20		
IPS Classic Add-On Powder	12.6 ± 0.5	13.3 ± 0.5	13.0 ± 0.5		440 ± 20		
					25–400		

¹Applicable to IPS Opaque Dentin V, IPS Classic V Dentin, IPS Classic V Incisal and IPS Classic Opal Incisal; all others: n.a.

Firing table

Notice

These firing parameters are guidance values. They are valid for the Programat® furnaces from Ivoclar Vivadent AG.

Deviations (approx. +/- 10 °C/18 °F) may occur:

- depending on the furnace generation
- when using ceramic furnaces from other manufacturers
- due to regional differences in the power supply
- if several electrical devices are operated on the same circuit

IPS Classic	Firing temperature T [°C/°F]	Stand-by temperature B [°C/°F]	Closing time S [min]	Heating rate t [°C/°F/min]	Holding time H [min]	Vacuum on V1 [°C/°F]	Vacuum off V2 [°C/°F]
1 st paste opaquer firing	980/1796	403/757	06:00	80/144	01:00	550/1022	979/1794
2 nd paste opaquer firing	970/1778	403/757	06:00	80/144	01:00	550/1022	969/1776
1 st /2 nd powder opaquer firing	960/176	403/757	04:00	100/180	02:00	450/842	959/1758
1 st Dentin / Incisal / Gingiva firing	920/1688	403/757	04:00	60/108	01:00	580/1076	919/1686
2 nd Dentin / Incisal / Gingiva firing	910/1670	403/757	04:00	60/108	01:00	580/1076	909/1668
Stains firing with IPS Shade V / IPS Stains P	900/1652	403/757	04:00	60/108	01:00	0	0
Glaze firing without glazing paste	920/1688	403/757	04:00	60/108	01:00	0	0
Glaze firing with IPS Classic glazing paste	900/1652	403/757	04:00	60/108	01:00-02:00	0	0
Add-On after glaze firing (1:1)	810/1490	403/757	04:00	60/108	01:00	420/788	809/1488
Add-On after glaze firing (pure)	690/1274	403/757	04:00	60/108	00:50	420/788	689/1274

4 Safety information

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, www.ivoclar.com and your responsible competent authority.
- The current Instructions for Use and the explanation of symbols are available on the website: www.ivoclar.com/IFU
- The Summary of Safety and Clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed>.
- Basic UDI-DI: 76152082ACERA003EU

Warnings

- Observe the Safety Data Sheet (SDS) (available at Ivoclar Vivadent AG www.ivoclar.com).
- Do not inhale ceramic dust during finishing. Use suction equipment, a face mask and protective goggles.
- If applied frequently or for a prolonged period of time, peroxide-based professional bleaching agents (carbamide peroxide; hydrogen peroxide) as well as acidulated phosphate fluorides, which are used for caries prevention, can make the surface of existing restorations rough and matte.

Disposal information

Remaining stocks or removed restorations must be disposed of according to the corresponding national legal requirements.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks.

The following known residual clinical risks exist:

- Chipping / fracture / debonding of the restorative material may lead to accidental swallowing or inhalation and dental retreatment

5 Shelf life and storage

- Storage temperature for ceramic pastes and mixing liquids: 2–28 °C (36–82 °F).
- Storage temperature for ceramic powders: No special storage conditions are required.
- Expiry date: see note on bottles and packaging.
- Do not use the product after the indicated expiry date.
- Store in a vibration-free environment.
- Keep dry.
- Keep out of sunlight.
- Check to make sure that the packaging and the product are intact and undamaged before use. If in doubt, contact Ivoclar Vivadent AG or your local sales partner.

6 Additional information

Keep material out of the reach of children!

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of use. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions.

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