EIGENSCHAFTEN

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PROPERTIES

InovaPrint sep
ALL purpose 3D print separator

- ▶ For orthodontics & ZT
- ▶ Extremely clean separation
- ▶ Ingenious fit
- ▶ Minimal application – minimal consumption
- ▶ Dry in a flash
- ▶ Washable without residue

Scope of delivery

- ▶ Bottle of InovaPrint sep – 100ml
- ▶ Special brush
- ▶ Instruction manual

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USAGE

InovaPrint sep is an insulation that insulates the 3D printed model optimally against dental resins (autopolymer and composite) with a wafer-thin layer. Whether in orthodontics or dental technology – InovaPrint sep is used wherever 3D printed models are used. Perfect results on 3D printed resin and filament models.

Application

- ▶ „Undercutting“ areas of the model are blocked out as usual with wax or already in the CAD program.
- ▶ Shake the material briefly, remove the screw cap and apply with pinpoint accuracy using the fine application tip.
- ▶ Use the brush supplied to apply the insulation in exactly the right place.
- ▶ Only apply InovaPrint Sep once to the model surface – that’s enough.
- ▶ Avoid puddling with the brush.
- ▶ The light coloring provides visual support and you retain control.
- ▶ After only 5 – 10 minutes drying time, you can get started.
- ▶ Insulation residues can be rinsed off in no time at all under warm running water or simply briefly in the ultra-sonic device.

Important: Close the bottle when not in use!!!

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