


# Technical Guide

IPS e.max<sup>®</sup>  
Ceram Art

Each smile  
a piece of art





# IPS e.max<sup>®</sup> Ceram Art

IPS e.max Ceram Art is an innovative line of universal ready-to-use stain, glaze and structure pastes for the individualized design of ceramic restorations. Staining, characterizing and structuring dental and gingival surfaces can be performed at the highest esthetic level – quickly, efficiently and with the proven quality of IPS e.max. The materials are compatible with zirconium oxide, lithium disilicate, glass-ceramics and veneering ceramics. They can be used safely and flexibly on a wide range of surfaces.

IPS e.max Ceram Art has been developed for professional dental design, combining the highest esthetic standards with economic efficiency – especially in the processing of monolithic restorations.

All-ceramics	IPS e.max Ceram Art	Metal-ceramics
Zirconium oxide IPS e.max ZirCAD	Stains and glazes IPS e.max Ceram Art <ul style="list-style-type: none"><li>• Magic Glaze</li><li>• Shade</li><li>• Essence</li><li>• Illusion</li><li>• Illusion Gingiva</li><li>• Structure</li><li>• Structure Gingiva</li><li>• Universal Liquid</li></ul>	Metal-supported veneering ceramics
Lithium disilicate IPS e.max CAD IPS e.max Press		IPS Style Ceram IPS InLine
Leucite-reinforced glass-ceramics IPS Empress CAD		
Metal-free veneering ceramics IPS e.max Ceram		






Please observe the regulatory Instructions for Use ([www.ivoclar.com/eIFU](http://www.ivoclar.com/eIFU)).



# Table of Contents

<b>Explanation of symbols</b>	<b>5</b>
<b>Product information</b>	<b>5</b>
<b>IPS e.max Ceram Art assortment</b>	<b>6</b>
<b>Application</b>	<b>16</b>
Preparing the ceramic surface	16
Artistry with brushes and instruments	18
<b>Characterization, glazing and structuring</b>	<b>20</b>
<b>White / red esthetics for monolithic restorations</b>	<b>22</b>
Wet-on-wet staining / one-shot	
Conventional, individual staining technique	
Structuring technique (3D effect)	
<b>Monolithic restorations characterized with IPS e.max Ceram Art</b>	<b>32</b>
Individual staining technique and structuring technique	
<b>White / red esthetics for veneering materials</b>	<b>34</b>
Layering technique	36
<b>Art gallery</b>	<b>38</b>
<b>Fundamentals of gingival anatomy</b>	<b>42</b>
<b>General information</b>	<b>44</b>
<b>Firing parameters</b>	<b>46</b>

# Explanation of symbols

Symbol	Note
	Information
	Note on firing
	See Instructions for Use

## Product information

### Product description

IPS e.max Ceram Art pastes are optimally coordinated with the layering, press and CAD ceramics<sup>1</sup> from Ivoclar. They are suitable for use on ceramic materials in the CTE range of  $9.4\text{--}17.5 \times 10^{-6}/\text{K}$  (25–500 °C).

The ready-to-use IPS e.max Ceram Art pastes can be applied to ceramic surfaces undiluted, diluted or custom-mixed. They allow for time-efficient and highly esthetic customization of both monolithic and layered ceramic restorations through various application techniques:

#### **Wet-on-wet technique / one-shot:**

simultaneous glazing, staining and characterization in a single firing cycle.

= efficient for streamlined and time-saving procedures

#### **Conventional individual staining technique:**

more intense, individual staining and characterization with a separate Stain and Glaze firing cycle.

= intense, individual and creative

#### **Structuring technique (3D effect):**

applying structure paste to monolithic restorations to create a three-dimensional surface texture.

= effective for creating a natural depth effect with refined esthetics

#### **Staining, glazing and structuring of tooth- and gingiva-coloured ceramics**

- IPS e.max® ZirCAD (sintered)
- IPS e.max® CAD (crystallized)
- IPS e.max® Press
- IPS Empress® CAD<sup>2</sup>
- IPS e.max® Ceram
- IPS Style® Ceram
- IPS InLine®

<sup>1</sup> Except for partially crystallized IPS e.max CAD restorations.

<sup>2</sup> IPS Empress CAD restorations must not be veneered with IPS e.max Ceram Art Structure paste.


# IPS e.max Ceram Art assortment

## White esthetics

IPS e.max Ceram Art Magic Glaze FLUO, 3 g / 9 g

IPS e.max Ceram Art Universal Liquid, 15 ml / 60 ml

IPS e.max Ceram Art Shade, 3 g



IPS e.max Ceram Art Essence, 3 g



IPS e.max Ceram Art Illusion, 5 g



IPS e.max Ceram Art Structure, 5 g



## Red esthetics

IPS e.max Ceram Art Magic Glaze, 3 g / 9 g

IPS e.max Ceram Art Illusion Gingiva Base, 5 g



rose



pink



orange-pink

IPS e.max Ceram Art Illusion Gingiva, 5 g



soft-mauve



intense-mauve



soft-bone



intense-bone



light-margin



dark-margin



capillary



fibers



intense-red

IPS e.max Ceram Art Structure Gingiva, 5 g



soft-pink



pink

## Glaze pastes

The two new IPS e.max Ceram Art Magic Glaze pastes, available with or without fluorescence, are used for glazing, as carrier medium for staining, and for wetting fully ground ceramic surfaces or those that are air-abraded or partially pre-polished with a rubber polisher. They can also be used to reduce surface tension. They do not produce a beading effect.



### IPS e.max Ceram Art Magic Glaze FLUO

IPS e.max Ceram Art Magic Glaze FLUO paste is designed for use on monolithic tooth-coloured restorations made of zirconium oxide ( $ZrO_2$ ) and lithium disilicate ( $LS_2$ ).

- Undiluted glaze paste for conventional glazing and the one-shot technique.
- Diluted glaze paste for the conventional technique.



### IPS e.max Ceram Art Magic Glaze

IPS e.max Ceram Art Magic Glaze paste is designed for use on gingival areas made of zirconium oxide ( $ZrO_2$ ) or lithium disilicate ( $LS_2$ ) and on ceramic veneering materials.

- Glaze paste for the layering technique with veneering materials.
- Adjusting the colour intensity of the gingiva pastes.



### IPS e.max Ceram Art Universal Liquid

IPS e.max Ceram Art Universal Liquid is a fluid for diluting the following pastes: IPS e.max Ceram Art Magic Glaze, Shade, Essence, Illusion and Illusion Gingiva.

---

**Note:** IPS e.max Ceram Structure pastes must not be diluted with IPS e.max Ceram Art Universal Liquid.

## IPS e.max Ceram Art pastes

IPS e.max Ceram Art pastes can be integrated into the glaze over the entire surface or selectively applied to create individualizations and subtle surface effects. The appearance before firing matches the final results, providing more control and confidence during the design process.

The ready-to-use pastes have a **form-stable** consistency that allows for easy and precise fabrication of individualized, modern restorations.



### Shade

#### IPS e.max Ceram Art Shade

IPS e.max Ceram Art Shade pastes are used for adjusting the shade of ceramic restorations to match A-D and Bleach shades.

IPS e.max Ceram Art Shade Incisal 1 and Incisal 2 are incisal materials designed to enhance depth and translucency.



#### Shade correspondence with A-D and BL shades

##### IPS e.max Ceram Art Shade

Dentin 1	Dentin 2	Dentin 3	Dentin 4	Dentin 5

##### VITA® classical shades

BL	A	B	C	D

\* no registered trademark of Ivoclar Vivadent AG.

##### IPS e.max Ceram Art Shade

Incisal 1	Incisal 2

### Essence

#### IPS e.max Ceram Art Essence

IPS e.max Ceram Art Essence pastes are highly pigmented and can be applied to the restoration for individualized surface characterization.



Essence white	Essence cream	Essence lemon
Essence sunset	Essence copper	Essence khaki
Essence mahogany	Essence anthracite	

## IPS e.max Ceram Art Illusion

### Special enamel

IPS e.max Ceram Art Illusion Opal are effect enamel pastes in a variety of nuances to create specific effects.



### IPS e.max Ceram Art Illusion Opal

IPS e.max Ceram Art Illusion Opal are specially shaded pastes designed to reproduce the dynamic light-optical properties of natural teeth in the enamel area.



Opal light-blue



Opal sky-blue



Opal night-blue



### IPS e.max Ceram Art Illusion Enamel

IPS e.max Ceram Art Illusion Enamel paste is applied to increase the value (brightness) in the incisal third and along the marginal ridges to create a butterfly effect.



Enamel pearl



### Light reflectors

IPS e.max Ceram Art Illusion Mamelon pastes are light-reflecting effect materials.



### IPS e.max Ceram Art Illusion Mamelon

IPS e.max Ceram Art Illusion Mamelon pastes are intensely coloured pastes designed for the representation of mamelon structures. They are applied to accentuate the dentin structure in the incisal third.



Mamelon reddish-orange



Mamelon light-salmon



### Light absorbers

#### IPS e.max Ceram Art Illusion Absorber

IPS e.max Ceram Art Illusion Absorber pastes are light-absorbing materials used to lower the value (brightness) in the incisal region and to create areas of light absorption in the incisal and cervical thirds as well as along the proximal surfaces.



### Special Effect Halo

#### IPS e.max Ceram Art Illusion Halo


IPS e.max Ceram Art Illusion Halo pastes are used to reproduce the halo effect, a natural optical phenomenon created by light reflection at the incisal edge of natural teeth.



### Special Effect Cervical

#### IPS e.max Ceram Art Illusion Cervical

IPS e.max Ceram Art Illusion Cervical pastes help create a natural transition between the restoration and gingiva.



### IPS e.max Ceram Art Illusion Gingiva

IPS e.max Ceram Art Illusion Gingiva pastes offer a new dimension in the reproduction of lifelike artificial gingiva.

The IPS e.max Ceram Art Illusion Gingiva pastes allow for the true-to-nature replication of various gingival phenotypes, ranging from pale pink and orange-yellow to intense purple-brown tones. Three newly designed, especially intense cadmium-free red shades further enhance the options for the individualized design of natural-looking artificial gingiva.

Due to their self-glazing properties, the IPS e.max Ceram Art Illusion Gingiva pastes do not require additional glaze firing.

#### Gingival phenotypes



pink



red-blue



red-yellow



purple

#### Shade correspondence with gingival phenotypes

Gingival phenotype	pink	red-blue	red-yellow	purple
IPS e.max Ceram Art Illusion Gingiva*	Base rose	Base pink	Base orange-pink	Base pink

\* IPS e.max Ceram Art Illusion Gingiva materials can be used for mixing purposes. Red-shaded IPS e.max Ceram Art Illusion Gingiva intense-red is not intended for mixing purposes with other shades.



### Gingiva base

#### IPS e.max Ceram Art Illusion Gingiva Base


IPS e.max Ceram Art Illusion Gingiva Base pastes are used to build up the bulk of the gingival tissue, with the shade selected to match the particular gingival phenotype.



### Gingiva effects


#### IPS e.max Ceram Art Illusion Gingiva mauve

IPS e.max Ceram Art Illusion Gingiva pastes in mauve are primarily used for effect characterization in the reproduction of dark gingival tissue (phenotype – purple/brown). In addition, they can be employed to simulate the intense purple vascular structure at the end of the alveolar mucosa.



#### IPS e.max Ceram Art Illusion Gingiva margin

IPS e.max Ceram Art Illusion Gingiva margin pastes can be used to characterize the marginal gingiva and interdental papillae.





**IPS e.max Ceram Art Illusion Gingiva bone**

IPS e.max Ceram Art Illusion Gingiva bone pastes are applied in the root region of the gingiva to simulate the appearance of the underlying bone showing through the attached gingiva.



Gingiva soft-bone



Gingiva intense-bone



**IPS e.max Ceram Art Illusion Gingiva capillary, fibers, intense-red**

IPS e.max Ceram Art Illusion Gingiva pastes in intense red are used in gingival areas requiring an increased reddish tone, such as areas with vascularization. The red stain is applied to the free gingiva and distributed along the border between the attached gingiva and the mucogingival line. In addition, the stains can be used to intensify the representation of the alveolar mucosa.



Gingiva capillary



Gingiva fibers



Gingiva intense-red

## IPS e.max Ceram Art Structure

When applied in a targeted fashion, IPS e.max Ceram Art Structure pastes intensify the optical depth and internal light refraction, especially in monolithic restorations, producing a lifelike 3D effect. Morphology and surface structure provide restorations with an especially natural appearance.



### IPS e.max Ceram Art Structure

Fluorescent IPS e.max Ceram Art Structure pastes can be easily layered over characterized and glazed restorations to create subtle three-dimensional shade effects and achieve natural macro- and microstructures as well as surface texture.



Structure neutral



Structure cloud



### IPS e.max Ceram Art Structure Gingiva

Non-fluorescent IPS e.max Ceram Art Structure Gingiva pastes can be layered over characterized and glazed gingival areas to create a subtle three-dimensional shade effect and to reproduce natural macro- and microstructures as well as texture on the surface of the gingiva.



Structure Gingiva soft-pink



Structure Gingiva pink

**i** After each use, opened products must be immediately resealed to preserve their material properties.

**i** An overview of the complete assortment can be found on pages 6/7 and at [www.ivoclar.com](http://www.ivoclar.com)

# Application

## Preparing the ceramic surface

**NOTICE!** Depending on the material, the restoration surface must be prepared in varying ways. Observe the Instructions for Use of the respective material.

**NOTICE!** Smooth or prepolish the areas for which a higher gloss is desired after glaze firing using silicone polishers.

### Monolithic restorations made of zirconium oxide (ZrO<sub>2</sub>) – IPS e.max® ZirCAD

**NOTICE!** Make changes to sintered restorations only if absolutely necessary.

1. If absolutely necessary: Prepare the restoration with suitable grinding instruments (see Instructions for Use) at low pressure and low speed.
2. Optional: Blast the restoration with aluminium oxide (Al<sub>2</sub>O<sub>3</sub>) 25–70 µm at 1 bar (15 psi) pressure or 70–110 µm at 1.5 bar (22 psi) pressure.
3. Smooth the incisal and occlusal contact points and the basal side of the bridge connectors using rubber polishers.
4. Thoroughly clean the restoration with running water or a steam jet.
5. Dry the restoration.

### Monolithic restorations made of lithium disilicate (LS<sub>2</sub>) CAD – IPS e.max® CAD

**NOTICE!** Do not blast lithium disilicate restorations with Al<sub>2</sub>O<sub>3</sub> or polishing beads.

✓ The restoration is in its precrystallized (blue) state.

1. Prepare the restoration using diamond-bonded grinding instruments at low pressure and low speed. Observe the Instructions for Use.
2. Clean the restoration with ultrasound in a water bath and / or with a steam jet.
3. Dry the restoration.
4. Crystallize the restoration.

### Monolithic restorations made of lithium disilicate (LS<sub>2</sub>) Press – IPS e.max® Press

1. Prepare the restoration using ceramic- and / or diamond-bonded grinding instruments at low pressure and low speed. Observe the Instructions for Use.
2. Blast the restoration with Al<sub>2</sub>O<sub>3</sub>, 100 µm at 1 bar (15 psi) pressure.
3. Thoroughly clean the restoration with running water or a steam jet.
4. Dry the restoration.

### Monolithic restorations made of leucite-reinforced glass-ceramic CAD – IPS Empress® CAD

1. Prepare the restoration using ceramic- and / or diamond-bonded grinding instruments at low pressure and low speed. Observe the Instructions for Use.
2. Carefully prepare the restoration margins using rubber and silicone polishers. Observe the Instructions for Use.
3. Thoroughly clean the restoration with running water or a steam jet.
4. Dry the restoration.

### Veneering ceramics – IPS e.max® Ceram, IPS Style® Ceram, IPS InLine®

1. Prepare the restoration using ceramic- and / or diamond-bonded grinding instruments at low pressure and low speed. Observe the Instructions for Use.
2. Thoroughly clean the restoration with running water or a steam jet.
3. Dry the restoration.


**NOTICE!** Select the appropriate grinding tools according to the recommendations provided in the Instructions for Use of the material.

# Preparing the ceramic surface



## **Monolithic restorations**

Standard surface finish of the restoration after processing



## **Veneered restorations**

Surface of the restoration after processing  
Individual natural macro- and microtextures with smoothing using silicone polishers

## Artistry with brushes and instruments

Standard dental instruments and brushes can be used to process the new generation of viscous IPS e.max Ceram Art pastes. The new pastes offer versatile options for finalizing restorations and adding an artistic touch. Embrace your creativity and turn each restoration into a small work of art – true 'ART'.

### Medium viscosity

Medium-viscosity IPS e.max Ceram Art Magic Glaze, Shade, Essence, Illusion and Illusion Gingiva pastes can be applied with modern synthetic brushes or brushes made with natural bristles. Medium-sized brushes are recommended for applying Magic Glaze and for incorporating stains, e.g. Shade Dentin, on larger areas. Fine-tipped brushes in the smallest size are ideal for artistic characterization.



### High viscosity

High-viscosity **IPS e.max Ceram Art Structure** and **Structure Gingiva** pastes should first be lightly spatulated on a mixing pad using a metal-free spatula to achieve a workable consistency. Very stiff brushes and metal-free instruments are ideal for applying and shaping the highly viscous structure paste over large areas on dental and gingival surfaces. Microtexture can be added to the surface using a fine brush with a firm tip or a fine instrument. To subtly replicate the natural stippling effect, or 'orange-peel' texture, on the gingival surface, a stiff, flat brush is recommended.



# Characterization, glazing and structuring

- ✓ The restoration is free of dirt and grease residues.
- ✓ The restoration demonstrates a smooth surface without sharp edges / ridges.

## White / red esthetics for monolithic restorations

IPS e.max® ZirCAD, IPS e.max® CAD, IPS e.max® Press, IPS Empress® CAD

### Wet-on-wet technique / one-shot

simultaneous glazing, staining and characterization in a single firing cycle.  
= efficient for streamlined and time-saving procedures

### IPS e.max Ceram Art pastes:

- before use, mix thoroughly in the tub using a clean, metal-free spatula.
- can also be mixed together on a mixing pad to further expand the variety of shades.

**Note:** IPS e.max Ceram Art Illusion Gingiva pastes are applied directly to the surface of the restoration (no wetting with glaze paste is required).





## White / red esthetics for monolithic restorations

### Wet-on-wet staining / one-shot



Prepared surface / restoration, see page 16.



1. Apply **undiluted** IPS e.max Ceram Art Magic Glaze FLUO on the teeth in an even layer using a brush.



2. Stain the teeth with IPS e.max Ceram Art Shade Dentin/Incisal and adjust the shade in accordance with the shade guide.
3. Characterize the teeth with IPS e.max Ceram Art Illusion and IPS e.max Ceram Art Essence.
4. Conduct Stain and Glaze firing.
5. Repeat Steps 1–4, if required.

 Firing parameters, see pages 46/47

### Final restoration



- IPS e.max Ceram Art**
- Magic Glaze FLUO
  - Shade Dentin 2
  - Shade Incisal 1
  - Illusion Mamelon light-salmon
  - Illusion Absorber fog



Prepared surface / restoration,  
see page 16.



1. Apply **undiluted** IPS e.max Ceram Art Glaze FLUO on the teeth in an even layer using a brush.

## Final restoration



2. Stain the teeth with IPS e.max Ceram Art Shade Dentin/ Incisal and adjust the shade in accordance with the shade guide.
3. Characterize the teeth with IPS e.max Ceram Art Illusion and IPS e.max Ceram Art Essence.
4. Apply undiluted IPS e.max Ceram Art Illusion Gingiva Base on the gingival area (self glazing).
5. Characterize the gingival area with IPS e.max Ceram Art Illusion Gingiva (self glazing).
6. Conduct Stain and Glaze firing.
7. Repeat Steps 1–6, if required.



- IPS e.max Ceram Art**
- Magic Glaze FLUO
  - Shade Dentin 2
  - Shade Incisal 1
  - Illusion Mamelon light-salmon
  - Illusion Absorber fog
  - Illusion Gingiva Base rose
  - Illusion Gingiva soft-bone
  - Illusion Gingiva capillary

 Firing parameters, see pages 46/47



## Conventional, individual staining technique

more intense, individual staining and characterization with a separate Stain and Glaze firing cycle.

= intense, individual and creative

### IPS e.max Ceram Art pastes:

- before use, mix thoroughly in the tub using a clean, metal-free spatula.
- can also be mixed together on a mixing pad to further expand the variety of shades.
- the viscosity of the pastes can be diluted and individually adjusted using IPS e.max Ceram Universal Liquid.

**NOTICE!** Avoid pooling and overly thick stain layers.


More intensive shades are achieved by repeating the staining procedure and firing, not by applying thicker layers.

**NOTE:** IPS e.max Ceram Art Illusion Gingiva pastes are directly applied on the surface of the restoration (no wetting with glaze paste required).



## White / red esthetics for monolithic restorations

### Conventional, individual staining technique



Prepared surface / restoration, see page 16.



1. Apply IPS e.max Ceram Art Glaze FLUO alone or slightly diluted with IPS e.max Ceram Art Universal Liquid on the teeth in a **thin** layer using a brush.



2. Stain the teeth with IPS e.max Ceram Art Shade Dentin/Incisal and adjust the shade in accordance with the shade guide.  
3. Characterize the teeth with IPS e.max Ceram Art Illusion and IPS e.max Ceram Art Essence.  
4. Conduct stain firing.  
5. Repeat Steps 1–4, if required.

### Final restoration



6. Apply undiluted IPS e.max Ceram Art Glaze FLUO on the teeth in an even layer using a brush.  
7. Conduct Stain and Glaze firing.  
8. Repeat Steps 6–7, if required.



#### IPS e.max Ceram Art

Magic Glaze FLUO

- Shade Dentin 2
- Shade Incisal 1
- Illusion Absorber fog
- Illusion Mamelon reddish-orange
- Essence sunset
- Illusion Enamel pearl
- Illusion Opal night-blue
- Illusion Halo ocher

 Firing parameters, see pages 46/47



Prepared surface / restoration, see page 16.



1. Apply IPS e.max Ceram Art Glaze FLUO alone or slightly diluted with IPS e.max Ceram Art Universal Liquid on the teeth in a thin layer using a brush.



2. Apply undiluted IPS e.max Ceram Art Illusion Gingiva Base on the gingival area (self glazing).

Optional: When applying the first base layer, selectively leave out the lighter zones at the tooth root and marginal gingiva.



3. Stain the teeth with IPS e.max Ceram Art Shade Dentin/Incisal and adjust the shade in accordance with the shade guide.
4. Characterize the teeth with IPS e.max Ceram Art Illusion and IPS e.max Ceram Art Essence.
5. Conduct stain firing.
6. Repeat Steps 1–5, if required.

## Final restoration



7. Apply undiluted IPS e.max Ceram Art Glaze FLUO on the teeth in an even layer using a brush.
8. Characterize the gingival area with IPS e.max Ceram Art Illusion Gingiva (self glazing).
9. Conduct Stain and Glaze firing.



Firing parameters, see pages 46/47

### IPS e.max Ceram Art

- Magic Glaze FLUO
- Shade Dentin 2
- Shade Incisal 1
- Illusion Absorber fog
- Illusion Mamelon reddish-orange
- Essence sunset
- Illusion Enamel pearl
- Illusion Opal night-blue
- Illusion Halo ocher
- Illusion Gingiva Base pink
- Illusion Gingiva intensive-bone
- Illusion Gingiva light-margin
- Illusion Gingiva capillary

## Structuring technique (3D effect)

applying structure paste to monolithic restorations to create a three-dimensional surface texture.  
= effective for creating a natural depth effect with refined esthetics

### IPS e.max Ceram Art Structure:

- the surface of the restoration has been characterized and glazed using the wet-on-wet / one-shot or conventional technique.
- lightly mix the high-viscosity IPS e.max Ceram Art Structure and IPS e.max Ceram Art Structure Gingiva pastes on a mixing pad using a clean, metal-free spatula to achieve a smooth viscosity.
- for better sculpting, lightly moisten the brush / instrument with IPS e.max Ceram Art Universal Liquid.
- use IPS e.max Art Magic Glaze to spread IPS e.max Ceram Art Structure paste thinly along the transitions in the proximal, incisal and tapered areas of the restoration.

**IMPORTANT!** IPS e.max Ceram Art Structure & Structure Gingiva pastes must not be diluted with IPS e.max Ceram Art Universal Liquid.

**NOTICE!** After the structure firing of IPS e.max Ceram Art Structure, check the contact area to the antagonist and the contact points.




### Layer thicknesses of IPS e.max Ceram Art Structure and IPS e.max Ceram Art Structure Gingiva

A maximum layer thickness of approx. 0.8 mm is recommended for the application of IPS e.max Ceram Art Structure pastes. After the firing process, this results in an a maximum layer thickness of approx. 0.5 mm, creating a 3D effect.

### Adjusting the gloss

The restoration's surface gloss can be adjusted with additional mechanical polishing, if desired. By creating shiny and matte areas, the surface can be made to look even more natural. This works especially well on restorations fired with IPS e.max Ceram Art Structure or IPS e.max Ceram Art Structure Gingiva. With the polishing instruments of the OptraGloss® Extra Oral range, mechanical pre- and high-gloss polishing can be achieved quickly and easily.



OptraGloss® Extra Oral from Ivoclar comprises universal polishers for the high-gloss polishing of ceramic restorations in two steps and a polishing paste for extraoral use.



## White / red esthetics for monolithic restorations

### Structuring technique (3D effect)



Application in the wet-on-wet technique or conventional technique has been performed, see pages 20–22.



1. Apply undiluted IPS e.max Ceram Art Structure to the glazed surface of the restoration (self glazing).
2. Sculpt the morphology of the restoration.
3. Incorporate the microtexture in the surface using instruments and brushes.
4. Conduct structure firing.



Firing parameters, see pages 46/47

### Final restoration



IPS e.max Ceram Art  
● Structure neutral  
● Structure cloud



Application in the wet-on-wet technique or conventional technique has been performed, see pages 20–27.

1. Apply undiluted IPS e.max Ceram Art Structure and Structure Gingiva to the glazed surface of the restoration (self glazing).
2. Sculpt the morphology of the restoration.
3. Incorporate the microtexture in the surface using instruments and brushes.
4. Conduct structure firing.



Firing parameters, see pages 46/47

## Final restoration



### IPS e.max Ceram Art

- Structure neutral
- Structure cloud
- Structure Gingiva soft-pink
- Structure Gingiva pink

**Monolithic restorations characterized with IPS e.max Ceram Art**

Individual staining technique and structuring technique



Phenotype: Gingiva pink

Phenotype: Gingiva red-blue





Phenotype: Gingiva red-yellow

Phenotype: Gingiva purple



## White / red esthetics for veneering materials

IPS e.max® Ceram, IPS Style® Ceram, IPS InLine®

### IPS e.max Ceram Art pastes:

- the restoration has been partially or completely veneered with all-ceramic or metal-ceramic veneering materials.
- before use, mix thoroughly in the tub using a clean, metal-free spatula.
- can also be mixed together on a mixing pad to further expand the variety of shades.

**TIP!** To achieve a higher gloss effect with diluted IPS e.max Ceram Magic Glaze, use a firing temperature 10–20 °C (50–68 °F) below the last firing of the veneering ceramic.





## White / red esthetics for veneering materials

### Layering technique



Prepared surface / restoration,  
see page 16



1. Apply IPS e.max Ceram Art Glaze in a **thin**  
even layer to the teeth using a brush.



2. Stain the teeth with IPS e.max Ceram Art Shade Dentin/Incisal and adjust the shade in accordance with the shade guide.
3. Characterize the teeth with IPS e.max Ceram Art Illusion and IPS e.max Ceram Art Essence.
4. Conduct Stain and Glaze firing.
5. Repeat Steps 1–4, if required.



Firing parameters, see pages 46/47

### Final restoration



- IPS e.max Ceram Art
- Magic Glaze
  - Shade Dentin 2
  - Essence creme
  - Essence white



Prepared surface / restoration,  
see page 16



1. Apply IPS e.max Ceram Art Glaze in a **thin**  
even layer to the teeth using a brush.

## Final restoration



2. Stain the teeth with IPS e.max Ceram Art Shade Dentin/Incisal and adjust the shade in accordance with the shade guide.
3. Characterize the teeth with IPS e.max Ceram Art Illusion and IPS e.max Ceram Art Essence.
4. Apply undiluted IPS e.max Ceram Art Illusion Gingiva Base to the gingival area (self glazing).
5. Characterize the gingival area with IPS e.max Ceram Art Illusion Gingiva and IPS e.max Ceram Art Essence.
6. Conduct Stain and Glaze firing.
7. Repeat Steps 1–6, if required.



- IPS e.max Ceram Art**
- Magic Glaze
  - Shade Dentin 2
  - Essence creme
  - Essence white
  - Illusion Gingiva fibers
  - Illusion Gingiva capillary


# Art gallery

Artistic characterization using IPS e.max Ceram Art



White esthetics – anterior

Shade Incisal 1



- Essence sunset
- Essence copper
- Essence mahogany
- Essence white

Shade Incisal 1

Shade Dentin 2



Essence kaki


White esthetics – posterior

# Art gallery


Artistic characterization using IPS e.max Ceram Art



White / red esthetics – anterior and gingiva



Full surface application with  
IPS e.max Ceram Art Structure Gingiva pink



Illusion Gingiva  
intense-red



Illusion Gingiva  
light-margin

Illusion Gingiva  
dark-margin

# Fundamentals of gingival anatomy

## Gingival morphology

Replicating gingiva with lifelike esthetics requires a fundamental understanding of the natural gingival morphology.



1. In the cervical area, you will find the marginal gingiva (=sulcus gingivalis) and the free gingiva, which contain microscopic spaces. The tissue in this region is flexible and is partially compressed or stretched during procedures such as implant placement, for example, when performing a sinus lift.
2. The attached gingiva is a firm tissue anchored to the alveolar bone, creating a connection to the teeth.
3. Below this, the gingiva extends along a wavy line, the mucogingival line. This is followed by the alveolar mucosa, the movable oral mucosa.

The areas of the marginal gingiva (sulcus gingiva) and the free gingiva, which surround the tooth like a band, appear slightly curved and light reddish in colour. The slightly curved, rounded ridge between the teeth, known as the interdental papilla, may exhibit varying degrees of recession depending on the patient's gingival condition. It is essential to adequately recreate the interdental papillae and preserve a healthy shape to prevent food particles from getting trapped.



The underlying attached gingiva appears as if the bone is pushing it outward. The tissue appears slightly tense and curved and it looks somewhat whitish. The alveolar mucosa extends below the gingival line (mucogingival line) of the attached mucosa. Its movable tissue has a slightly increased colour saturation in a reddish tone. Due to the thin tissue, fine blood vessels or violet colourations may also become visible.



Pictures courtesy of DT Ilhwan Jang

# General information

## Frequently asked questions

### Can I use IPS e.max Ceram Art for staining ceramic materials from other manufacturers?

**Yes.** IPS e.max Ceram Art can be used for all ceramic materials that fall within the specified CTE range. IPS e.max Ceram Art pastes are compatible with ceramics in the CTE range of  $9.4\text{--}17.5 \times 10^{-6}/\text{K}$  (25 - 500 °C).  
ⓘ Important: The use of IPS e.max Ceram Art on ceramic materials from other manufacturers is the sole responsibility of the user.

### Can I use IPS e.max Ceram Art for the crystallization firing of IPS e.max CAD restorations?

**No.** IPS e.max Ceram Art is not suitable for staining and glazing IPS e.max CAD restorations in the 'blue' state. After crystallization firing, IPS e.max Ceram Art can be used for characterization and glazing.

### Is it possible to mix the IPS e.max Ceram Art Shade, Essence, Illusion and Illusion Gingiva pastes with each other to further expand the variety of shades?

**Yes.** The IPS e.max Ceram Art pastes can be mixed with each other.  
ⓘ Important: Red-shaded IPS e.max Ceram Art Illusion Gingiva intense-red is not intended for mixing with gingiva shades.

### Can I dilute IPS e.max Ceram Art pastes with IPS e.max Ceram Art Universal Liquid?

**Yes.**  
ⓘ Important: IPS e.max Ceram Art Structure and Structure Gingiva pastes must not be diluted with IPS e.max Ceram Art Universal Liquid.

### Can I fire IPS e.max Ceram Art in furnaces from other manufacturers?

**Yes.** When using furnaces other than the Programat® furnaces from Ivoclar Vivadent AG, the specified firing parameters should be adopted. However, these firing parameters should be considered as guidelines only, as every furnace type may produce different firing results. We recommend carrying out a test firing first.

### Can I use the IPS Ivocolor Shade Dentin pastes, IPS Ivocolor Shade Incisal pastes and IPS Ivocolor Essence powders together with IPS e.max Ceram Art?

**Yes.** IPS Ivocolor Shade pastes and IPS Ivocolor Essence powders can be used together with the IPS e.max Ceram Art pastes.

## Can I carry out multiple Stain and Glaze firing cycles with the IPS e.max Ceram Art pastes?

**Yes.** When using the conventional staining technique in conjunction with IPS e.max Ceram Art Shade, Essence, Illusion and Illusion Gingiva, stain firing can be carried out multiple times.

**Yes.** For the one-shot technique using IPS e.max Ceram Art Magic Glaze FLUO, a second Stain and Glaze firing can be performed as an add-on.

## Can I carry out a corrective (add-on) firing cycle after Stain-Glaze firing?

**Yes.** A corrective firing can be performed with the following add-on materials: IPS e.max Ceram Add-On 690 °C, IPS Style Ceram Add-On 690 °C and IPS InLine Add-On 690 °C.

⊖ Important: The firing temperature of the add-on materials must not exceed the firing temperature of the IPS e.max Ceram Art pastes.

## As an alternative, can I perform a corrective (add-on) firing with the IPS e.max Ceram Art Structure pastes after the Stain-Glaze firing?

**Yes.** It is possible to perform a corrective firing with the IPS Art Structure pastes (structure firing).

⊖ Important: Only one corrective firing may be performed with the IPS e.max Ceram Art Structure pastes.

## Is a final glaze firing required after the structure firing and should it be performed at all?

**No.** An additional glaze firing cycle is not required, because the IPS e.max Ceram Art Structure pastes are self-glazing.

⊖ Note: If a higher level of gloss is required after the structure firing, high-gloss polishing can easily be achieved using OptraGloss® Extra Oral.

## Can I use rotary instruments to rework a fired IPS e.max Ceram Art Structure surface?

**Yes.** Reworking of the fired IPS e.max Ceram Art Structure surface is possible.

⊖ Note: Subsequently, the fired surface can be pre-polished and high-gloss polished using OptraGloss Extra Oral or, alternatively, a glaze firing can be performed without glaze paste at 20 °C (68 °F) below the structure firing temperature, with a reduced holding time of max. 1 minute.

# Firing parameters

## Stain and Glaze firing when using the staining technique (monolithic)

IPS e.max Ceram Art Magic Glaze FLUO, Shade, Essence, Illusion, Illusion Gingiva

IPS e.max Ceram Art	Stand-by temperature B [°C/°F]	Closing time S* [min]	Heating rate t [°C/°F/min]	Firing temperature T [°C/°F]	Holding time H [min]	Vacuum on V1 [°C/°F]	Vacuum off V2 [°C/°F]	Long-term cooling L [°C/°F]
IPS e.max® ZirCAD 1-4 units	403/757	IRT/6	45/81	710/1310	2:00	-	-	-
IPS e.max® ZirCAD 5-7 units	403/757	IRT/6	40/72	720/1328	2:00	-	-	**
IPS e.max® ZirCAD 8-14 units	403/757	IRT/6	30/54	730/1346	2:00	-	-	**
IPS e.max® CAD	403/757	IRT/6	45/81	710/1310	2:00	-	-	***
IPS e.max® Press	403/757	IRT/6	45/81	710/1310	2:00	-	-	-
IPS Empress® CAD	403/757	IRT/6	45/81	710/1310	2:00	-	-	-

\* IRT standard mode

\*\* In case of significant differences in the cross-sections of individual units within a restoration, long-term cooling L down to 500 °C/932 °F is recommended for each firing cycle.

\*\*\* For layer thicknesses exceeding 2 mm, long-term cooling L down to 500 °C/932 °F is required.

## Structure firing in the structuring technique (monolithic)

IPS e.max Ceram Art Structure / IPS e.max Ceram Art Structure Gingiva

IPS e.max Ceram Art	Stand-by temperature B [°C/°F]	Closing time S* [min]	Heating rate t [°C/°F/min]	Firing temperature T [°C/°F]	Holding time H [min]	Vacuum on V1 [°C/°F]	Vacuum off V2 [°C/°F]	Long-term cooling L [°C/°F]
IPS e.max® ZirCAD 1 unit	403/757	IRT/6	45/81	710/1310	2:00	550/1022	710/1310	-
IPS e.max® ZirCAD 2-4 units	403/757	IRT/6	40/72	710/1310	2:00	550/1022	710/1310	-
IPS e.max® ZirCAD 5-7 units	403/757	IRT/6	45/81	720/1328	2:00	550/1022	720/1328	500/932
IPS e.max® ZirCAD 8-14 units	403/757	IRT/6	30/54	720/1328	2:00	550/1022	720/1328	500/932
IPS e.max® CAD	403/757	IRT/6	45/81	710/1310	2:00	550/1022	710/1310	***
IPS e.max® Press	403/757	IRT/6	45/81	710/1310	2:00	550/1022	710/1310	-

\* IRT standard mode

\*\*\* For layer thicknesses exceeding 2 mm, long-term cooling L down to 500 °C/932 °F is required.

## Stain and Glaze firing in the layering technique (full veneering / partial veneering)

IPS e.max Ceram Art Magic Glaze Shade, Essence, Illusion, Illusion Gingiva

IPS e.max Ceram Art	Stand-by temperature B [°C/°F]	Closing time S* [min]	Heating rate t [°C/°F/min]	Firing temperature T [°C/°F]	Holding time H [min]	Vacuum on V1 [°C/°F]	Vacuum off V2 [°C/°F]	Long-term cooling L [°C/°F]
IPS e.max® Ceram / IPS e.max® ZirCAD 1-4 units	403/757	IRT/6	45/81	710/1310	1:00	-	-	450/842
IPS e.max® Ceram / IPS e.max® ZirCAD 5-7 units	403/757	IRT/6	40/72	720/1328	1:00	-	-	450/842
IPS e.max® Ceram / IPS e.max® ZirCAD 8-14 units	403/757	IRT/6	30/54	720/1328	1:00	-	-	450/842
IPS e.max® Ceram / IPS e.max® CAD	403/757	IRT/6	60/108	710/1310	1:00	-	-	-
IPS e.max® Ceram / IPS e.max® Press	403/757	IRT/6	60/108	710/1310	1:00	-	-	-
IPS Style® Ceram	403/757	IRT/6	60/108	750/1382	1:00	-	-	-
IPS InLine®	403/757	IRT/6	60/108	830/1526	1:00	-	-	-

\* IRT standard mode

These firing parameters are guidance values and they are valid for the Programat® furnaces from Ivoclar Vivadent AG.

Deviations (approx. ± 10 °C/50 °F) may occur:

- Depending on the furnace generation.
- Ceramic furnaces of other manufacturers.
- Regional differences in the supply voltage.
- Operating several electrical devices in the same circuit.



**Manufacturer**

Ivoclar Vivadent AG  
Bendererstrasse 2  
9494 Schaan/Liechtenstein  
[www.ivoclar.com](http://www.ivoclar.com)

Date information prepared: 2025-07 / Rev.0  
768998 / EN