Cercon®ht, Cercon®xt

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Instructions for use – multilingual

THE DENTAL
SOLUTIONS
COMPANY™

 Dentsply
Sirona

Cercon®ht, Cercon®xt

CZ	Návod k použití Cercon® ht, Cercon® xt	對	11
DA	Brugsanvisning for Cercon® ht, Cercon® xt	對	21
DE	Gebrauchsanweisung Cercon®ht, Cercon® xt		31
EN	Instructions for Use Cercon®ht, Cercon® xt	對	41
ES	Instrucciones de uso Cercon®ht, Cercon® xt		51
FR	Notice d'utilisation Cercon®ht, Cercon® xt		61
GR	Cercon® ht, Cercon® xt, Οδηγίες Χρήσης		71
HR	Cercon® ht, Cercon® xt Upute za upotrebu	對	81
HU	Cercon® ht, Cercon® xt használati útmutató		91
IT	Istruzioni per l'uso Cercon®ht, Cercon® xt		101
JA	Cercon ht, Cercon xt 取扱説明書	對	111

LT	„Cercon®ht“, „Cercon®xt“ naudojimo instrukcija	121
LV	Cercon® ht, Cercon® xt lietošanas instrukcija	131
NL	Gebruiksaanwijzing Cercon® ht, Cercon® xt	141
NO	Bruksanvisning for Cercon® ht, Cercon® xt	151
PL	Instrukcja postępowania - Cercon® ht, Cercon® xt	161
PT	Cercon® ht, -Cercon® xt Instruções de uso	對 171
RO	Instrucțiuni de utilizare Cercon®ht, Cercon® xt	181
RU	Инструкция по применению Cercon® ht, Cercon® xt	191
SK	Cercon® ht, Cercon® xt – Návod na používanie	201
SL	Navodila za uporabo Cercon® ht, Cercon® xt	211
SV	Cercon® ht, Cercon xt® användningsinstruktioner	221
TR	Cercon® ht, Cercon® xt için kullanım talimatları	231
ZH	泽康高透和超透氧化锆的使用说明	241....

Instructions for use Cercon®ht, Cercon®xt

EN

Product description:

-Cercon® ht and -Cercon® xt blanks are made of yttrium oxide- (yttria-) stabilized zirconium oxide (zirconia) (Y-TZP).

They are used in fabricating frameworks for fixed prosthetic restorations.

Depending on the framework design, -Cercon® ht and Cercon® xt frameworks can be ceramically veneered or deli-vered as fully contoured restorations.- Which blank is selected will depend on the tooth shade to be -reproduced and the space available for the veneer. With fully contoured restorations, no space is required for the ceramic -veneer, which may allow the dentist to preserve more of the tooth substance during preparation.-

Framework material	Zirconium oxide (Y-TZP)
Temporary cementation	Possible for Cercon®ht Not possible for Cercon®- xt
Definitive cementation	Adhesive cementing Conventional cementing

Objects are individually fabricated to your digital design specifications such as anatomic contour, framework and wall thickness, connector diameter and cementing gap.

Technical specifications Cercon®ht:

- Type II, class 5
(pursuant to DIN EN ISO 6872:2015)
- CTE: 10.5 µm/m·K (25–500 °C)
- Modulus of elasticity: 210 GPa
- Flexural strength: approx. 1 200 MPa
(three-point flexural testing)

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Cercon®ht:

Cer-con® ht is indicated in the anterior and posterior segments for:

- Crowns
- Telescopic primary crowns
- Multi-unit bridges (with no more than two pontics between- abutment crowns; with no more than 6 units*)
- Two piece abutments**

Technical specifications Cercon®xt:

- Type II, class 4
(pursuant to DIN EN ISO 6872:2015)
- CTE: 10.1 µm/m·K (25–500 °C)
- Modulus of elasticity: 210 GPa
- Flexural Strength: approx. 750 MPa
(three-point flexural testing)

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Cercon®xt:

Cer-con® xt is indicated for all ceramic restorations for anterior and posterior locations:

- Crowns
- 3-unit bridges (up to the 2nd premolar)

Cer-con® ht can be used as a substructure (framework) which is then veneered with a dental veneering ceramic or can be used for full-contour application (without veneering) as well. In the case of telescopic primary crowns, the substructure is not veneered.

*for Canada only ** not valid for USA

Contraindications:

This medical product may not be used in patient hypersensitivity to zirconia (Y-TZP) or one - of the other ingredients
Bruxism or recalcitrant parafunctional habits (for ceramically veneered frameworks)
Insufficient available space
Endodontic posts
Endosseous implants
Inlay bridges
3-unit bridges in the molar region (for Cercon® xt only)

Warnings:

Possible cross-reactions or interactions of this product with other products or material- -already present in the oral environment must be taken into consideration by the dentist when -selecting this product.

Precautions:

Please note:

Keep product dust away from eyes.
Avoid any contact with mucosa.
After use, wash your hands and apply a hand cream.
Do not smoke, eat or drink while handling the product.
Do not swallow the product.
Do not inhale dust particles during grinding.
Use local vacuum suction and suitable mouth/face protection during manual machining at the workplace.

The safety and warning notes listed here describe how to use our product in a safe and risk-free manner-.

Notify- the dentist in charge of all factors described above if you use this product for a custom design and make sure to comply with the pertinent Material Safety Data Sheets (MSDS).

Adverse effects:

If properly processed and used, adverse effects of this product are highly unlikely. However, -reactions of the immune system (such as allergies) to substances contained in the material or localized paraesthesia (such as taste disturbances or irritation of the oral mucosa) cannot be completely ruled out as a matter of principle. Should you hear or be informed of any adverse effects – even when doubtful – we would like to -request notification.

Special notes on framework design

Framework design for fully contoured restorations:

Indicated especially where limited occlusal space is available, for frameworks that are not to be veneer-ed or for frameworks with body stains.

The surface of fully contoured frameworks can carefully be optimized with rotary tools such as fine -cutters prior- to sintering.

Make sure that the occlusal surface is not changed by subsequent deepening of fissures, as notching action- may compromise the strength of the material. Please note that flat occlusal reliefs can extend the live expec- tancy of fully contoured restorations. When making manual adjustments, make sure nev-er to separate the interdental- spaces of the frameworks with cutting discs or other rotary -instruments, Doing so might damage the framework and compromise the strength of the material!

Important note:

Please ensure that the minimum wall thickness of the framework in the area of the occlusal surface is -respected even after occlusal adjustments.

Framework design for veneered restorations:

Frameworks to be ceramically veneered are designed to reduced anatomical contour to provide maximum- support for the veneer. The frameworks can be veneered using the press-on technique or the build-up technique.

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Wall and border thickness:	Cercon®ht	Cercon®xt
Wall thickness, single crowns	0.4 mm	0.7 mm
Margin thickness, single crowns	0.2 mm	—
Wall thickness, bridges	0.5 mm	0.7 mm
Margin thickness, bridges	0.2 mm	—

Additional dimensional requirements for the anterior region:	Cercon®ht	Cercon®xt
Number of pontics	2	1
Connector cross-section	6 mm ₂	12 mm ₂

Additional dimensional requirements for the posterior region:	Cercon®ht	Cercon®xt
Number of pontics	2	1
Connector cross-section	9 mm ₂	16 mm ₂
Cantilever pontic at tooth position (only one pontic, up to one premolar width)	Up to the second premolar	—
Connector cross-section for this cantilever- pontic	12 mm ₂	—

Nesting

Special nesting notes for large-span bridges \9 units or more, only for Cercon® ht\

For distortion-free sintering of long-span bridges (9 units or more) when nesting the object in a Cer-con® ht 98 disk, make sure to add a sintering reinforcement in the form of a "tongue" to create.

Always place objects (single copings, three-unit frameworks) within the sintering reinforcement ("tongue")- in order to achieve a homogeneous shrinkage of large-span objects during sintering.

Entering the shrinkage factor for the milling unit:

In CAM software allowing the entry of 3 dimensions, enter X, Y and Z values.

In CAM software allowing the entry of 2 dimensions, enter X or Y and Z values.

In CAM software allowing the entry of 1 dimension, enter the X value.

Finishing

Special notes on finishing:

Please read the respective Instructions for Use for your device regarding further processing.

We recommend the following milling strategies for secure processing:

Procedure	Tool, spherical	Feed rateZ mm/min	Feed rateF mm/min	Speed rpm	a ₁ mm	a ₂ mm	Strategy	Dimension
Contour rough-machining, occlusal- side (OS)	HM Ø2	800	1800	22 000	0.8	1	Contour-parallel rough-machining	0.3
Contour roughmachining, cavity- side (CS)	HM Ø2	800	1800	22 000	0.8	1	Contour-parallel rough-machining	0.3
Pre-finishing (CS)	HM Ø2	400	1200	22 000		0.2	Residue material rough-machining	0.1
Pre-finishing (CS)	HM Ø2	400	1200	22 000		0.2	Residue material rough-machining	0.1
Finishing (OS)	HM Ø1	400	1200	27 000		0.15	3D offset	0
Finishing without cavity (CS)	HM Ø1	400	1200	27 000		0.15	3D offset from inside	0
Finishing cavities	HM Ø1	250	1000	27 000		0.1	3D offset	0
Fissures	HM Ø 0.5 - 0.6	250	1000	35 000		0.1	Partial finishing, 3D	0

The milling strategies are recommendations. If necessary, carry out trial millings and adapt the parameters as needed.

Separating

Notes on separating the objects:

Separate the objects from the disk by sandblasting with aluminium oxide (50 μm , max. 1.5 bar). The sandblasting aid facilitates the removal of objects from the blank. It helps prevent framework fracture or other damage to the objects. For large-span bridges (9 units or more), separate only the labial and buccal sprues of the objects and the "tongue" connector, because the objects must be sintered together- with that "tongue". Any protruding ridges at the bottom of the "tongue" must be removed in order to assure that the objects will stay firmly put on the sintering tray. The smaller objects nested inside the "tongue" are detached completely and sintered separately.

Sandblasting aid

Milling and sandblasting the object

Sintering

Sintering in the Cercon®- heat plus P8:

1 500 °C in the -Cercon® heat plus P8

– Program #4 for bridges up to 8 units (Cercon® ht), up to 3 units (Cercon® xt), $T_{\text{max}} = 1\,500\,^{\circ}\text{C}$

– Program #5 for bridges of 9 units (Cercon® ht) or more, $T_{\text{max}} = 1\,500\,^{\circ}\text{C}$

Special sintering notes for large-span bridges (9 units or more, Cercon®- ht)

Two wide-span bridges (9 or more units) can be sintered in the Cercon® heat plus P8 at one time. Place the objects- on the sintering block with due regard to the internal vertical clearance of the Cercon® heat plus P8 (130 mm) and the need to facilitate contraction without any mechanical constraints. The objects must not touch the sintering block during the sintering.

Sintertray with sintering block

Correct position on the sintering tray

Sintering in the heat DUO/Multimat2Sinter:

- Program #6: Speed program for bridges up to 6 units (Cer-con® ht), for bridges up to 3 units (Cer-con® xt), $T_{\max} = 1\,540\text{ °C}$
- Program #7: Standard sintering for bridge frameworks up to 8 units (Cer-con® ht); for bridges up to 3 units (Cer-con® xt), $T_{\max} = 1\,520\text{ °C}$
- Program #8: Sintering program for bridges of 9 units (Cer-con® ht) or more, $T_{\max} = 1\,520\text{ °C}$

Special sintering notes for large-span bridges 9 units or more Cer-con® ht

Two wide-span bridges (9 or more units) can be sintered in the heat DUO or Multimat2Sinter at one time. Place the objects on the special sintering bar for large-span bridges and keep in mind that the objects must not touch the sintering block during the sintering.

Note the maximum dimensions for milling:

Height: 65 mm

Width: 90 mm

Placement of bridges with sintering bar (9 or more units)

Sintering in third-party furnaces:

Sintering results may be adversely affected by e.g.:

- Incorrect sintering temperatures
- Insufficient heating power
- Incorrect temperature curves
- Incorrect object placement
- Insufficient heat-storing capacity of the furnace over the time of the sintering cycle
- Manufacturer-related or age-related variations in heating performance
- Object contamination by oxidation products emitted by non-enclosed heating elements

Any of these factors by itself or in combination may reduce the maximum strength of our above-mentioned zirconium dioxide materials and compromise the life expectancy of the frameworks.

For these reasons, we cannot grant general approval for the use of third-party furnaces for sintering two-piece abutments (mesostructures) and crown and bridge frameworks made of Cer-con® ht and Cer-con® xt. We will, however, open the system technically for the use of furnaces of third-party furnaces only on condition that the following requirement are met:

IMPORTANT!

The sintering of Cercon®ht and Cercon®xt in third-party sintering furnaces is done at the user's own responsibility and risk. Dentsply Sirona can accept no responsibility for any kind of damage caused by the sintering of zirconia- materials in third-party furnaces, including but not limited to damage to or damage caused by the-sintered objects such as abutments, crowns or bridge frameworks.

Additional notes:

The programming you use for your sintering furnace should be analogous to the DeguDent -sintering programs.

General Sinter programs Cercon® ht and Cercon® xt all shades

Material:	Start-- Temp.	Ramp time	Temp. 1	Holding time	Ramp time	Temp. 2	Holding time	Cooling
	°C	min	°C	min	min	°C	min	
Sinter program for bridge frameworks up to 8 units (Cercon®- ht) and bridge frameworks up to 3 units (Cercon®- xt)	RT 1)	40	900	0	55	1500	145	With closed furnace cooling down to 200 °C
Sinter program for bridge frameworks for 9 or more units (Cercon®ht)	RT 1)	120	860	0	320	1500	120	With closed furnace cooling down to 200 °C
Speed-Sinter program for bridge frameworks up to 6 units (Cercon® ht) and for bridge frameworks- up to 3 units (Cercon® xt)	RT 1)	90	1540 1)	35	20	1150	0	Gradual opening of the furnace within 35 Min down to 200 °C

Material:	Start-- Temp.	Ramp time	Temp. 1	Holding time	Ramp time	Temp. 2	Holding time	Cooling
	°C	°C / min	°C	min	°C / min	°C	min	
Sinter program for bridge frameworks up to 8 units (Cercon®- ht) and bridge frameworks up to 3 units (Cercon®- xt)	RT 1)	22	900	0	1)	1500	145	With closed furnace cooling down to 200 °C
Sinter program for bridge frameworks for 9 or more units (Cercon®ht)	RT 1)	7	860	0	2	1500	120	With closed furnace cooling down to 200 °C
Speed-Sinter program for bridge frameworks up to 6 units (Cercon® ht) and for bridge frameworks- up to 3 units (Cercon® xt)	RT 1)	17	1540 1)	35	18	1150	0	Gradual opening of the furnace within 35 Min down to 200 °C

Material:	Start-- Temp.	Ramp time	Temp. 1	Holding time	Ramp time	Temp. 2	Holding time	Cooling
	°C	°C / h	°C	h : min	°C / h	°C	h : min	
Sinter program for bridge frameworks up to 8 units (Cercon®- ht) and bridge frameworks up to 3 units (Cercon®- xt)	RT 1)	1320	900	00 : 00	660	1500	02 : 15	With closed furnace cooling down to 200 °C
Sinter program for bridge frameworks for 9 or more units (Cercon®ht)	RT 1)	420	860	00 : 00	20	1500	02 : 00	With closed furnace cooling down to 200 °C
Speed-Sinter program for bridge frameworks up to 6 units (Cercon® ht) and for bridge frameworks- up to 3 units (Cercon® xt)	RT 1)	1020	1540 1)	00 : 35	1080	1150	00 : 00	Gradual opening of the furnace within 35 Min down to 200 °C

1) Room temperature 2) valid for closed sinter bowls, otherwise 1520 °C

Sintering programs, Multimat2Sinter / heat DUO / Sirona HTC-speed sintering furnace:

Speed sintering of bridge frameworks with up to 6 units \Cercon® ht\ and bridge frameworks up to 3 units \Cercon® xt\

Step	Heating rate	Temperature	Holding time
	°C / min	°C	min
S4	70	750	0
S3	70	1100	0
S2	70	1540	35
S1	70	0	0

Cercon® ht and Cercon®- xt standard sintering:

Step	Heating rate	Temperature	Holding time
	°C / min	°C	min
S4	31	300	0
S3	30	750	0
S2	11	1520	130
S1	22	880	0

Cercon® ht sintering of bridge frameworks with 8 or more units:

Step	Heating rate	Temperature	Holding time
	°C / min	°C	min
S4	11	170	0
S3	11	1100	0
S2	2	1520	125
S1	7	820	0

Sintering temperatures are recommendations. If necessary, carry out a trial sintering cycle and adapt the sintering temperatures or times as needed.

Separating the sintering reinforcement in the case of large-span bridges:
Objects are separated from the "tongue" after sintering using irrigated rotary diamond cutters.

Manual finishing after sintering:

Sandblast the inside and outside of the framework with aluminium oxide (110–125 µm, max. 2–3 bar, 45 ° angle).
Eliminate premature contacts one by one until the framework has reached its final position on the die(s).
During the try-on and fitting of the framework, keep the dies on the cast and try on the framework as a whole.
Once the try-on and fitting has been completed, do not perform any additional adjustments such as finishing- the entire framework.

Note: Zirconia crown or bridge abutments should fit passively, without friction. The rationale for requiring a passive- fit if zirconia frameworks lies in the physical properties of the material itself: Ceramics will tolerate compressive stress but not tensile stress. In crowns that exhibit a friction fit, this friction is generated- by certain "bearing" parts of the framework, as the crown itself never rests completely on the die- due to processing-- related surface roughness (the same actually holds for cast-metal frameworks). So it is only the peaks of the surface roughness (the named "bearing" parts) that are in contact with the die surface. This converts the -compressive forces generated by masticatory pressure into tensile forces, which may damage the crown.

Adapt the framework margins exactly to the preparation margins.

Note: Sintered zirconia should be finished using diamond instruments under proper irrigation only. - Keep the pressure on the framework material to a minimum and work only in one direction.

Sandblast the finished areas once again with aluminium oxide (110–125 µm, max. 2-3 bar, 45 ° angle).
Finally, clean the framework using a steam cleaner.

Veneering/staining:

Adjust your build-up technique to the respective -Cercon® ht/Cercon® xt shade.

Note that the light transmission of -Cercon® ht/Cercon®- xt is very high due to its translucency.

General

The reproduction of the individual patient's tooth colour can be significantly influenced by:

Shade of the die
Shade of the material used for temporary or definitive cementing
Framework wall thickness
Liner, if applied

Staining technique

For tooth-coloured staining of fully contoured Cer-con® ht and Cer-con® xt restorations we recommend Cer-con® ceram stains, Cer-con® TCT Stains, Celtra® Stains and Dentsply Sirona Universal Stains & Glaze.

Layering technique

We recommend -Cercon® ceram Kiss / Cercon® ceram press / -Cercon® ceram love/Celtra® Ceram ceramic -veneering -materials for veneering zirconia frameworks. Make sure to follow the applicable Instructions for Use.

Annealing:

Based on the results of our scientific testing of zirconia frameworks, we consider a separate annealing ("healing")- step both unnecessary and inappropriate.

Polishing in the laboratory:

Unveneered -Cercon® ht and -Cercon® xt frameworks should be highly polished or glazed with a high-gloss -ceramic -glaze. This also makes the restoration more amenable to proper oral hygiene procedures.

Polishing in the dental office:

Extensive studies have shown that the abrasive action of Cer-con® ht and Cercon® xt on antagonists is less - than that of con-ventional ceramic veneers and not more than standard lithium disilicate ceramics even after finishing and polishing.

Important: After making minor final occlusal adjustments at chairside, we recommend prophylactic polishing- of the adjusted sites to a high lustre or adding glaze with a smooth surface prior to temporary or definitive cementing- to protect the antagonists from possible abrasion. This also makes the restoration- more amenable to proper oral hygiene procedures.

Tooth surfaces whose shade was obtained by staining may present brighter spots where occlusal adjustments- were made.

Transport and storage:

No particular requirements.

Shelf Life:

Shelf life of Cercon® ht and Cercon® xt is 7 years from date of manufacture.

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DeguDent GmbH