

Bluephase[®] Curing light

Product care tips



Dos & Don'ts of Bluephase® Care

For the safety of all users and patients, proper cleaning is a top priority!

Cleaning the Bluephase curing light



✓ **DO** wipe the handpiece and charging base.



✓ **DO** use recommended disinfectants (see chart). Wring out extra moisture from the wipe before disinfecting the handpiece & charging base (excess moisture can lead to damage).

✓ **DO** disconnect the charging base from the power source when cleaning it.



✓ **DO** remove the light guide from the handpiece and place both horizontally for one minute to dry, after disinfecting the handpiece.



✗ **DON'T** remove the light guide before cleaning the handpiece (the light guide must be attached to the handpiece during cleaning so that no liquid enters the device).

✗ **DON'T** clean with highly aggressive disinfecting solutions (see chart).

Cleaning the light guide

✓ **DO** clean the light guide with a disinfectant and a soft cloth. Carefully remove any possible material residue. Partially polymerized composite can be removed with alcohol and a plastic spatula.

✓ **DO** sterilize the light guide before each use if disposable protective sleeves are not used. See Instructions for Use for further directions.

✗ **DON'T** use sharp or pointed objects, as they may scratch the surface.

Cleaning & Disinfecting agents

Acceptable cleaners:	DO NOT use agents containing:
The pH value of the cleaning and disinfecting solution should be between pH 7 and pH 11. (e.g. FD 366/Dürr Dental, Incidin Liquid/Eco-lab, CaviCideTM/Metrex Research)	Organic, mineral and oxidizing acids (minimum acceptable pH value is 5.5)
Isopropyl alcohol-based cleaners (<40% alcohol)	Alkaline solutions (maximum acceptable pH value is 11)
Ethyl alcohol-based cleaners (<40% alcohol)	Oxidizing agent (e.g. hydrogen peroxide)

*Damage to the Bluephase curing light caused by the use of cleaning products, including discolouration and corrosion, is not covered under the Bluephase curing light warranty.



WATCH THE VIDEO
for more product care tips.