



### CUTS A REALLY FINE FIGURE



**BEGO SEMADOS® RI-IMPLANTS**  
 QUICK & FUNCTIONAL · ECONOMICAL & MODERN  
 RELIABLE & PROVEN · VARIED & SOPHISTICATED  
 100% GERMAN QUALITY PRODUCT  
 INTUITIVE & SIMPLE · QUICK & FUNCTIONAL  
 ECONOMICAL & MODERN · RELIABLE & PROVEN  
 VARIED & SOPHISTICATED · INTUITIVE & SIMPLE

## The system for primary and long-term stability

The BEGO Semados® RI implant is an effective system when primary stability is vital. Supplemented by individual CAD/CAM-fabricated prosthetic components and the BEGO MultiPlus concept, a system that provides patients with a perfectly fitting denture in one day, the implants ensure standards ensure a high degree of satisfaction—for implant dentists and patients in equal measure.

#### Your benefits

- Microthread in the neck region—for improved force transfer to the crestal bone
- Increase in the primary stability due to the condensing thread design for lateral bone compression
- Well thought-out—cover screw included
- Integrated platform switch supports the preservation of crestal bone
- For the most stringent quality requirements —100% German development and manufacturing
- Stable connection with no micro-movements during functional masticatory loading (up to 150 N)

#### Technical details

- Condensing conical implant with platform switch and 0.5 mm machined implant shoulder
- Diameter: 3.75/4.1/4.5/5.5 mm
- Length: 8.5/10/11.5/13 mm
- Fabricated from grade 4 commercially pure titanium with a TiPure<sup>Plus</sup> surface
- Internal taper ensures no gaps
- Internal hex prevents rotation of abutments

# PRODUCT OVERVIEW

Distributor contact data

## BEGO Semados® RI-Implants

Art. no	Colour code	Product
---------	-------------	---------

### BEGO Semados® RI-Implants 3.75

57708	■	RI 3.75 L8.5
57709	■	RI 3.75 L10
57710	■	RI 3.75 L11.5
57711	■	RI 3.75 L13

### BEGO Semados® RI-Implants 4.1

57718	■	RI 4.1 L8.5
57719	■	RI 4.1 L10
57720	■	RI 4.1 L11.5
57721	■	RI 4.1 L13

### BEGO Semados® RI-Implants 4.5

57728	■	RI 4.5 L8.5
57729	■	RI 4.5 L10
57730	■	RI 4.5 L11.5
57731	■	RI 4.5 L13

### BEGO Semados® RI-Implants 5.5

57738	■	RI 5.5 L8.5
57739	■	RI 5.5 L10
57740	■	RI 5.5 L11.5
57741	■	RI 5.5 L13