

# READ ME!

## EXCITING CHANGES COMING

We are improving our current delivery system.

The new delivery system saves time and money while helping you achieve better results for your patients and your practice. The system is easier to use, provides less waste, and improves mixing.

### NEW IMPROVED DELIVERY SYSTEM available Spring 2021.

Look for updated material refill (no item number change) and NEW mix tip refills at your preferred authorized dealer.

Current Order Numbers for Mixing Tips		New Mixing Tips	
REF 72001	Green Mixing Tips (48 ct)	0072100FG	BLUE Mixing Tips (50 ct) <i>For Heavier Viscosity Materials (tray material)</i>
REF 72003	Yellow Mixing Tips (48 ct)	0072101FG	RED Mixing Tips (50 ct) <i>for Lower Viscosity Materials (wash material)</i>
REF 72007	Yellow Intraoral Tips (100 ct)	0072008FG	RED Intraoral Tips (100 ct) <i>for Lower Viscosity Materials (wash material)</i>
REF 78840	Dynamic Mixing Tips (50 ct)	0072102FG	Dynamic Mixing Tips (50 ct) <i>use with Machine Mix Impression Material</i>

#### Important to note about your new delivery system:

- ✓ Can be used with your current dispensing gun.  
No special additional dispenser is required.
- ✓ New mixing tips will **ONLY** fit the new system
- ✓ New mixing tip will **NOT** fit your old cartridge system

#### INTRODUCTORY BONUS

20% more mixing tips  
in new mixing tip refill vs. original packaging




Pour lire ceci en français,  
veuillez numériser ce code QR

# NEW CARTRIDGE DELIVERY SYSTEM

Can eliminate **up to 24%** of the waste you're used to.\*


Plus, you only need to bleed once for the life of the cartridge.



**Improved Mix Quality, More Homogeneous Mix**  
with the Stonehenge Mixing Chamber Technology.  
Ensures proper mix from start to finish



**For 5:1 Cartridges**  
No Bayonet Ring Necessary



\*when comparing Genie impression material with the new blue mixing tip versus the old green mixing tip  
MLS10008F (7-15-21)