



EFFEGI BREGA®

“ COMPACT ”

- POLYMERIZING UNIT -

The great skilled polymerizing unit, modern, compact and functional. A long lasting product that respects the tradition of solidity and practicalness.

Outlet RS232

Connection to printer. It allows to connect the Polyprinter. The released printing shows in details the data used in the applied cycle. A control and a written certification for the most demanding customers.

Display

It shows:
Cycle start time, pressure, temperature, maintenance, cooling, water/air discharge, data print.

Chamber

- Big capacity: a normal articulator can be inserted.
- Equipped with basket, cup and support with tweezers.
- Equipped with a safety device that prevents the opening when the machine is in pressure.



Uses

- It polymerizes any kind of resine.
- Hydroscopic bath.
- Pressure hydroscopic bath.
- Removing wax from the cylinders.
- Steamer
- Plaster coating press.
- Resin cooking in slow-polymerization muffle.
- 10 programmable cycles.
- Programmable switch.
- Programmable timer.

Control buttons

They allow: :

- To program the different cycles.
- To visualize the set data..
- To print the report of the cycle.
- To control the appliance.