

Adapta deep-drawing system

- Simple and fast deep-drawing of copings
- Reasonably priced system which has been proven over many years with special plastic foils
- A uniform minimum wall thickness ensures a high level of stability in the copings
- Ideal for the double-crown technique; coated with milling wax, the Adapta coping offers a high level of protection against inadvertent milling through
- The thin 0.1 mm spacer foil, which is used as part of the system, frees up the necessary, defined space for the luting material

Product details				
Availability		Presentation	Qty	REF
Adapta deep drawing system comprising:		set	1	20500
1 Forming tub with Adapta mastic				
1 Spare pack Adapta mastic				
1 Foil holder				
100 Adapta foils, 0.6 mm in foil dispenser				
1 Pack, 100 Adapta foils, 0.6 mm				
200 Adapta foils, 0.1 mm red, in foil dispenser				
Adapta deep drawing system intro set comprising:		set	1	20520
1 Forming tub with Adapta mastic				
1 Foil holder				
50 Adapta foils, 0.6 mm				
50 Spacer foils, 0.1 mm				
Accessories	Contents	Presentation	Qty	REF
Adapta mastic (spare pack)		pack	1	20503
Forming tub with Adapta mastic, 1 Forming tub		piece	1	20504
Adapta foil holder		piece	1	20510
Adapta foil dispenser incl. 100 x 0.6 mm		pack	1	20519
Adapta foil dispenser incl. 200 x 0.1 mm		pack	1	20521
Adapta Spacer foils, 0.1 mm transparent	200 pcs.	pack	1	20517
Adapta Spacer foils 0.1 mm red	200 pcs.	pack	1	20502
Adapta foils 0.6 mm, transparent	100 pcs.	pack	1	20501



Rapidi modelling knife

- The Rapidi modelling knife is ideal for cutting, scraping and modelling
- Easy-to-change blade

Product details					
Availability		Contents	Presentation	Qty	REF
Rapidi modelling knife			piece	1	52270
Rapidi spare blades	40 pieces		set	1	52280

