

Trade name: BegoForm (REF 52785)

Current version: 2.0.3, issued: 19.10.2015 Replaced version: 2.0.2, issued: 16.12.2014 Region: GB

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1 Product identifier

Trade name

BegoForm (REF 52785)

1.2 Relevant identified uses of the substance or mixture and uses advised against

Relevant identified uses of the substance or mixture

Manufacturing of dental prosthesis in a dental laboratory

Uses advised against

No data available.

1.3 Details of the supplier of the safety data sheet

Address

BEGO Bremer Goldschlägerei Wilh. Herbst GmbH & Co. KG Wilhelm-Herbst-Str. 1

28359 Bremen
Telephone no. +49/ 421/ 2028 – 0

Fax no. +49/ 421/ 2028 – 115 e-mail msds@bego.com

Information provided by / telephone

Research & Development Department - Materials, alloys and ceramics; +49/ 421/ 2028 – 130 (Chief Development Officer alloys)

Advice on Safety Data Sheet

msds@bego.com

1.4 Emergency telephone number

For medical advice (in German and English): +49 (0)551 192 40 (Giftinformationszentrum Nord)

SECTION 2: Hazards identification

2.1 Classification of the substance or mixture

Classification in accordance with Regulation (EC) No 1272/2008 (CLP)

STOT RE 2: H373i

2.2 Label elements

<u>Labelling according to Regulation (EC) No 1272/2008 (CLP Regulation)</u>

Hazard pictograms



Signal word

Warning

Hazardous component(s) to be indicated on label:

Quartz (SiO2)

Hazard statements

H373i May cause damage to organs through prolonged or repeated exposure if inhaled.

Precautionary statements

P260 Do not breathe dust.

P314 Get medical advice/attention if you feel unwell.

P403+P233 Store in a well-ventilated place. Keep container tightly closed.



Trade name: BegoForm (REF 52785)

Current version: 2.0.3, issued: 19.10.2015 Replaced version: 2.0.2, issued: 16.12.2014 Region: GB

2.3 Other hazards

PBT assessment No data available.

vPvB assessment No data available.

SECTION 3: Composition/information on ingredients

3.1 Substances

Not applicable. The product is not a substance.

3.2 Mixtures

Hazardous ingredients

No	Substance name		Additi	onal information	1	
	CAS / EC / Index /	Classification (EC) 1272/2008 (CLP)	Conce	entration		%
	REACH no					
1	Quartz (SiO2)					
	14808-60-7	STOT RE 2; H373i	>	10.00 - <	50.00	%-b.w.
	238-878-4					
	-					
	-					
2	cristobalite					
	14464-46-1	STOT RE 2; H373i	>	5.00 - <	20.00	%-b.w.
	238-455-4					
	-					
	-					
3	zinc oxide					
	1314-13-2	Aquatic Acute 1; H400	<	0.50		%-b.w.
	215-222-5	Aquatic Chronic 1; H410				
	030-013-00-7					
	01-2119463881-32					
4	aluminium oxide					
	1344-28-1	-	>	10.00 - <	30.00	%-b.w.
	215-691-6					
	-					
	01-2119529248-35					
5	MAGNESIUM OXID	E				
	1309-48-4	-	>	10.00 - <	30.00	%-b.w.
	215-171-9					
	-					
	-					

Full Text for all H-phrases and EUH-phrases: pls. see section 16

No	Note	Specific concentration limits	M-factor (acute)	M-factor (chronic)
3	-	-	M = 1	-

No	Route, target organ, concrete effect
1	H373i
	inhalational; -; -
2	H373i
	inhalational; -; -

SECTION 4: First aid measures

4.1 Description of first aid measures

General information

If the patient is likely to become unconscious, place and transport in stable sideways position. In case of persisting adverse effects, consult a physician. Remove contaminated clothing and shoes immediately, and launder thoroughly before reusing.



Trade name: BegoForm (REF 52785)

Current version: 2.0.3, issued: 19.10.2015 Replaced version: 2.0.2, issued: 16.12.2014 Region: GB

After inhalation

Remove affected person from the immediate area. Ensure supply of fresh air. Irregular breathing/no breathing: artificial respiration. Call a doctor immediately.

After skin contact

In case of contact with skin wash off with water.

After eye contact

Remove contact lenses, irrigate copiously with clean, fresh water for at least 15 minutes, holding the eyelids apart and seek medical advice.

After ingestion

Rinse mouth thoroughly with water. Never give anything by mouth to an unconscious person. Call a doctor immediately.

4.2 Most important symptoms and effects, both acute and delayed

No data available.

4.3 Indication of any immediate medical attention and special treatment needed

No data available.

SECTION 5: Firefighting measures

5.1 Extinguishing media

Suitable extinguishing media

Carbon dioxide; Extinguishing powder; Water spray jet; Foam

Unsuitable extinguishing media

High power water jet

5.2 Special hazards arising from the substance or mixture

No data available.

5.3 Advice for firefighters

Product itself does not burn. Adapt extinguisher and fire-fighting measures to fire in the environment. Use self-contained breathing apparatus. Wear protective clothing.

SECTION 6: Accidental release measures

6.1 Personal precautions, protective equipment and emergency procedures

For non-emergency personnel

Avoid dust formation. Refer to protective measures listed in sections 7 and 8.

For emergency responders

No data available. Personal protective equipment (PPE) - see Section 8.

6.2 Environmental precautions

Do not discharge into the drains/surface waters/groundwater. Do not discharge into the subsoil/soil.

6.3 Methods and material for containment and cleaning up

Take up mechanically. Avoid raising dust. When picked up, treat material as prescribed under heading "Disposal considerations".

6.4 Reference to other sections

No data available.

SECTION 7: Handling and storage

7.1 Precautions for safe handling

Advice on safe handling

Keep container tightly closed.

General protective and hygiene measures

Wash hands before breaks and after work. Do not eat, drink or smoke during work time. Remove soiled or soaked clothing immediately. Do not inhale dust. Keep away from foodstuffs and beverages.



Trade name: BegoForm (REF 52785)

Current version: 2.0.3, issued: 19.10.2015 Replaced version: 2.0.2, issued: 16.12.2014 Region: GB

Advice on protection against fire and explosion

No special measures necessary.

7.2 Conditions for safe storage, including any incompatibilities

Technical measures and storage conditions

Keep container tightly closed and dry in a cool, well-ventilated place.

Requirements for storage rooms and vessels

Keep only in the original container.

Advice on storage assembly

Protect against humidity; Do not store together with: explosive substances

7.3 Specific end use(s)

No data available.



Trade name: BegoForm (REF 52785)

Current version: 2.0.3, issued: 19.10.2015 Replaced version: 2.0.2, issued: 16.12.2014 Region: GB

SECTION 8: Exposure controls/personal protection

8.1 Control parameters

Occupational exposure limit values

No	Substance name	CAS no.		EC no.	
1	Dust				
_	List of approved workplace exposure limits (WELs) / I	EH40			
	total inhalable dust				
	TWA	10	mg/m³		
	List of approved workplace exposure limits (WELs) / I				
	respirable dust				
	TWA	4	mg/m³		
2	aluminium oxide	1344-28-1	<i>-</i>	215-691-6	
	List of approved workplace exposure limits (WELs) / I	EH40			
	Aluminium oxides				
	total inhalable dust				
	TWA	10	mg/m³		
	List of approved workplace exposure limits (WELs) / I	EH40	•		
	Aluminium oxides				
	respirable dust				
	TWA	4	mg/m³		
3	MAGNESIUM OXIDE	1309-48-4	_	215-171-9	
	List of approved workplace exposure limits (WELs) / I	EH40			
	Magnesium oxide				
	inhalable dust				
	Mg				
	TWA	10	mg/m³		
	List of approved workplace exposure limits (WELs) / I	EH40			
	Magnesium oxide				
	fume and respirable dust				
	Mg	Τ.			
4	TWA	4	mg/m³	204 205 2	
4	ammonia, anhydrous	7664-41-7		231-635-3	
	List of approved workplace exposure limits (WELs) / I	EH40			
	Ammonia STEL	25	m a /m 3	25	ml/m³
	TWA	18	mg/m³ mg/m³	35 25	ml/m³
	2000/39/EWG	10	mg/m²	25	mi/m²
	Ammonia, anhydrous				
	STEL	36	mg/m³	50	ml/m³
	TWA	14	mg/m³	20	ml/m³
5	zinc oxide	1314-13-2	mg/m²	215-222-5	1111/1111-
3	List of approved workplace exposure limits (WELs) / I			£13-222-3	
	Zinc oxide, fume				
	STEL	10	mg/m³		
	TWA	5	mg/m³		
	IVA	J	mg/m ⁻		



Trade name: BegoForm (REF 52785)

Current version: 2.0.3, issued: 19.10.2015 Replaced version: 2.0.2, issued: 16.12.2014 Region: GB

DNEL and PNEC values

DNEL values (worker)

No	Substance name	Substance name			C no
	Route of exposure	Exposure time	Effect	Value	
1	zinc oxide			1314-13 215-222	
	dermal	Long term (chronic)	systemic	83	mg/kg/day
	with reference to: Zn Remarks: insoluble				
	inhalative	Long term (chronic)	systemic	5	mg/m³
	with reference to: Zn Remarks: insoluble				

DNEL value (consumer)

No	Substance name			CAS / EC	no
	Route of exposure	Exposure time	Effect	Value	
1	zinc oxide			1314-13- 215-222-	—
	oral	Long term (chronic)	systemic	0.83	mg/kg/day
	with reference to: Zn Remarks: insoluble				
	dermal	Long term (chronic)	systemic	83	mg/kg/day
	with reference to: Zn Remarks: insoluble				
	inhalative	Long term (chronic)	systemic	2.5	mg/m³
	with reference to: Zn Remarks: insoluble				

PNEC values

No	Substance name		CAS / EC no	
	ecological compartment	Туре	Value	
1	zinc oxide		1314-13-2	
			215-222-5	
	water	fresh water	20.6	μg/l
	with reference to: Zn			
	water	marine water	6.1	μg/l
	with reference to: Zn			
	water	fresh water sediment	117.8	mg/kg
	water	marine water sediment	56.5	mg/kg
	with reference to: Zn, dry weight			
	soil	-	35.6	mg/kg
	with reference to: Zn, dry weight			
	sewage treatment plant	-	100	μg/l

8.2 Exposure controls

Appropriate engineering controls

No data available.

Personal protective equipment

Respiratory protection

If workplace exposure limits are exceeded, a respiration protection approved for this particular job must be worn. In case of dust formation, take appropriate measures for breathing protection in the event workplace threshold values are not specified.

Eye / face protection

Safety glasses with side protection shield (EN 166)



Trade name: BegoForm (REF 52785)

Current version: 2.0.3, issued: 19.10.2015 Replaced version: 2.0.2, issued: 16.12.2014 Region: GB

Hand protection

Sufficient protection is given wearing suitable protective gloves checked according to i.e. EN 374, in the event of risk of skin contact with the product. Before use, the protective gloves should be tested in any case for its specific workstation suitability (i.e. mechanical resistance, product compatibility and antistatic properties). Adhere to the manufacturer's instructions and information relating to the use, storage, care and replacement of protective gloves. Protective gloves shall be replaced immediately when physically damaged or worn. Design operations thus to avoid permanent use of protective gloves.

Other

Normal chemical work clothing.

Environmental exposure controls

No data available.

SECTION 9: Physical and chemical properties

9.1 Information on basic physical and chemical properties

Form/Colour				
Powder				
white				
Odour				
odourless				
Odour threshold				
No data available				
pH value				
Value	5	- 7		
Concentration		80	% (water)	



Trade name: BegoForm (REF 52785)

Current version: 2.0.3, issued: 19.10.2015 Replaced version: 2.0.2, issued: 16.12.2014 Region: GB

Boiling point / boiling range

No data available

Melting point / melting range

No data available

Decomposition point / decomposition range

No data available

Flash point

No data available

Auto-ignition temperature

No data available

Oxidising properties

No data available

Explosive properties

No data available

Flammability (solid, gas)

No data available

Lower flammability or explosive limits

No data available

Upper flammability or explosive limits

No data available

Vapour pressure

No data available

Vapour density

No data available

Evaporation rate

No data available

Relative density

No data available

Density

No data available

Bulk density

Value 1.2 - 1.5 g/cm³ Reference temperature 21 °C

Solubility in water

No data available

Solubility(ies)

No data available

Partition coefficient: n-octanol/water

No data available

Viscosity

No data available

9.2 Other information

Other information

The product is hygroscopic.

Product hardens with the addition of water.



Trade name: BegoForm (REF 52785)

Current version: 2.0.3, issued: 19.10.2015 Replaced version: 2.0.2, issued: 16.12.2014 Region: GB

SECTION 10: Stability and reactivity

10.1 Reactivity

Dangerous reactions are not expected if the product is handled according to its intended use.

10.2 Chemical stability

Stable under recommended storage and handling conditions (See section 7).

10.3 Possibility of hazardous reactions

Dangerous reactions are not to be expected when handling product according to its intended use.

10.4 Conditions to avoid

Protect from humid air and water.

10.5 Incompatible materials

None known.

10.6 Hazardous decomposition products

Under normal heating of the hardened ground-mass in the pre-heat furnace ammonia is released with a concentration of 1 percent by weight. Provide ventilation to exhaust vapours and conduct them directly to the outside.

SECTION 11: Toxicological information

11.1 Information on toxicological effects

Acute oral toxicity	
No data available	

Acute dermal toxicity

No data available

Acute inhalational toxicity

No data available

Skin corrosion/irritation

No data available

Serious eye damage/irritation

No data available

Respiratory or skin sensitisation

No data available

Germ cell mutagenicity

No data available

Reproduction toxicity

No data available

Carcinogenicity

No data available

STOT-single exposure

No data available

STOT-repeated exposure

No data available

Aspiration hazard

No data available

Delayed and immediate effects as well as chronic effects from short and long-term exposure

Contact with the skin and eyes may cause mechanical irritation. Inhalation of dusts may irritate the respiratory tract. Can cause silicosis.



Trade name: BegoForm (REF 52785)

Current version: 2.0.3, issued: 19.10.2015 Replaced version: 2.0.2, issued: 16.12.2014 Region: GB

SECTION 12: Ecological information

12.1 Toxicity

No data available

Toxicity to fish (chronic)

No data available

Toxicity to Daphnia (acute)

No data available

Toxicity to Daphnia (chronic)

No data available

Toxicity to algae (acute)

No data available

Toxicity to algae (chronic)

No data available

Bacteria toxicity

No data available

12.2 Persistence and degradability

No data available.

12.3 Bioaccumulative potential

No data available.

12.4 Mobility in soil

No data available.

12.5 Results of PBT and vPvB assessment

Results of PBT and vPvB assessment		
PBT assessment	No data available.	
vPvB assessment	No data available.	

12.6 Other adverse effects

No data available.

12.7 Other information

Other information

Do not discharge product unmonitored into the environment.

SECTION 13: Disposal considerations

13.1 Waste treatment methods

Product

Allocation of a waste code number, according to the European Waste Catalogue, should be carried out in agreement with the regional waste disposal company.

Packaging

Residuals must be removed from packaging and when emptied completely disposed of in accordance with the regulations for waste removal. Incompletely emptied packaging must be disposed of in the form of disposal specified by the regional disposer.



Trade name: BegoForm (REF 52785)

Current version: 2.0.3, issued: 19.10.2015 Replaced version: 2.0.2, issued: 16.12.2014 Region: GB

SECTION 14: Transport information

14.1 Transport ADR/RID/ADN

The product is not subject to ADR/RID/ADN regulations.

14.2 Transport IMDG

The product is not subject to IMDG regulations.

14.3 Transport ICAO-TI / IATA

The product is not subject to ICAO-TI / IATA regulations.

14.4 Other information

No data available.

14.5 Environmental hazards

Information on environmental hazards, if relevant, please see 14.1 - 14.3.

14.6 Special precautions for user

No data available.

14.7 Transport in bulk according to Annex II of Marpol and the IBC Code

Not relevant

SECTION 15: Regulatory information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture EU regulations

Regulation (EC) No 1907/2006 (REACH) Annex XIV (List of substances subject to authorisation)

According to the data available and/or specifications supplied by upstream suppliers, this product does not contain any substances considered as substances requiring authorisation as listed on Annex XIV of the REACH regulation (EC) 1907/2006.

REACH candidate list of substances of very high concern (SVHC) for authorisation

According to available data and the information provided by preliminary suppliers, the product does not contain substances that are considered substances meeting the criteria for inclusion in annex XIV (List of Substances Subject to Authorisation) as laid down in Article 57 and article 59 of REACH (EC) 1907/2006.

Regulation (EC) No 1907/2006 (REACH) Annex XVII: RESTRICTIONS ON THE MANUFACTURE, PLACING ON THE MARKET AND USE OF CERTAIN DANGEROUS SUBSTANCES, PREPARATIONS AND ARTICLES

According to the data available and/or specifications supplied by upstream suppliers, this product does not contain any substances requiring authorisation as listed on Annex XIV of the REACH regulation (EC) 1907/2006.

Directive 2012/18/EU on the control of major-accident hazards involving dangerous substances

This product is not subject to Part 1 or 2 of Annex I.

Other regulations

Adhere to the national sanitary and occupational safety regulations when using this product.

15.2 Chemical safety assessment

No data available.

SECTION 16: Other information

Sources of key data used to compile the data sheet:

Regulation (EC) No 1907/2006 (REACH), 1272/2008 (CLP) as amended in each case.

EC Directives 2000/39/EC, 2006/15/EC, 2009/161/EU

National Threshold Limit Values of the corresponding countries as amended in each case.

Transport regulations according to ADR, RID, IMDG, IATA as amended in each case.

The data sources used to determine physical, toxic and ecotoxic data, are indicated directly in the corresponding chapter.



Trade name: BegoForm (REF 52785)

Current version: 2.0.3, issued: 19.10.2015 Replaced version: 2.0.2, issued: 16.12.2014 Region: GB

Full text of the H- and EUH- phrases drawn up in sections 2 and 3 (provided not already drawn up in these sections)

H400 Very toxic to aquatic life.

H410 Very toxic to aquatic life with long lasting effects.

Department issuing safety data sheet

UMCO Umwelt Consult GmbH - D-21107 Hamburg, Georg-Wilhelm-Strasse 183, Tel.: +49(40)79 02 36 300, Fax: +49(40)79 02 36 357, e-mail: umco@umco.de

This information is based on our present state of knowledge and experience.

The security data sheet describes products with a view to the security requirements.

However, it should not constitute a guarantee for any specific product properties and shall not establish a legally valid relationship.

Document protected by copyright. Alterations or reproductions require the express written permission of UMCO Umwelt Consult GmbH.