🔿 BEGO

Current version : 2.0.3, issued: 19.10.2015

Replaced version: 2.0.2, issued: 16.12.2014

Region: GB

SECTION 1: Identification of the substance/mixture and of the company/undertaking

#### 1.1 Product identifier

Trade name

### BegoStone plus (REF 54811)

 Identification numbers

 CAS no.
 10034-76-1

 EC no.
 231-900-3

#### 1.2 Relevant identified uses of the substance or mixture and uses advised against

Relevant identified uses of the substance or mixture Manufacturing of dental prosthesis in a dental laboratory Uses advised against

No data available.

#### **1.3** Details of the supplier of the safety data sheet

#### Address

BEGO Bremer Goldschlägerei Wilh. Herbst GmbH & Co. KG Wilhelm-Herbst-Str. 1 28359 Bremen

 Telephone no.
 +49/421/2028 - 0

 Fax no.
 +49/421/2028 - 115

 e-mail
 msds@bego.com

#### Information provided by / telephone

Research & Development Department - Materials, alloys and ceramics; +49/ 421/ 2028 – 130 (Chief Development Officer alloys)

# Advice on Safety Data Sheet msds@bego.com

#### 1.4 Emergency telephone number

For medical advice (in German and English): +49 (0)551 192 40 (Giftinformationszentrum Nord)

#### **SECTION 2: Hazards identification**

#### 2.1 Classification of the substance or mixture

#### **Classification information**

This product does not meet the classification and labelling criteria given in the Regulation (EC) No 1272/2008 (CLP). This product does not meet the classification and labelling criteria given in the Dangerous Substances Directive 67/548/EC (DSD)

#### 2.2 Label elements Not relevant

# 2.3 Other hazards

PBT assessment No data available. vPvB assessment

No data available.

### **SECTION 3: Composition/information on ingredients**

#### 3.1 Substances

#### Chemical characterization

Substance name Calciumsulfathemihydrat

Current version : 2.0.3, issued: 19.10.2015

**Replaced version:** 2.0.2, issued: 16.12.2014

Region: GB

BEG

#### Identification numbers

CAS no.	10034-76-1	
EC no.	231-900-3	

#### 3.2 Mixtures

Not applicable. The product is not a mixture.

#### **SECTION 4: First aid measures**

#### 4.1 Description of first aid measures

#### **General information**

If the patient is likely to become unconscious, place and transport in stable sideways position. In case of persisting adverse effects, consult a physician. Remove contaminated clothing and shoes immediately, and launder thoroughly before reusing.

#### After inhalation

Remove affected person from the immediate area. Ensure supply of fresh air. Irregular breathing/no breathing: artificial respiration. Call a doctor immediately.

#### After skin contact

In case of contact with skin wash off with water.

#### After eve contact

Remove contact lenses, irrigate copiously with clean, fresh water for at least 15 minutes, holding the eyelids apart and seek medical advice.

#### After ingestion

Rinse mouth thoroughly with water. Never give anything by mouth to an unconscious person. Let plenty of water be drunk in small gulps. Call a doctor immediately.

#### **4.2 Most important symptoms and effects, both acute and delayed** No data available.

# **4.3** Indication of any immediate medical attention and special treatment needed No data available.

#### **SECTION 5: Firefighting measures**

#### 5.1 Extinguishing media

#### Suitable extinguishing media Carbon dioxide; Extinguishing powder; Water spray jet; Foam

Unsuitable extinguishing media High power water jet

5.2 Special hazards arising from the substance or mixture No data available.

#### 5.3 Advice for firefighters

Product itself does not burn. Adapt extinguisher and fire-fighting measures to fire in the environment. Use selfcontained breathing apparatus. Wear protective clothing.

#### **SECTION 6: Accidental release measures**

#### 6.1 Personal precautions, protective equipment and emergency procedures

#### For non-emergency personnel

Avoid dust formation. Refer to protective measures listed in sections 7 and 8.

#### For emergency responders

No data available. Personal protective equipment (PPE) - see Section 8.

#### 6.2 Environmental precautions

Do not discharge into the drains/surface waters/groundwater. Do not discharge into the subsoil/soil.

Current version : 2.0.3. issued: 19.10.2015



Trade name: BegoStone plus (REF 54811)

Replaced version: 2.0.2, issued: 16.12.2014

Region: GB

# 6.3 Methods and material for containment and cleaning up

Take up mechanically. Avoid raising dust. When picked up, treat material as prescribed under heading "Disposal considerations".

6.4 Reference to other sections No data available.

#### **SECTION 7: Handling and storage**

#### 7.1 Precautions for safe handling

Advice on safe handling Keep container tightly closed.

#### General protective and hygiene measures

Wash hands before breaks and after work. Do not eat, drink or smoke during work time. Remove soiled or soaked clothing immediately. Do not inhale dust. Keep away from foodstuffs and beverages.

Advice on protection against fire and explosion

No special measures necessary.

#### 7.2 Conditions for safe storage, including any incompatibilities

#### Technical measures and storage conditions

Keep container tightly closed and dry in a cool, well-ventilated place.

**Requirements for storage rooms and vessels** Keep only in the original container.

Advice on storage assembly Protect against humidity; Do not store together with: explosive substances

#### 7.3 Specific end use(s)

No data available.

#### **SECTION 8: Exposure controls/personal protection**

#### 8.1 Control parameters

#### **Occupational exposure limit values**

No	Substance name	CAS no.		EC no.
1	Dust			
	List of approved workplace exposure limits (WELs) / E	EH40		
	total inhalable dust			
	TWA	10	mg/m³	
	List of approved workplace exposure limits (WELs) / E	EH40		
	respirable dust			
	TWA	4	mg/m³	

#### 8.2 Exposure controls

Appropriate engineering controls No data available.


#### Personal protective equipment

#### **Respiratory protection**

If workplace exposure limits are exceeded, a respiration protection approved for this particular job must be worn. In case of dust formation, take appropriate measures for breathing protection in the event workplace threshold values are not specified.

#### Eye / face protection

Safety glasses with side protection shield (EN 166)



Current version : 2.0.3, issued: 19.10.2015

Replaced version: 2.0.2, issued: 16.12.2014

Region: GB

Hand protection

Sufficient protection is given wearing suitable protective gloves checked according to i.e. EN 374, in the event of risk of skin contact with the product . Before use, the protective gloves should be tested in any case for its specific work-station suitability (i.e. mechanical resistance, product compatibility and antistatic properties). Adhere to the manufacturer's instructions and information relating to the use, storage, care and replacement of protective gloves. Protective gloves shall be replaced immediately when physically damaged or worn. Design operations thus to avoid permanent use of protective gloves.

#### Other

Normal chemical work clothing.

Environmental exposure controls

No data available.

**SECTION 9: Physical and chemical properties** 

#### 9.1 Information on basic physical and chemical properties

Form/Colour		
Powder		
ivory colours		
Odour		
odourless		
Odour threshold		
No data available		
pH value		
No data available		
Boiling point / boiling range		
No data available		

Melting point / melting range

No data available



Replaced version: 2.0.2, issued: 16.12.2014

Region: GB

BEG

Decomposition point / decomposition range
No data available
Flash point
No data available
Auto-ignition temperature
No data available
Oxidising properties
No data available
Explosive properties
No data available
Flammability (solid, gas)
No data available
Lower flammability or explosive limits
No data available
Upper flammability or explosive limits
No data available
Vapour pressure
No data available
Vapour density
No data available
Evaporation rate
No data available
Relative density
No data available
Density
No data available
Solubility in water
Remarks insoluble
Solubility(ies)
No data available
Partition coefficient: n-octanol/water
No data available
Viscosity
No data available
2 Other information

## 9

Other information No data available.

#### **SECTION 10: Stability and reactivity**

#### 10.1 Reactivity

Dangerous reactions are not expected if the product is handled according to its intended use.

#### 10.2 Chemical stability

Stable under recommended storage and handling conditions (See section 7).

#### 10.3 Possibility of hazardous reactions None, if handled according to order.

Conditions to avoid 10.4

Protect from humid air and water.



Current version : 2.0.3, issued: 19.10.2015

Replaced version: 2.0.2, issued: 16.12.2014

Region: GB

#### 10.5 Incompatible materials

None known.

10.6 Hazardous decomposition products No hazardous decomposition products known.

#### **SECTION 11: Toxicological information**

#### 11.1 Information on toxicological effects

Acute oral toxicity
No data available
Acute dermal toxicity
No data available
Acute inhalational toxicity
No data available
Skin corrosion/irritation
No data available
Serious eye damage/irritation
No data available
Respiratory or skin sensitisation
No data available
Germ cell mutagenicity
No data available
Reproduction toxicity
No data available
Carcinogenicity
No data available
STOT-single exposure
No data available
STOT-repeated exposure
No data available
Aspiration hazard
No data available
Delayed and immediate effects as well as chronic effects from short and long-term exposure
Contact with the skin and eyes may cause mechanical irritation. Inhalation of dusts may irritate the respiratory tract.

Current version : 2.0.3, issued: 19.10.2015

#### **SECTION 12: Ecological information**

#### 12.1 Toxicity

Toxicity to fi	sh (acute)
----------------	------------

No data available

Toxicity to fish (chronic) No data available

## Toxicity to Daphnia (acute)

No data available

### Toxicity to Daphnia (chronic)

No data available

### Toxicity to algae (acute)

No data available

#### Toxicity to algae (chronic)

No data available

#### Bacteria toxicity

No data available

#### 12.2 Persistence and degradability

No data available.

#### **12.3 Bioaccumulative potential** No data available.

#### 12.4 Mobility in soil

No data available.

#### 12.5 Results of PBT and vPvB assessment

Results of PBT and vPvB assessment	
PBT assessment	No data available.
vPvB assessment	No data available.

#### 12.6 Other adverse effects

No data available.

### 12.7 Other information

Other information

Do not discharge product unmonitored into the environment.

#### **SECTION 13: Disposal considerations**

#### 13.1 Waste treatment methods

#### Product

Allocation of a waste code number, according to the European Waste Catalogue, should be carried out in agreement with the regional waste disposal company.

#### Packaging

Residuals must be removed from packaging and when emptied completely disposed of in accordance with the regulations for waste removal. Incompletely emptied packaging must be disposed of in the form of disposal specified by the regional disposer.





Current version : 2.0.3, issued: 19.10.2015

Replaced version: 2.0.2, issued: 16.12.2014

Region: GB

#### **SECTION 14: Transport information**

**14.1 Transport ADR/RID/ADN** The product is not subject to ADR/RID/ADN regulations.

#### **14.2 Transport IMDG** The product is not subject to IMDG regulations.

- **14.3 Transport ICAO-TI / IATA** The product is not subject to ICAO-TI / IATA regulations.
- **14.4** Other information No data available.
- **14.5** Environmental hazards Information on environmental hazards, if relevant, please see 14.1 - 14.3.
- **14.6** Special precautions for user No data available.
- 14.7 Transport in bulk according to Annex II of Marpol and the IBC Code Not relevant

#### **SECTION 15: Regulatory information**

# 15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture EU regulations

#### Regulation (EC) No 1907/2006 (REACH) Annex XIV (List of substances subject to authorisation)

In accordance with the Reach regulation (EC) 1907/2006, the product does not contain any substances that are considered as subject to listing in annex XIV, inventory of substances requiring authorisation.

REACH candidate list of substances of very high concern (SVHC) for authorisation

In accordance with article 57 and article 59 of the Reach regulation (EC) 1907/2006, this substance is not considered as subject to listing in annex XIV, inventory of substances requiring authorisation ("Authorization list").

Regulation (EC) No 1907/2006 (REACH) Annex XVII: RESTRICTIONS ON THE MANUFACTURE, PLACING ON THE MARKET AND USE OF CERTAIN DANGEROUS SUBSTANCES, PREPARATIONS AND ARTICLES The substance is not subject to the provisions of annex XVII (restriction entries) of the Reach regulation (EC) 1907/2006.

Directive 2012/18/EU on the control of major-accident hazards involving dangerous substances This substance is not subject to Part I or 2 of Annex I

#### Other regulations

Adhere to the national sanitary and occupational safety regulations when using this product.

#### 15.2 Chemical safety assessment

No data available.

#### **SECTION 16: Other information**

#### Sources of key data used to compile the data sheet:

Regulation (EC) No 1907/2006 (REACH), 1272/2008 (CLP) as amended in each case.

EC Directives 2000/39/EC, 2006/15/EC, 2009/161/EU

National Threshold Limit Values of the corresponding countries as amended in each case.

Transport regulations according to ADR, RID, IMDG, IATA as amended in each case.

The data sources used to determine physical, toxic and ecotoxic data, are indicated directly in the corresponding chapter.

#### Department issuing safety data sheet

UMCO Umwelt Consult GmbH - D-21107 Hamburg, Georg-Wilhelm-Strasse 183, Tel.: +49(40)79 02 36 300, Fax: +49(40)79 02 36 357, e-mail: umco@umco.de



Current version : 2.0.3, issued: 19.10.2015

Replaced version: 2.0.2, issued: 16.12.2014

Region: GB

This information is based on our present state of knowledge and experience. The security data sheet describes products with a view to the security requirements. However, it should not constitute a guarantee for any specific product properties and shall not establish a legally valid relationship.

Document protected by copyright. Alterations or reproductions require the express written permission of UMCO Umwelt Consult GmbH.