

Trade name: Wiropaint plus (REF 51100)

Current version: 2.0.3, issued: 19.10.2015 Replaced version: 2.0.2, issued: 16.12.2014 Region: GB

# SECTION 1: Identification of the substance/mixture and of the company/undertaking

#### 1.1 Product identifier

Trade name

# Wiropaint plus (REF 51100)

## 1.2 Relevant identified uses of the substance or mixture and uses advised against

#### Relevant identified uses of the substance or mixture

Manufacturing of dental prosthesis in a dental laboratory

#### Uses advised against

No data available.

#### 1.3 Details of the supplier of the safety data sheet

#### Address

BEGO Bremer Goldschlägerei Wilh. Herbst GmbH & Co. KG Wilhelm-Herbst-Str. 1

28359 Bremen
Telephone no. +49/ 42

Telephone no. +49/ 421/ 2028 - 0 Fax no. +49/ 421/ 2028 - 115 e-mail msds@bego.com

#### Information provided by / telephone

Research & Development Department - Materials, alloys and ceramics; +49/ 421/ 2028 - 130 (Chief Development Officer alloys)

#### **Advice on Safety Data Sheet**

msds@bego.com

#### 1.4 Emergency telephone number

For medical advice (in German and English): +49 (0)551 192 40 (Giftinformationszentrum Nord)

# **SECTION 2: Hazards identification**

#### 2.1 Classification of the substance or mixture

#### **Classification information**

This product does not meet the classification criteria given in the Regulation (EC) No 1272/2008 (CLP).

This product does not meet the labelling criteria given in the Directives 67/548/EWG respectively 1999/45/EC (DPD)

#### 2.2 Label elements

# Labelling according to Regulation (EC) No 1272/2008 (CLP Regulation)

# **Hazard pictograms**

-

# Signal word

-

#### **Hazard statements**

-

#### Hazard statements (EU)

EUH208 Contains TRIISOBUTYLPHOSPHAT. May produce an allergic reaction.

EUH210 Safety data sheet available on request.

#### Labelling information

The labelling (EU hazard statements) meets the criteria of annex II of Directive (EC) Nr. 1272/2008 (CLP).

## 2.3 Other hazards

PBT assessment

No data available.



Trade name: Wiropaint plus (REF 51100)

Current version: 2.0.3, issued: 19.10.2015 Replaced version: 2.0.2, issued: 16.12.2014 Region: GB

vPvB assessment No data available.

# **SECTION 3: Composition/information on ingredients**

#### 3.1 Substances

Not applicable. The product is not a substance.

#### 3.2 Mixtures

**Hazardous ingredients** 

No	o Substance name Addition		onal informati	nal information		
	CAS / EC / Index / REACH no	Classification (EC) 1272/2008 (CLP)	Conce	entration		%
1	cristobalite					
	14464-46-1	STOT RE 2; H373i	>	40.00 - <	70.00	%-b.w.
	238-455-4					
	-					
	-					
2	Silicic acid, sodium					
	1344-09-8	Eye Irrit. 2; H319	>	5.00 - <	20.00	%-b.w.
	215-687-4	Skin Irrit. 2; H315				
	-	STOT SE 3; H335				
	-					
3	TRIISOBUTYLPHO	SPHAT				
	126-71-6	Skin Sens. 1; H317		<	1.00	%-b.w.
	204-798-3					
	-					
	-					
4	titanium-dioxide					
	13463-67-7	-	>	5.00 - <	10.00	%-b.w.
	236-675-5					
	-					
	01-2119489379-17					

Full Text for all H-phrases and EUH-phrases: pls. see section 16

No	Route, target organ, concrete effect
1	H373i
	inhalational; -; -

## 3.3 Other information

According to the manufacturer's information, all harmful ingredients are physically bound in this mixture. Therefore they do not represent a potential hazard to humans and thus it is not subject to mandatory labelling.

# **SECTION 4: First aid measures**

# 4.1 Description of first aid measures

#### **General information**

If the patient is likely to become unconscious, place and transport in stable sideways position. In case of persisting adverse effects, consult a physician. Remove contaminated clothing and shoes immediately, and launder thoroughly before reusing.

#### After inhalation

Remove affected person from the immediate area. Ensure supply of fresh air. Irregular breathing/no breathing: artificial respiration. Call a doctor immediately.

#### After skin contact

In case of contact with skin wash off with water.

#### After eye contact

Remove contact lenses, irrigate copiously with clean, fresh water for at least 15 minutes, holding the eyelids apart and seek medical advice.



Trade name: Wiropaint plus (REF 51100)

Current version: 2.0.3, issued: 19.10.2015 Replaced version: 2.0.2, issued: 16.12.2014 Region: GB

#### After ingestion

Rinse mouth thoroughly with water. Never give anything by mouth to an unconscious person. Call a doctor immediately.

#### 4.2 Most important symptoms and effects, both acute and delayed

No data available.

#### 4.3 Indication of any immediate medical attention and special treatment needed

No data available.

# **SECTION 5: Firefighting measures**

#### 5.1 Extinguishing media

#### Suitable extinguishing media

Carbon dioxide; Extinguishing powder; Water spray jet; Foam

#### Unsuitable extinguishing media

High power water jet

#### 5.2 Special hazards arising from the substance or mixture

In the event of fire, the following can be released: Carbon monoxide and carbon dioxide

#### 5.3 Advice for firefighters

Product itself does not burn. Adapt extinguisher and fire-fighting measures to fire in the environment. Use self-contained breathing apparatus. Wear protective clothing.

# **SECTION 6: Accidental release measures**

#### 6.1 Personal precautions, protective equipment and emergency procedures

# For non-emergency personnel

Ensure adequate ventilation. Refer to protective measures listed in sections 7 and 8.

#### For emergency responders

No data available. Personal protective equipment (PPE) - see Section 8.

# 6.2 Environmental precautions

Do not discharge into the drains/surface waters/groundwater. Do not discharge into the subsoil/soil.

#### 6.3 Methods and material for containment and cleaning up

Pick up with absorbent material (e.g., sand, kieselguhr, acid binder, universal binder, sawdust). When picked up, treat material as prescribed under heading "Disposal considerations".

#### 6.4 Reference to other sections

No data available.

# **SECTION 7: Handling and storage**

## 7.1 Precautions for safe handling

#### Advice on safe handling

Keep container tightly closed.

#### General protective and hygiene measures

Wash hands before breaks and after work. Do not eat, drink or smoke during work time. Remove soiled or soaked clothing immediately. Keep away from foodstuffs and beverages.

#### Advice on protection against fire and explosion

No special measures necessary.

#### 7.2 Conditions for safe storage, including any incompatibilities

#### Technical measures and storage conditions

Keep container tightly closed and dry in a cool, well-ventilated place.

#### Requirements for storage rooms and vessels

Keep only in the original container. Containers which are opened must be carefully closed and kept upright to prevent leakage.



Trade name: Wiropaint plus (REF 51100)

Current version: 2.0.3, issued: 19.10.2015 Replaced version: 2.0.2, issued: 16.12.2014 Region: GB

#### Advice on storage assembly

Protect against humidity; Do not store together with: explosive substances

#### 7.3 Specific end use(s)

No data available.

# **SECTION 8: Exposure controls/personal protection**

#### 8.1 Control parameters

#### Occupational exposure limit values

No	Substance name	CAS no.		EC no.
1	titanium-dioxide	13463-67-7		236-675-5
	List of approved workplace exposure limits (WELs) / E	H40		
	Titanium dioxide			
	total inhalable dust			
	TWA	10	mg/m³	
	List of approved workplace exposure limits (WELs) / E	H40		
	Titanium dioxide			
	respirable dust	•	•	
	TWA	4	mg/m³	

#### 8.2 Exposure controls

#### Appropriate engineering controls

No data available.

#### Personal protective equipment

#### Respiratory protection

If workplace exposure limits are exceeded, a respiration protection approved for this particular job must be worn. In case of aerosol and mist formation, take appropriate measures for breathing protection in the event workplace threshold values are not specified.

#### Eye / face protection

Safety glasses with side protection shield (EN 166)

## Hand protection

Sufficient protection is given wearing suitable protective gloves checked according to i.e. EN 374, in the event of risk of skin contact with the product . Before use, the protective gloves should be tested in any case for its specific workstation suitability (i.e. mechanical resistance, product compatibility and antistatic properties). Adhere to the manufacturer's instructions and information relating to the use, storage, care and replacement of protective gloves. Protective gloves shall be replaced immediately when physically damaged or worn. Design operations thus to avoid permanent use of protective gloves.

#### Other

Normal chemical work clothing.

# **Environmental exposure controls**

No data available.

# **SECTION 9: Physical and chemical properties**

#### 9.1 Information on basic physical and chemical properties

Form/Colour
liquid to pasty; Suspension
white

Odour	
odourless	



Trade name: Wiropaint plus (REF 51100)

Current version: 2.0.3, issued: 19.10.2015 Replaced version: 2.0.2, issued: 16.12.2014 Region: GB

Odour threshold	
No data available	

 pH value
 10

 Value Concentration
 500 g/l H2O

Boiling point / boiling range
No data available

Melting point / melting range

No data available

Decomposition point / decomposition range

No data available

Flash point

No data available

**Auto-ignition temperature** 

No data available

Oxidising properties

No data available

**Explosive properties** 

No data available

Flammability (solid, gas)

No data available

Lower flammability or explosive limits

No data available

**Upper flammability or explosive limits** 

No data available

Vapour pressure

No data available

Vapour density
No data available

Evaporation rate

No data available

Relative density

No data available

Density

Value 1.8 g/cm³

Solubility in water

No data available

Solubility(ies)

No data available

Partition coefficient: n-octanol/water

No data available

**Viscosity** 

No data available

9.2 Other information

Other information

No data available.



Trade name: Wiropaint plus (REF 51100)

Current version: 2.0.3, issued: 19.10.2015 Replaced version: 2.0.2, issued: 16.12.2014 Region: GB

# **SECTION 10: Stability and reactivity**

#### 10.1 Reactivity

Dangerous reactions are not expected if the product is handled according to its intended use.

#### 10.2 Chemical stability

Stable under recommended storage and handling conditions (See section 7).

#### 10.3 Possibility of hazardous reactions

None, if handled according to order.

#### 10.4 Conditions to avoid

Protect from humid air and water.

#### 10.5 Incompatible materials

None known.

#### 10.6 Hazardous decomposition products

No hazardous decomposition products known.

# **SECTION 11: Toxicological information**

# 11.1 Information on toxicological effects

Acute oral toy	1 - 14

No data available

#### Acute dermal toxicity

No data available

#### Acute inhalational toxicity

No data available

# Skin corrosion/irritation

No data available

#### Serious eye damage/irritation

No data available

# Respiratory or skin sensitisation

No data available

# Germ cell mutagenicity

No data available

# Reproduction toxicity

No data available

#### Carcinogenicity

No data available

# STOT-single exposure

No data available

#### STOT-repeated exposure

No data available

#### **Aspiration hazard**

No data available

### Delayed and immediate effects as well as chronic effects from short and long-term exposure

Contact with the skin and eyes may cause mechanical irritation. Inhalation of dusts may irritate the respiratory tract. Can cause silicosis.



Trade name: Wiropaint plus (REF 51100)

Current version: 2.0.3, issued: 19.10.2015 Replaced version: 2.0.2, issued: 16.12.2014 Region: GB

# **SECTION 12: Ecological information**

#### 12.1 Toxicity

Toxicity to fish (acute
-------------------------

No data available

## Toxicity to fish (chronic)

No data available

#### **Toxicity to Daphnia (acute)**

No data available

#### Toxicity to Daphnia (chronic)

No data available

# Toxicity to algae (acute)

No data available

#### Toxicity to algae (chronic)

No data available

#### **Bacteria toxicity**

No data available

#### 12.2 Persistence and degradability

No data available.

#### 12.3 Bioaccumulative potential

No data available.

## 12.4 Mobility in soil

No data available.

### 12.5 Results of PBT and vPvB assessment

Results of PBT and vPvB assessment		
PBT assessment	No data available.	
vPvB assessment	No data available.	

#### 12.6 Other adverse effects

No data available.

#### 12.7 Other information

## Other information

Do not discharge product unmonitored into the environment.

# **SECTION 13: Disposal considerations**

#### 13.1 Waste treatment methods

#### Product

Allocation of a waste code number, according to the European Waste Catalogue, should be carried out in agreement with the regional waste disposal company.

# **Packaging**

Residuals must be removed from packaging and when emptied completely disposed of in accordance with the regulations for waste removal. Incompletely emptied packaging must be disposed of in the form of disposal specified by the regional disposer.



Trade name: Wiropaint plus (REF 51100)

Current version: 2.0.3, issued: 19.10.2015 Replaced version: 2.0.2, issued: 16.12.2014 Region: GB

# **SECTION 14: Transport information**

#### 14.1 Transport ADR/RID/ADN

The product is not subject to ADR/RID/ADN regulations.

#### 14.2 Transport IMDG

The product is not subject to IMDG regulations.

#### 14.3 Transport ICAO-TI / IATA

The product is not subject to ICAO-TI / IATA regulations.

## 14.4 Other information

No data available.

#### 14.5 Environmental hazards

Information on environmental hazards, if relevant, please see 14.1 - 14.3.

#### 14.6 Special precautions for user

No data available.

#### 14.7 Transport in bulk according to Annex II of Marpol and the IBC Code

Not relevant

# **SECTION 15: Regulatory information**

# 15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture EU regulations

#### Regulation (EC) No 1907/2006 (REACH) Annex XIV (List of substances subject to authorisation)

According to the data available and/or specifications supplied by upstream suppliers, this product does not contain any substances considered as substances requiring authorisation as listed on Annex XIV of the REACH regulation (EC) 1907/2006.

#### REACH candidate list of substances of very high concern (SVHC) for authorisation

According to available data and the information provided by preliminary suppliers, the product does not contain substances that are considered substances meeting the criteria for inclusion in annex XIV (List of Substances Subject to Authorisation) as laid down in Article 57 and article 59 of REACH (EC) 1907/2006.

# Regulation (EC) No 1907/2006 (REACH) Annex XVII: RESTRICTIONS ON THE MANUFACTURE, PLACING ON THE MARKET AND USE OF CERTAIN DANGEROUS SUBSTANCES, PREPARATIONS AND ARTICLES

According to the data available and/or specifications supplied by upstream suppliers, this product does not contain any substances requiring authorisation as listed on Annex XIV of the REACH regulation (EC) 1907/2006.

#### Directive 2012/18/EU on the control of major-accident hazards involving dangerous substances

This product is not subject to Part 1 or 2 of Annex I.

#### Other regulations

Adhere to the national sanitary and occupational safety regulations when using this product.

#### 15.2 Chemical safety assessment

No data available.

#### **SECTION 16: Other information**

## Sources of key data used to compile the data sheet:

Regulation (EC) No 1907/2006 (REACH), 1272/2008 (CLP) as amended in each case.

EC Directives 2000/39/EC, 2006/15/EC, 2009/161/EU

National Threshold Limit Values of the corresponding countries as amended in each case.

Transport regulations according to ADR, RID, IMDG, IATA as amended in each case.

The data sources used to determine physical, toxic and ecotoxic data, are indicated directly in the corresponding chapter.



Trade name: Wiropaint plus (REF 51100)

Current version: 2.0.3, issued: 19.10.2015 Replaced version: 2.0.2, issued: 16.12.2014 Region: GB

# Full text of the H- and EUH- phrases drawn up in sections 2 and 3 (provided not already drawn up in these sections)

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.
H319 Causes serious eye irritation.
H335 May cause respiratory irritation.

H373i May cause damage to organs through prolonged or repeated exposure if inhaled.

#### Department issuing safety data sheet

UMCO Umwelt Consult GmbH - D-21107 Hamburg, Georg-Wilhelm-Strasse 183, Tel.: +49(40)79 02 36 300, Fax: +49(40)79 02 36 357, e-mail: umco@umco.de

This information is based on our present state of knowledge and experience.

The security data sheet describes products with a view to the security requirements.

However, it should not constitute a guarantee for any specific product properties and shall not establish a legally valid relationship.

Document protected by copyright. Alterations or reproductions require the express written permission of UMCO Umwelt Consult GmbH.