



EFFEGI BREGA®

- P O L Y M E R I Z A T I O N U N I T -

P O L Y A



Item code: IV800428

TECHNICAL DATA

Chamber capacity	2.0	dm ³
Chamber diameter	123	mm
Chamber height	180	mm
Max water height	160	mm
Temperature range	0÷120	°C
Max. working pressure	0.4 (4)	MPa (bar)
Safety valve set pressure	0.8 (8)	MPa (bar)
Sound pressure level	66.8	dB (A)
Input power	800	W
Voltage	230	V
Frequency	50	Hz
Electric plug	10÷16	A
Depth	280	mm
Length	210	mm
Height	280	mm
Weight	13	kg
Weight packaging included	16	kg
Dimension of packaging	59x29x44	cm

IV800428-00 - EN - EFFEGI BREGA® s.r.l. reserves the right to make changes or amendments without further notice.