



EFFEGI BREGA®

« AURORE »

- ROTARY DRUM DENTAL CLEANER -

The time spent cleaning impression trays in the lab can now be dedicated to other activities thanks to the new "Aurore". By means of a strong jet of cold water, this rotary drum dental cleaner can remove alginates from the impression trays and make them as good as new. Convenient, fast and simple: an appliance that will help you save time and money by returning clean and reusable impression trays to your customers.

Control panel

Functions:

- Start
- Stop
- 3 Leds
- 3 different colours for each led
- Combination of colours and positions to signal, by intuition, the operating conditions of the appliance.

Impression tray cleaner

- Simple: load and press START.
- Easy to use: an electronic card manages all the services and the intuitive led signals.
- Professional: designed for dental offices.



Load hatch

The operator can load up to a maximum of 20 impression trays.

Signals:

The leds can signal the following conditions:

- Appliance ready
- Single cycle
- Double cycle
- Triple cycle
- Upper door not closed
- Lower tray not closed
- Inlet water failure
- Water drainage failure
- Pre-warning: service within 100 hours
- Warning: call for service
- Number of working hours.

Debris collection tray

- To prevent clogging and pollution of the drainage system, the alginates removed from the impression trays are piled up in the collection tray.