



Fact 1

... is a gradiation of color-translucancy without visible layers attainable at all?

Yes, with InfinityLayer® Technology it is!

Fact 2

... do highly-esthetic PMMA's necessarily have to be that expensive?

No, utilizing InfinityLayer® technology we even reduced production cost!

Fact 3

... is there really no PMMA, allowing restorations with a wear time bigger than 1 year?

Sure there is, restorations made of InfinityTemp® are approved for a 3 year wear time for now!

Say Hello to ...

InfinityTemp®





Trust your eyes!

InfinityTemp® the new generation ultra-esthetic PMMA

image – Plaster Dental, Nürnberg



InfinityTemp®

Premium InfinityLayer® PMMA

- InfinityLayer® – Translucency-gradation without visible layers
- Superior esthetic in a heartbeat
- Dead easy to mill and polish
- Approved for 3 year wear time
- Ø 98,5 und Ø 95 – 16mm / 20mm height
- Fair pricing – **MSRP € 85,00** no matter which height



Elasticity Modulus	> 2,3 MPa
Bending Resistance	≥ 120 MPa
Density	≥ 1,17 (g/cm³)
Vickers hardness (0,05/10)	≥ 200
Water absorption	≤ 20 (µg/mm³)
Solubility	≤ 0,8(µg/mm³)
Residual monomer	< 0,5%
Impact Resistance (Charpy)	> 13 KJ/m²
Impact Resistance (Izod)	> 1,3 KJ/m²

Why??

InfinityTemp®

image – Plaster Dental, Nürnberg



... because a beautiful smile tells more than a thousand words!

Want to spread smiles now?

o.mannella@hp-dent.com

+49 7731 38 11 044

www.hp-dent.com

ClearTemp

Premium clear PMMA

- Cristal Clear PMMA
- 1-4-all -> Splints, surgical guides, castables/pressables
- Dead easy to mill and polish
- Ø 98,5 und Ø 95 - 16mm / 20mm height
- Fair pricing - **MSRP € 38,00** no matter which height



Elasticity Modulus	> 2,3 MPa
Bending Resistance	≥ 120 MPa
Density	≥ 1,17 (g/cm³)
Vickers hardness (0,05/10)	≥ 200
Water absorption	≤ 20 (µg/mm³)
Solubility	≤ 0,8(µg/mm³)
Residual monomer	< 0,5%
Impact Resistance (Charpy)	> 13 KJ/m²
Impact Resistance (Izod)	> 1,3 KJ/m²

