

# Safety Data Sheet

according to HPR, Schedule 1

Printing date 07/01/2016


Version number 8

Reviewed on 06/03/2016

## 1 Identification

- **Product identifier**
- **Trade name:** ***AESTHETIC BLUE Polymer***
- **Application of the substance / the mixture** *Denture base material*
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**  
CANDULOR AG  
Boulevard Lilienthal 8  
CH-8152 Glattpark (Opfikon)  
SWITZERLAND
- 
- **Tel.** +41 (0) 44 805 9000  
**Fax** +41 (0) 44 805 9090
- **Information department:**  
Regulatory Affairs  
info@candulor.ch
- **Emergency telephone number:**  
24 Hour Emergency Assistance:  
Emergency-Call USA - Infotrac: 1-800-535-5053  
Emergency-Call Canada - Canutec: 1-613-996-6666
- 
- **General SDS Assistance:**  
US: 1-800-533-6825  
Canada: 1-800-263-8182

## 2 Hazard identification

- **Classification of the substance or mixture**  
*Skin Sensitizer - Category 1 H317 May cause an allergic skin reaction.*
- **Label elements**
- **GHS label elements**  
*The product is classified and labeled according to the Globally Harmonized System (GHS).*
- **Hazard pictograms**  
  
GHS07
- **Signal word** *Warning*
- **Hazard-determining components of labeling:**  
*dibenzoyl peroxide*
- **Hazard statements**  
*May cause an allergic skin reaction.*
- **Precautionary statements**  
*Avoid breathing dust/fume/gas/mist/vapours/spray.  
Wear protective gloves/protective clothing/eye protection/face protection.  
Contaminated work clothing should not be allowed out of the workplace.  
Specific treatment (see on this label).  
IF ON SKIN: Wash with plenty of water.  
If skin irritation or rash occurs: Get medical advice/attention.  
Take off contaminated clothing and wash it before reuse.  
Dispose of contents/container in accordance with local/regional/national/international regulations.*

(Contd. on page 2)

## Safety Data Sheet according to HPR, Schedule 1

Printing date 07/01/2016

Version number 8

Reviewed on 06/03/2016

**Trade name: AESTHETIC BLUE Polymer**

(Contd. of page 1)

- **Classification system:**
- **NFPA ratings (scale 0 - 4)**



Health = 0  
Fire = 1  
Reactivity = 0

- **HMIS-ratings (scale 0 - 4)**



Health = 0  
Fire = 1  
Reactivity = 0

- **Other hazards**

Dust generation

Particular danger of slipping on leaked/spilled product.

- **Results of PBT and vPvB assessment**

· **PBT:** Not applicable.

· **vPvB:** Not applicable.

### 3 Composition/Information on ingredients

- **Chemical characterization: Mixtures**
- **Description:** > 95% Polymethylmethacrylate

- **Dangerous components:**

CAS: 94-36-0 dibenzoyl peroxide

≤1% w/w

### 4 First aid measures

- **Description of first aid measures**
- **General information:** No special measures required.
- **After inhalation:** Supply fresh air and to be sure call for a doctor.
- **After skin contact:** Rinse with water.
- **After eye contact:** Rinse opened eye for several minutes under running water. Then consult a doctor.
- **After swallowing:**  
Rinse out mouth and then drink plenty of water.  
If symptoms persist consult doctor.
- **Information for doctor:**
- **Most important symptoms and effects, both acute and delayed** No further relevant information available.
- **Indication of any immediate medical attention and special treatment needed**  
No further relevant information available.

### 5 Firefighting measures

- **Extinguishing media**
- **Suitable extinguishing agents:**  
CO<sub>2</sub>, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
- **Special hazards arising from the substance or mixture** No further relevant information available.
- **Advice for firefighters**
- **Protective equipment:** No special measures required.

CA

(Contd. on page 3)

## Safety Data Sheet according to HPR, Schedule 1

Printing date 07/01/2016

Version number 8

Reviewed on 06/03/2016

**Trade name: AESTHETIC BLUE Polymer**

(Contd. of page 2)

### 6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures** Not required.
- **Environmental precautions:** Do not allow to enter sewers/ surface or ground water.
- **Methods and material for containment and cleaning up:** Pick up mechanically.
- **Reference to other sections**  
See Section 7 for information on safe handling.  
See Section 8 for information on personal protection equipment.  
See Section 13 for disposal information.

### 7 Handling and storage

- **Handling:**
- **Precautions for safe handling**  
Only adequately trained personnel should handle this product.  
For use in dentistry only.
- **Information about protection against explosions and fires:** No special measures required.
- **Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** Store only in the original receptacle.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:** Protect from heat and direct sunlight.
- **Specific end use(s)** No further relevant information available.

### 8 Exposure controls/ Personal protection

- **Additional information about design of technical systems:** No further data; see item 7.
- **Control parameters**

- **Components with limit values that require monitoring at the workplace:**

**CAS: 94-36-0 dibenzoyl peroxide**

EL Long-term value: 5 mg/m<sup>3</sup>

EV Long-term value: 5 mg/m<sup>3</sup>

- **Additional information:** The lists that were valid during the creation were used as basis.

- **Exposure controls**
- **Personal protective equipment:**
- **General protective and hygienic measures:**  
Usual hygienic measures for dental practice and dental laboratories.  
Keep away from foodstuffs, beverages and feed.  
Wash hands before breaks and at the end of work.  
Avoid contact with the skin.
- **Breathing equipment:** Use respiratory protective device against the effects of dust.
- **Protection of hands:**



Protective gloves

After use of gloves apply skin-cleaning agents and skin cosmetics.

- **Material of gloves**  
Butyl rubber, BR  
Chloroprene rubber, CR

(Contd. on page 4)

CA

## Safety Data Sheet

### according to HPR, Schedule 1

Printing date 07/01/2016

Version number 8

Reviewed on 06/03/2016

**Trade name: AESTHETIC BLUE Polymer**

(Contd. of page 3)

Natural rubber, NR

Nitrile rubber, NBR

Fluorocarbon rubber (Viton)

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

· **Penetration time of glove material**

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

· **Eye protection:** Safety glasses

### 9 Physical and chemical properties

· **Information on basic physical and chemical properties**

· **General Information**

· **Appearance:**

Form: Powder  
Color: According to product specification

· **Odor:** Characteristic

· **Odor threshold:** Not determined.

· **pH-value:** Not applicable.

· **Change in condition**

Melting point/Melting range: 150 °C

Boiling point/Boiling range: Undetermined.

· **Flash point:** 250 °C

· **Flammability (solid, gaseous):** Not determined.

· **Auto igniting:** Product is not selfigniting.

· **Danger of explosion:** Product does not present an explosion hazard.

· **Explosion limits:**

Lower: Not determined.

Upper: Not determined.

· **Vapor pressure:** Not applicable.

· **Density at 20 °C:** 1.2 g/cm<sup>3</sup>

· **Relative density** Not determined.

· **Vapor density** Not applicable.

· **Evaporation rate** Not applicable.

· **Solubility in / Miscibility with**

Water: Nearly insoluble.

· **Partition coefficient (n-octanol/water):** Not determined.

· **Viscosity:**

Dynamic: Not applicable.

Kinematic: Not applicable.

· **Other information** No further relevant information available.

### 10 Stability and reactivity

· **Reactivity** No further relevant information available.

· **Chemical stability** Stable under normal handling and storage conditions.

(Contd. on page 5)

## Safety Data Sheet according to HPR, Schedule 1

Printing date 07/01/2016

Version number 8

Reviewed on 06/03/2016

**Trade name: AESTHETIC BLUE Polymer**

(Contd. of page 4)

- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **Possibility of hazardous reactions** No dangerous reactions known.
- **Conditions to avoid** No further relevant information available.
- **Incompatible materials:** No further relevant information available.
- **Hazardous decomposition products:** None under normal conditions of storage and use.

### 11 Toxicological information

- **Information on toxicological effects**
- **Acute toxicity:**
- **on the skin:** No irritant effect.
- **on the eye:** No irritating effect.
- **Sensitization:** Sensitization possible through skin contact.
- **Additional toxicological information:** No further relevant information available.

- **Carcinogenic categories**

- **IARC (International Agency for Research on Cancer)**

CAS: 9011-14-7	polymethylmethacrylate	3
CAS: 94-36-0	dibenzoyl peroxide	3

- **NTP (National Toxicology Program)**

None of the ingredients is listed.

- **OSHA-Ca (Occupational Safety & Health Administration)**

None of the ingredients is listed.

### 12 Ecological information

- **Toxicity**
- **Aquatic toxicity:** No further relevant information available.
- **Persistence and degradability** No further relevant information available.
- **Behavior in environmental systems:**
- **Bioaccumulative potential** No further relevant information available.
- **Mobility in soil** No further relevant information available.
- **Additional ecological information:**
- **General notes:**
- *Water hazard class 1 (Self-assessment): slightly hazardous for water*
- *Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.*
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.
- **Other adverse effects** No further relevant information available.

### 13 Disposal considerations

- **Waste treatment methods**
- **Recommendation:**
- *Take to an approved landfill or a waste incineration plant, under conditions approved by the local authority.*

(Contd. on page 6)

CA

## Safety Data Sheet

### according to HPR, Schedule 1

Printing date 07/01/2016

Version number 8

Reviewed on 06/03/2016

**Trade name: AESTHETIC BLUE Polymer**

(Contd. of page 5)

- **Uncleaned packagings:**
- **Recommendation:** Disposal must be made according to official regulations.

#### 14 Transport information

- |  |  |                 |
|--|--|-----------------|
| · <b>UN-Number</b>   |  | Void            |
| · <b>DOT, TDG, ADN, IMDG, IATA</b>   |  | Void            |
| · <b>UN proper shipping name</b>   |  | Void            |
| · <b>DOT, TDG, ADN, IMDG, IATA</b>   |  | Void            |
| · <b>Transport hazard class(es)</b>  |  | Void            |
| · <b>DOT, TDG, ADN, IMDG, IATA</b>   |  | Void            |
| · <b>Class</b>   |  | Void            |
| · <b>Packing group</b>   |  | Void            |
| · <b>DOT, TDG, IMDG, IATA</b>  |  | Void            |
| · <b>Environmental hazards:</b>  |  |                 |
| · <b>Marine pollutant:</b>   |  | No              |
| · <b>Special precautions for user</b>  |  | Not applicable. |
| · <b>Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code</b> |  | Not applicable. |
| · <b>Transport/Additional information:</b>                                       | Product is not classified as a dangerous good for transport (ADR, IMDG, IATA). |                 |
| · <b>UN "Model Regulation":</b>  |  | Void            |

#### 15 Regulatory information

- **Safety, health and environmental regulations/legislation specific for the substance or mixture**

- **Sara**

- **Section 355 (extremely hazardous substances):**

None of the ingredients is listed.

- **Section 313 (Specific toxic chemical listings):**

CAS: 94-36-0 dibenzoyl peroxide

- **TSCA (Toxic Substances Control Act):**

All ingredients are listed.

- **Proposition 65**

- **Chemicals known to cause cancer:**

None of the ingredients is listed.

- **Chemicals known to cause reproductive toxicity for females:**

None of the ingredients is listed.

- **Chemicals known to cause reproductive toxicity for males:**

None of the ingredients is listed.

- **Chemicals known to cause developmental toxicity:**

None of the ingredients is listed.

(Contd. on page 7)

## Safety Data Sheet

### according to HPR, Schedule 1

Printing date 07/01/2016

Version number 8

Reviewed on 06/03/2016

**Trade name: AESTHETIC BLUE Polymer**

(Contd. of page 6)

- **Carcinogenic categories**

- **EPA (Environmental Protection Agency)**

None of the ingredients is listed.

- **TLV (Threshold Limit Value established by ACGIH)**

CAS: 94-36-0 dibenzoyl peroxide

A4

- **NIOSH-Ca (National Institute for Occupational Safety and Health)**

None of the ingredients is listed.

- **GHS label elements**

The product is classified and labeled according to the Globally Harmonized System (GHS).

- **Hazard pictograms**



GHS07

- **Signal word** Warning

- **Hazard-determining components of labeling:**

dibenzoyl peroxide

- **Hazard statements**

May cause an allergic skin reaction.

- **Precautionary statements**

Avoid breathing dust/fume/gas/mist/vapours/spray.

Wear protective gloves/protective clothing/eye protection/face protection.

Contaminated work clothing should not be allowed out of the workplace.

Specific treatment (see on this label).

IF ON SKIN: Wash with plenty of water.

If skin irritation or rash occurs: Get medical advice/attention.

Take off contaminated clothing and wash it before reuse.

Dispose of contents/container in accordance with local/regional/national/international regulations.

- **National regulations:**

- **Other regulations, limitations and prohibitive regulations**

The product is a medical device according to the Directive 93/42/EEC.

This product is classified as a medical device under US and Canadian regulations and has been reviewed by the US Food and Drug Administration and Health Canada.

- **Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

## 16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

- **Date of preparation / last revision** 07/01/2016 / 7

- **Abbreviations and acronyms:**

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation

IATA: International Air Transport Association

ACGIH: American Conference of Governmental Industrial Hygienists

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

NFPA: National Fire Protection Association (USA)

HMIS: Hazardous Materials Identification System (USA)

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

(Contd. on page 8)

CA

**Safety Data Sheet**  
**according to HPR, Schedule 1**

Printing date 07/01/2016

Version number 8

Reviewed on 06/03/2016

Trade name: **AESTHETIC BLUE Polymer**

NIOSH: National Institute for Occupational Safety  
OSHA: Occupational Safety & Health  
Skin Sensitizer - Category 1: Sensitisation - Skin, Hazard Category 1  
· **\* Data compared to the previous version altered.**

(Contd. of page 7)

CA