

## NEW POWER steamer descaler

Concentrated descaling solution  
(for POWER steamer 1 and 2)



POWER steamer descaler is a concentrated descaling solution for effective as well as gentle descaling of both POWER steamers.

### Advantages

- Descaling that is both effective and gentle on material – the formulation is especially adapted for the aluminum pressure tank (incl. corrosion protection).
- Considerably faster descaling than with standard commercial descaling solutions.

### Details

- 1 liter of concentrate: practical handling, effective cleaning.
- Can be poured in directly and only requires dilution with water.

### Ordering information

POWER steamer descaler

No. 18450100



## NEW POWER steamer water softener

Water softening system for POWER steamer 2



silent



compact



Reliable water softening and filter system, especially adapted for POWER steamer 2. The system provides optimal protection of the POWER steamer, increasing its efficiency, reliability, and durability.

### Advantages

- Inexpensive and advanced water softening / filter technology that is usually only available in large expensive systems.
- Provides the POWER steamer 2 with clean water with reduced calcium content.
- Easy, compact installation (free-standing or wall-mounted) that can be completed within minutes.
- Simple operation and monitoring (visual and acoustic signal) via a flow sensor unit with an OLED display.

### Details

- Fast replacement of filter cartridges within seconds – remove, insert, that's it.
- Long-lasting filter cartridges.
- Venting and rinsing directly at the filter head, without disconnecting the tubes (via automated inlet and outlet shut-off valves).
- Automatic replacement of all seals between the head and cartridge each time a cartridge is replaced.

### Technical data

Weight (cartridge) ~1,8 kg // ~3.96 lbs

### Ordering information

POWER steamer water softener	No. 18460100
POWER steamer water softener Flowmeter	No. 18460101
POWER steamer water softener Filter head	No. 18460102
POWER steamer water softener Filter cartridge M	No. 18460103
POWER steamer water softener Filter cartridge L	No. 18460104

