

PC minimal requirements

for individual setup

We strongly recommend to use one of the systems provided by CADstar Technology to guarantee full compatibility. These computers will be delivered with all necessary software installed, configured and tested with the system purchased. Keep in mind that components like USB version, drivers, other applications, place of installation ... could cause serious compatibility problems or a reduction of the workflow speed.

min. requirements for CS.Neo2, CS.Neo2 Pro, CS.Neo2 Max

CPU: i7 (min. Gen. 10) or better

Mainstorage: SSD

RAM: 16GB or better

OS: min. Windows 10, 64Bit - licensed

GPU: RTX3060 (min. 12GB VRAM) or better

Monitor: we suggest 24" - full HD resolution

Netzteil: 750W or better (depends on other components)

Connection: min. 1x USB 3.0 or better (2x for CS.Neo2 Max)

CADstar provided workstations



Mainboard: Intel Z690 Chipset

SSD: Samsung 970 Evo Plus 500GB M.2 PCIe 3.0x4 NVMe 1.3

CPU: Intel Core i7-12700 2.50GHZ

GPU: MSI RTX3060 Ventus 2X 12G OC

RAM: 2x DRAM 8GB DDR4 3200MHz

OS: Windows 10 Professional, 64Bit

EQUIPMENT: Basic mouse and keyboard

Warranty: 2 years Bring-In basic warranty

POWER: 550W