

PrograScan® PS7



Simultaneous high speed for perfection

The high-end lab scanner PrograScan PS7 scans upper and lower jaw models in one single process: simultaneously, fully automatically and in just 10 seconds.

Technical Data

Camera	2x (4x8Mpix)
Scan tables	2
Auto-Start	yes
Accuracy (ISO 12836)	5 µm
Repeatability	<2 µm
Scanning speed (tooth model)	10 sec
Scanning speed (impression)	45 sec
Texture	Colour
Scan strategy	All-in-one /die in model
Format of output	STL, PLY
Mesh voltage	90 - 264VAC
Mains voltage	90 Wh
Scanner voltage	24V DC +-10%
Type of USB port	USB 3.0 TYP C (5Gbit/s)
Number of USB ports	2
Weight	13.6 kg
Dimensions (W x H x D)	396.0 x 436.5 x 380.5 mm