



HARVARD®

Harvard **MultiChrome**

NEW



EN

Harvard MultiChrome

It's magic.

Shade-adaptive
composite for almost all
tooth shades.

- Hyper-Nano technology, spherical fillers
- Without any color pigments
- Light-optical chameleon effect for perfect match with natural tooth substance
- Simple storage logistics



The Next


Harvard MultiChrome

Shade-adaptive composite.


The innovative high-tech composite **Harvard MultiChrome** continuously absorbs the colour spectrum of the tooth to be treated and matches it to the existing tooth shade substance thanks to its special Hyper-Nano technology. The secret is the composition, shape and distribution of the fillers inside the material.

The pleasant consistency, excellent polishability and permanent high gloss complement the excellent physical properties of **Harvard MultiChrome**.

Harvard MultiChrome with round Hyper-Nano fillers. White light is refracted and reflected in cavity color.



Results that shine well.




“ Really great! ”


As a practitioner who places a composite filling, I need a material that is nice to handle and has a very smooth consistency. It should not be sticky on the instrument, have a good adaptation to the cavity walls and be thixotropic, it should stay where I want it. And it should be easy to polish. I find all these properties in Harvard MultiChrome. And on top of that, the big plus: with one syringe I can hit almost all tooth shades perfectly. Really great.




Amir Ibrik
D.D.S, N.D.B, E.R.B, MSc




■ Harvard MultiChrome
■ Tokuyama Omnichroma



■ Harvard MultiChrome
■ Tokuyama Omnichroma



■ Harvard MultiChrome
■ Tokuyama Omnichroma



■ Harvard MultiChrome
■ Tokuyama Omnichroma

Generation Composite.

Incredibly adaptable!

One composite – one shade fits all. Modelable and Flow.



Cavity – before



A1

A2

A3

A3.5

A4

B1

B2



Restoration – after

Generation Composite.





B3

B4

C1

C2

C3

C4

D3



Havard MultiChrome is a high-tech composite with active shade matching. Thanks to Hyper-Nano technology, **Havard MultiChrome** does not require shade determination. It is a good feeling for the user to know that he is always getting the right shade result.

Havard MultiChrome continuously absorbs the color spectrum of the tooth to be treated and returns it through its special light-optical effect.

- **Multi-Shade Restorative – shade-adaptive composite**
- Additionally available composite for blocking and masking

Harvard MultiChrome

Multi-Shade Restorative

Shade-adaptive composite

- For user who only want to use one single composite
- Shade-adaptive composite for almost all tooth shades
- Invisible transition between enamel, dentin and composite, creating margins that simply disappear
- Hyper-Nano technology, spherical fillers
- Without any color pigments
- Light-optical chameleon effect
- Superior handling properties
- Low shrinkage
- Permanent high gloss after polishing
- High abrasion resistance
- Excellent physical properties
- Radiopaque
- Opalescence
- Simple storage logistics


Harvard Block 'n Mask

Composite for blocking and masking

- Corrects shade matching in difficult situations (e.g. when the dark oral cavity shines through)



Optimal shade-adaptation for almost all tooth shades after light polymerization.



Posterior restoration A3

Restoration with
Harvard MultiChrome

Restoration of different tooth shades with Harvard MultiChrome

- Multi-Shade Restorative – shade-adaptive flowable composite
- Additionally available flowable composite for blocking and masking

Harvard MultiChrome Flow

Multi-Shade Flowable Restorative

Shade adaptive flowable composite

- For user who only want to use one single composite
- Shade-adaptive composite for almost all tooth shades
- Invisible transition between enamel, dentin and composite, creating margins that simply disappear
- Also for extended fissure sealing
- Hyper-Nano technology, spherical fillers
- Without any color pigments
- Light-optical chameleon effect
- Superior handling properties
- Optimal flowable consistency, thixotropic
- Low shrinkage
- Permanent high gloss after polishing
- High abrasion resistance
- Excellent physical properties
- Radiopaque
- Opalescence
- Simple storage logistics

Harvard Block 'n Mask Flow

Flowable composite for blocking and masking

- Corrects shade matching in difficult situations, severe discolorations (e.g. when the dark oral cavity shines through)



Optimal shade-adaptation for small enamel restorations.



before



after



Extended posterior
fissure sealing

Restoration with
Harvard MultiChrome Flow

Harvard MultiChrome at a glance:



Order information

Harvard MultiChrome		Harvard MultiChrome Flow	
3 g syringe	7100001	2 x 1 ml syringe, incl. 8 needle tips H18	7100003
Harvard MultiChrome		Harvard Block 'n Mask	
20 x 0.25 g OptiTips®	7100002	2 x 1 ml syringe, incl. 8 needle tips H18	7100006
Harvard Block 'n Mask		Harvard NeedleTips H18	7095158
3 g syringe	7100005	Refill bag with 50 needle tips	
Harvard Applier OptiTips®	7095200		
Kit			
Harvard MultiChrome Kit			
3 syringes e.a. 3 g Harvard MultiChrome, 1 syringe 3 g Harvard Block 'n Mask	7100009		

Liability is excluded for all printing errors and omissions. Before using our Harvard products, the respective directions for use should be noticed in every case. All measurements are internal measurements of Harvard Dental International.

Harvard Distribution Partner.

 **HARVARD®**
Marke und Qualität seit 1892

Harvard Dental International GmbH
Margaretenstr. 2 - 4, 15366 Hoppegarten, Germany
Phone: + 49 (0) 30/99 28 978-0
Fax: + 49 (0) 30/99 28 978-19
info@harvard-dental-international.de

www.harvard-dental.de