


TECHNICAL DATA SHEET

BEGO Semados® RS/RXS-Implants 3.0

BEGO Semados® RS/RXS^{Pro}-Implants 3.75–5.5



ON THE RIGHT TRACK



BEGO SEMADOS® IMPLANTS

QUICK & FUNCTIONAL · UNIVERSAL & SUCCESSFUL
RELIABLE & PROVEN · VARIED & SOPHISTICATED
100% GERMAN QUALITY PRODUCT
SIMPLE & INTUITIVE · QUICK & FUNCTIONAL
UNIVERSAL & SUCCESSFUL · RELIABLE & PROVEN
VARIED & SOPHISTICATED · SIMPLE & INTUITIVE

The right product for many indications in implant dentistry

BEGO Semados® RS/RXS or RS/RXS^{Pro} implants are individual and are available with machined shoulder (RS/RXS^{Pro}) or microstructured shoulder (RSX/RXS^{Pro}) and platform switching, depending on the clinician's preference and the requirements of the patient treatment. The range is supplemented by individual CAD/CAM-fabricated prosthetic components and surgical guides. The exceptionally high technical standards ensure a high degree of satisfaction—for implant dentists and patients in equal measure.

Your benefits

- Low stresses in the implant—self-tapping bionic thread designed to withstand stresses
- Integrated platform switching supports the preservation of crestal bone
- Fast—just a few rotations needed when inserting
- Simple to use—insertion with self-centering
- Well-thought out—cover screw included
- Optional—guided/template-guided preparation of the implants with diameters 3.0–4.5 can be carried out using the BEGO Guide system
- For the most stringent quality requirements—100% German development and production

Technical details

- Conical implant design with rounded apex
- Self-tapping thread design with optimal cutting angle
- Diameters: 3.0 / 3.75 / 4.1 / 4.5 / 5.5 mm
- Lengths: 7 / 8.5 / 10 / 11.5 / 13 / 15 mm
- Manufactured from grade 4 pure titanium with ultrapure homogeneous TiPure^{Plus} surface
- Internal taper on the abutments ensures no gaps
- RS/RXS 3.0 is suitable for use in narrow anterior gaps to replace teeth 12, 22, 32–42
- Diameter of 5.5 mm in length 7 mm, for cases with a relatively low bone height

Partners in Progress



PRODUCT OVERVIEW

Distributor contact data

BEGO Semados® RS/RS^{Pro}-Implants

Art. no.	Color code	Product	Quantity
BEGO Semados® RS-Implants 3.0			
57962	□	RS 3.0 L10	—
57963	□	RS 3.0 L11.5	—
57964	□	RS 3.0 L13	—

BEGO Semados® RS^{Pro}-Implants 3.75

58966	■	RS ^{Pro} 3.75 L8.5	—
58967	■	RS ^{Pro} 3.75 L10	—
58968	■	RS ^{Pro} 3.75 L11.5	—
58969	■	RS ^{Pro} 3.75 L13	—
58970	■	RS ^{Pro} 3.75 L15	—

BEGO Semados® RS^{Pro}-Implants 4.1

58971	■	RS ^{Pro} 4.1 L7	—
58972	■	RS ^{Pro} 4.1 L8.5	—
58973	■	RS ^{Pro} 4.1 L10	—
58974	■	RS ^{Pro} 4.1 L11.5	—
58975	■	RS ^{Pro} 4.1 L13	—
58976	■	RS ^{Pro} 4.1 L15	—

BEGO Semados® RS^{Pro}-Implants 4.5

58977	■	RS ^{Pro} 4.5 L7	—
58978	■	RS ^{Pro} 4.5 L8.5	—
58979	■	RS ^{Pro} 4.5 L10	—
58980	■	RS ^{Pro} 4.5 L11.5	—
58981	■	RS ^{Pro} 4.5 L13	—
58982	■	RS ^{Pro} 4.5 L15	—

BEGO Semados® RS^{Pro}-Implants 5.5

58983	■	RS ^{Pro} 5.5 L7	—
58984	■	RS ^{Pro} 5.5 L8.5	—
58985	■	RS ^{Pro} 5.5 L10	—
58986	■	RS ^{Pro} 5.5 L11.5	—
58987	■	RS ^{Pro} 5.5 L13	—

BEGO Semados® RSX/RSX^{Pro}-Implants

Art. no.	Color code	Product	Quantity
BEGO Semados® RSX-Implants 3.0			
57862	□	RSX 3.0 L10	—
57863	□	RSX 3.0 L11.5	—
57864	□	RSX 3.0 L13	—

BEGO Semados® RSX^{Pro}-Implants 3.75

58866	■	RSX ^{Pro} 3.75 L8.5	—
58867	■	RSX ^{Pro} 3.75 L10	—
58868	■	RSX ^{Pro} 3.75 L11.5	—
58869	■	RSX ^{Pro} 3.75 L13	—
58870	■	RSX ^{Pro} 3.75 L15	—

BEGO Semados® RSX^{Pro}-Implants 4.1

58871	■	RSX ^{Pro} 4.1 L7	—
58872	■	RSX ^{Pro} 4.1 L8.5	—
58873	■	RSX ^{Pro} 4.1 L10	—
58874	■	RSX ^{Pro} 4.1 L11.5	—
58875	■	RSX ^{Pro} 4.1 L13	—
58876	■	RSX ^{Pro} 4.1 L15	—

BEGO Semados® RSX^{Pro}-Implants 4.5

58877	■	RSX ^{Pro} 4.5 L7	—
58878	■	RSX ^{Pro} 4.5 L8.5	—
58879	■	RSX ^{Pro} 4.5 L10	—
58880	■	RSX ^{Pro} 4.5 L11.5	—
58881	■	RSX ^{Pro} 4.5 L13	—
58882	■	RSX ^{Pro} 4.5 L15	—

BEGO Semados® RSX^{Pro}-Implants 5.5

58883	■	RSX ^{Pro} 5.5 L7	—
58884	■	RSX ^{Pro} 5.5 L8.5	—
58885	■	RSX ^{Pro} 5.5 L10	—
58886	■	RSX ^{Pro} 5.5 L11.5	—
58887	■	RSX ^{Pro} 5.5 L13	—