


PREMIUM QUALITY

MADE IN GERMANY

5TH GENERATION
ZIRCONIA



CopraSupreme ML PRO

Discover yourself the

Smart Gradient Technology

available from November 4, 2024




without
loss of
translucency

Super-Speedsintering**

60 min* (up to 3 units)

** Super-Speedsintering under vacuum
in the Ivoclar© Programat CS6



Speedsintering*

2 hrs. (up to 4 units)

The Allrounder among Zirconia

The *Smart Gradient Technology* enables a stepless color gradient. The new CopraSupreme ML PRO is a further development of our successful CopraSupreme Symphony product range.

Isostatically pressed with a lifetime warranty.

* can lead to reduced bending strength

	CopraSupreme ML PRO
Zirconia	5th generation, 100% Tosoh
Isostatically pressed	yes
Bending strength	~ 1200 MPa
Class	Ila MDR
Aesthetics	layerless transition multilayer
Translucency	up to 45%
Veneering	all porcelains suitable for zirconia
Colours	A1-D4 and 0M1- 0M3
Heights	12mm, 14mm, 16mm, 18mm, 20mm, 25mm, 30mm in Ø 98 and Ø 95



Areas of application



veneer



inlay



anterior crown



posterior crown



full arch
crown bridge



full crown
bridge