

01 PC Configuration

Recommended PC Configuration		
	Desktop	Laptop
CPU	12th Generation Intel® Core™ i9-12900K AMD Ryzen™ 9 5900X	13th Generation Intel® Core™ i9-13900HX AMD Ryzen™ 9 7845HX
RAM	32 GB (For complex workflows or full-arch implant cases, 64 GB is recommended)	
Graphics*	NVIDIA GeForce RTX 3060 Ti 8GB	NVIDIA GeForce RTX 4060 Laptop 8GB
Disk	1 TB SSD	
OS	Win 10 (build 18362+) / Win 11, 64 bit	

Minimum PC Configuration		
	Desktop	Laptop
CPU	12th Generation Intel® Core™ i5-12600K AMD Ryzen™ 7 5700X	13th Generation Intel® Core™ i7-13700H Processor AMD Ryzen™ 7 7840H
RAM	32 GB	
Graphics*	NVIDIA GeForce RTX 2060 6GB	NVIDIA GeForce RTX 3060 Laptop 6GB
Disk	1 TB SSD	
OS	Win 10 (build 18362+) / Win 11, 64 bit	