



	Cercon <sup>®</sup> xt ML (5Y)	Cercon <sup>®</sup> yo ML (Hybrid 4Y/5Y)	Cercon <sup>®</sup> ht ML (Hybrid 3Y/5Y)
<b>STRENGTH</b>	750 MPa	1000 MPa (dentin) 900 MPa (incisal)	1200 MPa (dentin) 750 MPa (incisal)
<b>TRANSLUCENCY</b>	49%	49% (incisal) 45% (dentin)	49% (incisal) 41% (dentin)
<b>ADVANTAGES</b>	Xtra translucency for high esthetic anterior cases	Optimal balance between strength/translucency for all zirconia cases	High translucency for large posterior cases
<b>SINTERING</b>	1500°C; 4 sintering programs (90min, 3:00, 5:00 and 9:00h)	1500°C; 1550°C; 6 sintering programs (45min, 59min, 90min, 3:00, 5:00 and 9:00h)	1500°C; 4 sintering programs (90min, 3:00, 5:00 and 9:00h)
<b>AVAILABLE COLORS</b>	A1 – D4, BL2	A1 – D4, BL2, BL3	A1 – D4, BL2
<b>AVAILABLE SIZES</b>	14mm, 18mm, 25mm	14mm, 18mm, 25mm	14mm, 18mm, 25mm
	<b>High esthetic anterior cases</b>	<b>High esthetic + high bending strenght</b>	<b>Esthetic large span posterior cases</b>